

CALIBRATION PLASMA NORMAL

COAGULATION CALIBRATION PLASMA (LYOPHILISED)



CALIBRATION PLASMA NORMAL 10 x 1 ml Cat. No K-720

Calibration plasmas are processed from human plasma collected with <0.1% sodium citrate anticoagulant.

RECONSTITUTION

Open the vial carefully and add exactly 1 ml of distilled water, close bottle and let stand undisturbed for 15 minutes at room temperature. Do not invert vial or mix vigorously. Use only after complete reconstitution.

STABILITY AND STORAGE

1. Lyophilised calibration plasmas remain stable when stored at -25°C until expiry date given on the product label.
2. The reconstituted plasmas are stable for 2 hours when stored in capped vial at 2-8°C.

PROCEDURE

Use the CALIBRATION PLASMA NORMAL to establish standard curves in the performance of Factor assays or as a quality control check on standard curves established by other means. Use in accordance with instructions accompanying instrumentation and reagents for the performance of each coagulation Factor.

EXPECTED VALUES

Assay values are provided for each lot. However, actual results depend largely on the standard, reagents, instrument, and individual laboratory protocols in use. Each laboratory should evaluate these products with their reagent system and methodology to establish acceptable limits.

NOTES

- Product for in vitro diagnostic use only.
- The reagents must be used only for the purpose intended by suitably qualified laboratory personnel, under appropriate laboratory conditions.
- Each unit of source material used in the preparation of this product has been tested and found non-reactive for HBsAg and negative for antibodies to HIV and HCV. This product should be handled as potentially infectious biological material.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Thelin, G.M.: Preparation and Standardization of a Stable AHF Plasma, 1968. Thromb. Diath. Haemorrh. 19:423.

Date of issue: 01. 2008.

MANUFACTURER

PZ CORMAY S.A.
ul. Wiosenna 22,
05-092 Łomianki, POLAND
tel.: +48 (0) 22 751 79 10
fax: +48 (0) 22 751 79 14
<http://www.pzcormay.pl>

01/08/01/08