ORGENTEC Diagnostika GmbH

Carl-Zeiss-Straße 49-51

55129 Mainz - Germany

Phone: +49 (0) 61 31 / 92 58-0

Fax: +49 (0) 61 31 / 92 58-58

Internet: www.orgentec.com

Instruction For Use

2016-07

ORG 286 Alpha-1-Antitrypsin

INTENDED USE AND CLINICAL INDICATION

Alegria® Alpha-1-Antitrypsin is an ELISA-based test system for the quantitative measurement of Alpha-1-Antitrypsin in human stool samples. This product is intended for professional in vitro diagnostic use only.

Determination of the alpha-1-antitrypsin concentration in stool is used to evaluate the activity level of chronic inflammatory bowel diseases. Alpha-1-Antitrypsin can be found in stool during intestinal protein leakage, or increased permeability of the intestinal mucosa. Increased values can also occur during celiac disease, polyps, diverticulitis, colon cancer, food allergies etc.

SYMBOLS USED

| | |
|--|------------------------------------|
| | In vitro diagnostic medical device |
| | Manufacturer |
| | Catalogue number |
| | Sufficient for ... determinations |
| | Batch code |
| | Use by |
| | Temperature limitation |
| | Consult instructions for use |
| | Keep away from sunlight |

| | |
|--|-------------------------|
| | Alegria® Test Strips |
| | Wash Buffer |
| | System Fluid |
| | Stool Extraction Medium |
| | Ready to use |

SUMMARY AND EXPLANATION OF THE TEST

Determination of the alpha-1-antitrypsin concentration in stool is used to evaluate the activity level of chronic inflammatory bowel diseases. In addition, faecal alpha-1-antitrypsin is also a marker for protein leakage or increased permeability of the intestinal mucosa when the mucosal barrier is not intact.

This protein is a primary inhibitor that reversibly forms complexes with serine proteases in the digestive tract and with active inflammatory immune cells. In contrast to the digestive enzymes trypsin and chymotrypsin, the elastase of polymorphonuclear neutrophil granulocytes (PMN-elastase) is a proteolytic enzyme that is released during inflammatory processes. Alpha-1-antitrypsin binds to PMN-elastase and ensures that their proteolytic activity is limited to the inflamed area. Alpha-1-antitrypsin thus plays an important regulatory and anti-inflammatory role, preventing an attack on neighbouring healthy tissue during an inflammatory event. This function makes alpha-1-antitrypsin an indicator for active inflammatory processes.

Excretion of alpha-1-antitrypsin with stool is also a marker for the crossing of serum proteins into the intestinal lumen. Elevated values for alpha-1-antitrypsin in stool are primarily found in cases of ulcerative colitis and Crohn's disease. However, they are also found in cases of enteral protein leakage from other causes, such as Wipple's disease, necrotising enterocolitis, intestinal tuberculosis, celiac disease, intestinal lymphoma, or systemic lupus erythematosus. The levels of alpha-1-antitrypsin may also be elevated in patients with intestinal polyps, colon cancer, diverticulitis, food allergies, and pregnancy.

Alpha-1-antitrypsin is an independent indicator of activity for Crohn's disease. Studies show there is also good correlation with anatomical extent of the disease. Measurement of alpha-1-antitrypsin values in stool are also a good way to monitor the success of treatment, detect oncoming flare-ups, and determine the extend of residual inflammation after surgery.

PRINCIPLE OF THE TEST

Alpha-1-Antitrypsin contained in a stool sample has to be released from the stool with stool extraction medium. Then, extracted Alpha-1-Antitrypsin can be analyzed with the Alegria® Alpha-1-Antitrypsin assay.

The Alegria® assay features barcoded 8-well-microstrips, called Alegria® Test Strips. Each strip is designed for a single determination of one patient sample. The Alegria® Test Strip holds a complete set of reagents. Included are enzyme conjugate, enzyme substrate, sample buffer and a test specific control. Furthermore each strip has two antibody-coated wells which serve as reaction wells for one control and one patient sample.

The determination is based on an indirect enzyme linked immune reaction with the following steps: Alpha-1-Antitrypsin present in control / positive samples binds to the surface of the reaction wells forming an antibody-antigen-complex. After incubation, a first washing step removes unbound and unspecifically bound molecules. Subsequently added enzyme conjugate binds to the immobilized antibody-antigen-complex. After incubation, a second washing step removes unbound enzyme conjugate. Addition of enzyme substrate solution results in hydrolysis and colour development during incubation. The intensity of the blue colour correlates with the concentration of the antibody-antigen-complex and can be measured photometrically at 650 nm. The Alegria® Test Strip is based on the proprietary SMC®-Technology (Sensotronic Memorized Calibration): information about the assay, analysis and evaluation, and the lot-specific expiry date is contained on the barcode printed on each Alegria® Test Strip.

The Alegria® Test Strip can be used with the diagnostic instrument Alegria® - a fully automated Random Access Analyser. By means of SMC®-Technology data encoded on the barcode are transferred from the Alegria® Test Strip to the instrument and the assay is automatically processed and evaluated. The instrument reads the date of expiry and rejects further processing if the Alegria® Test Strip is out of date.

WARNINGS AND PRECAUTIONS

- All reagents of this kit are intended for professional in vitro diagnostic use only.
- The control contains material of human origin. Since no test can guarantee the absence of infective agents in human material, we recommend handling Alegria® Test Strips as potentially infective material.
- Bovine serum albumin (BSA) used in components has been tested for BSE and found negative.
- Avoid contact with the substrate TMB (3,3',5,5'-Tetramethyl-benzidine).
- System fluid contains acid, classification is non-hazardous. Avoid contact with skin.
- Control, sample buffer, stool extraction medium and wash buffer contain sodium azide 0.09% as preservative. This concentration is classified as non-hazardous.
- Enzyme conjugate, control and sample buffer contain ProClin 300 0.05% as preservative. This concentration is classified as non-hazardous.

During handling of all reagents, controls and serum samples observe the existing regulations for laboratory safety


regulations and good laboratory practice:

- First aid measures: In case of skin contact, immediately wash thoroughly with water and soap. Remove contaminated clothing and shoes and wash before reuse. If system fluid comes into contact with skin, wash thoroughly with water. After contact with the eyes carefully rinse the opened eye with running water for at least 10 minutes. Remove contact lenses if this can be done easily. Get medical attention if necessary. Respiratory tract: take person to the fresh air. Swallowing: Rinse the mouth and spit the fluids out. Drink 1 - 2 glasses of water immediately. During spontaneous vomiting hold the head of the casual low with the body in a prone position in order to avoid the penetration of vomit into the air tube.
 - Personal precautions, protective equipment and emergency procedures: Observe laboratory safety regulations. Avoid contact with skin and eyes. Do not swallow. Do not pipette by mouth. Do not eat, drink, smoke or apply makeup in areas where specimens or kit reagents are handled. When spilled, absorb with an inert material and put the spilled material in an appropriate waste disposal.
 - Exposure controls / personal protection: Wear protective gloves of nitril rubber or natural latex. Wear protective glasses. Used according to intended use no dangerous reactions known.
 - Conditions to avoid: Since substrate solution is light-sensitiv. Store Alegria® strips in the dark.
 - For disposal of laboratory waste the national or regional legislation has to be observed.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.

CONTENTS OF THE KIT

▽ 24 ORG 286

ALEGRIA TEST STRIPS

Sufficient for 24 determinations

Alegria® Test Strips are modules of 8 wells each composed of:
Wells 1 + 2: empty and not coated (wells for the sample dilution)
Wells 3 + 4: coated with antibody (reaction wells)

Well 5: **Control**: yellow; containing alpha-1-Antitrypsin, PBS, BSA, detergent, preservative sodium azide 0.09% and ProClin 300 0.05%.

Well 6: **Enzyme Conjugate**; light red; containing alpha1-Antitrypsin antibodies, HRP labelled; PBS, BSA, detergent, preservative ProClin 300 0.05%.

Well 7: **Sample Buffer**: yellow; containing TBS, stabilizing protein, detergent, preservative sodium azide 0.09%.

Well 8: **TMB Substrate**: clear; containing 3,3', 5,5'- Tetramethylbenzidin.

Alpha-1-Antitrypsin antibodies are bound to microwells.

Product code on barcode: **alpha1 Antitr Stool**

WASH

1x 20 ml **Wash Buffer**, containing Tris, detergent, preservative sodium azide 0.09%; 50 x concentrate.

SYSTEM FLUID

1x 2.5 ml **System Fluid**, contains acid; 1000 x concentrate.

EXTRACT

1x 20 ml **Stool Extraction Medium**; containing TBS, stabilizing protein, extraction reagent, preservative sodium azide 0.09%. Ready to use. (**F5126**)

1 Alegria® Instruction for Use: Alegria® Mini-DVD

1 Certificate of Analysis

MATERIALS REQUIRED

- Vortex mixer, rocking shaker, benchtop centrifuge, optional: precision scales
- Pipettes for 10 µl
- Measuring cylinder for 1000 ml and 2500 ml, microtubes
- Distilled or deionized water

AUXILIARY IVD REAGENTS

- F5126 **EXTRACT** Stool Extraction Medium; ready to use, 20 ml.
This kit component is also available separately.

FURTHER ACCESSORIES

- ORG 282 Stool Extraction Tubes; containing 100 tubes.
- SSCO Stool Sample Collector; containing 100 pieces.
These auxiliaries are available separately.

STORAGE AND STABILITY

- Store test kit at 2-8 °C in the dark.
- Do not expose reagents to heat, sun, or strong light during storage and usage.
- Store Alegria® Test Strips sealed and desiccated in the clip bag provided.
- Shelf life of the unopened test kit is 15 months from day of production.
Unopened reagents are stable until expiration of the kit. See labels for individual batch.
- Diluted Wash Buffer and System Fluid are stable for at least 30 days when stored at 2-8 °C.
Once transferred to the reagent container we recommend consumption on the same day.

PROCEDURAL NOTES

- Do not use kit components beyond their expiration dates.
- All materials must be at room temperature (20-28 °C) prior to use.
- To avoid carryover or contamination, change the pipette tip between samples.

SPECIMEN COLLECTION, PREPARATION, STORAGE AND HANDLING

Collecting stool samples

Stool Sample Collector (SSCO) may be used.

- Avoid contamination with toilet water containing disinfectants.
- Keep the collected stool sample no longer than 5 days at 2-8 °C or store at -20 °C no longer than 1 month.

Preparing stool sample

Extraction of Alpha-1-Antitrypsin using Stool Extraction Tube (ORG 282)

For the Alegria® Alpha-1-Antitrypsin assay the stool sample has to be diluted 1:50 in the **EXTRACT** provided. Inhomogeneous samples should be homogenised prior to sampling, e.g. by using an inoculation loop.

For liquid stool samples we recommend to weigh out sample (see below method without extraction tube).

1 Prepare Stool Extraction Tube (ORG 282)

Stool sample and **EXTRACT** have to be at room temperature.

Fill 750 µl **EXTRACT** into the transparent tube (part 1) and lock with the spatula (part 2) with the blue closure thread at the tube.

2 Pick stool sample

Open tube and pull out spatula at the orange closure thread.


Remove the stool sample by picking the spatula into the stool. The notches of the spatula need to be fully covered with stool sample. This results in picking 15 mg of stool.

3 Portion stool sample

Press the orange dosing spatula into the blue cone insert of the tube. The surplus stool is wiped off by the cone insert. **The resulting dilution is 1:50.**

4 Extraction

Vortex 30 sec at 1800 rpm to remove stool sample from the notches of the spatula. Further homogenize for 15 min at top speed on a rocking shaker.

Open the tube at the blue cone insert and transfer the homogenate to a microtube. Centrifuge for 2 min at 3000xg. Transfer clear fluid to another microtube and test for Alpha-1-Antitrypsin **immediately**.
Extracted stool sample may be stored at 2-8°C for up to 5 days, or at -20°C for up to 1 week.

Extraction without Stool Extraction Tube (ORG 282)

- Stool sample and **EXTRACT** have to be at room temperature.

• Weigh out empty tube. Pick stool sample and transfer to tube. Weigh out and determine net weight of sample.

• Fill **EXTRACT** into the tube: 49-times net weight of sample needed to gain a 1:50 dilution

| | |
|-------------------------------|-------------------------------|
| stool sample + EXTRACT | stool sample + EXTRACT |
| 15 mg + 0.75 ml | 60 mg + 2.9 ml |
| 20 mg + 1.0 ml | 70 mg + 3.4 ml |
| 30 mg + 1.5 ml | 80 mg + 3.9 ml |
| 40 mg + 2.0 ml | 90 mg + 4.4 ml |
| 50 mg + 2.5 ml | 100 mg + 4.9 ml |

Close tube. Vortex 30 sec at 1800 rpm. Further homogenize for 15 min at top speed on a rocking shaker. Open tube. Transfer the homogenate to a microtube. Centrifuge for 2 min at 3000xg. Transfer clear fluid to another microtube and test for Alpha-1-Antitrypsin **immediately**.
Extracted stool sample may be stored at 2-8°C for up to 5 days, or at -20°C for up to 1 week.

PREPARATION OF REAGENTS

WASH

Dilute the content of the Wash Buffer concentrate (50x) with distilled or deionized water to a final volume of 1000 ml prior to use. Transfer the diluted Wash Buffer into the instrument reagent container. If only one Alegria run is to be performed on one day we recommend transferring only 500 ml diluted Wash Buffer.

SYSTEM FLUID

Dilute the content of the System Fluid concentrate (1000x) with distilled or deionized water to a final volume of 2500 ml prior to use. Transfer the diluted System Fluid into the instrument reagent container.

ALEGRIA TEST STRIPS

Take the required number of Alegria® Test Strips out of the clip bag and let them reach room temperature (20-28 °C). Do not remove foil covering the empty wells until you are ready to start the assay.

EXTRACT

Ready to use


TEST PROCEDURE

Alegria® Test Strips with SMC® technology are used with the diagnostic instrument Alegria®. Detailed information about operating the instrument can be taken from the Instrument User Manual.

(1) Remove the foil from the empty wells 1 to 4 of the Alegria® Test Strip.

Do not remove foil with printed barcode, covering wells 5 to 8.

- (2) Pipette **10 µl extracted stool sample** at the bottom of well 1.
 (3) Insert the strip into the SysTray.
 (4) Place loaded SysTrays into the correct position in the Alegria® instrument and start run. All further steps will be done automatically. The test run is completed when the instrument starts printing the results.

CALIBRATION

The assay system is calibrated using a Alpha-1-Antitrypsin reference preparation.

CALCULATION OF RESULTS

By means of SMC® Technology (Sensotronic Memorized Calibration), all test data are transferred to the system through individual barcodes on the Alegria® Test Strip. Calculation and interpretation of results will be performed automatically.

PERFORMANCE CHARACTERISTICS

Measuring range

The calculation range of this Alegria® assay is 80 - 3600 µg/g

Expected values

In a normal range study with samples from healthy donors the following ranges have been established with this Alegria® assay: Cut-off 400 µg/g

Interpretation of results

| | |
|-----------------|------------|
| Normal range | < 400 µg/g |
| Elevated values | ≥ 400 µg/g |

Linearity

Samples containing Alpha-1-Antitrypsin were serially diluted in **Stool Extraction Medium** to demonstrate the dynamic range of the assay. Activity for each dilution was calculated by means of SMC® Technology.

| Sample | Factor | Observed | Expected | O/E |
|--------|--------|----------|----------|-----|
| | | [µg/g] | [µg/g] | [%] |
| 1 | 1 | 2443.0 | 2443.0 | 100 |
| | 2 | 1168.8 | 1221.5 | 96 |
| | 4 | 581.5 | 610.8 | 95 |
| | 8 | 275.0 | 305.4 | 95 |
| 2 | 16 | 156.1 | 152.7 | 90 |
| | 1 | 2000.0 | 2000.0 | 100 |
| | 2 | 1067.7 | 1000.0 | 107 |
| | 4 | 544.8 | 500.0 | 109 |
| 3 | 8 | 266.5 | 250.0 | 107 |
| | 1 | 1671.5 | 1671.5 | 100 |
| | 2 | 888.7 | 835.8 | 106 |
| | 4 | 445.2 | 417.9 | 107 |
| | 8 | 244.6 | 208.9 | 117 |

Sensitivity

Lowest detectable concentration: 139 µg/g

Reproducibility

Intra-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 20 determinations in a single run. Results for precision-within-assay are shown in the table below.

Inter-assay precision: Coefficient of variation (CV) was calculated for each of three samples from the results of 2 determinations in 5 different runs. Results for run-to-run precision are shown in the table below.

| Intra-Assay | | |
|-------------|-----------|------|
| Sample | Mean µg/g | % CV |
| 1 | 357.8 | 5.8 |
| 2 | 660.7 | 9.0 |
| 3 | 2777.0 | 6.6 |

| Inter-Assay | | |
|-------------|-----------|------|
| Sample | Mean µg/g | % CV |
| 1 | 433.7 | 5.0 |
| 2 | 2578.4 | 9.7 |
| 3 | 2736.2 | 8.7 |

Interfering substances

The following substances were tested per 15 mg of stool sample and were found to be non-interfering: Vancomycin hydrochlorid 0.201 mg, Ciprofloxacin 0.15 mg, Plevazid 0.006 mg, Azathioprine 0.021 mg, Prednisone 0.003 mg, Pentasa 0.399 mg, Vitamin A 2.4 UI, Vitamin C 0.015 mg, Vitamin E 0.03 mg, Vitamin D3 0.33 UI, Hemoglobin 1.749 mg.

The following microorganisms were tested at 1.1×10^7 CFU per 15 mg of stool sample and were found to be non-interfering: *E. coli*, *Klebsiella pneumoniae*, *Salmonella enterica*, *Shigella flexneri*, *Yersinia enterocolitica*.

Study results

171 samples were tested against a comparative method. The study population contained 60 samples from asymptomatic donors, 38 samples from hospitalized patients without further information and 74 samples from hospitalized patients with bowel disease or other abdominal disorders.

169 of 171 samples showed matching results which concludes in 98.8 % agreement.

| | | Reference method | | |
|---------------------|------|------------------|-----|-----|
| | | pos | neg | |
| ORG 286 | pos | 52 | 1 | 171 |
| | neg | 1 | 117 | |
| | | 53 | 118 | |
| Positive agreement: | 98.1 | % | | |
| Positive agreement: | 99.2 | % | | |
| Overall agreement | 98.8 | % | | |

LIMITATIONS OF THE PROCEDURE

This assay is a diagnostic aid. A definite clinical diagnosis should not be based on the results of a single test, but should be made by the physician after all clinical and laboratory findings have been evaluated concerning the entire clinical picture of the patient. Also every decision for therapy should be taken individually.

The above pathological and normal reference ranges in patient stool samples should be regarded as recommendations only. Each laboratory should establish its own ranges according to ISO 15189 or other applicable laboratory guidelines.


Results above 3600 µg/g measuring range: we recommend to predilute the extracted stool sample in Stool Extraction Medium e.g. 1:5 and retest using 10µl predilution in well 1 as described in section TEST PROCEDURE. Multiply the result by dilution factor e.g. factor 5.

REFERENCES

- Rodriguez-Otero P, Porcher R, Peffault de LR, Contreras M, Bouhnik Y, Xhaard A, Andreoli A, Ribaud P, Kapel N, Janin A, Socie G and Robin M.: "Fecal calprotectin and Alpha-1-Antitrypsin predict severity and response to corticosteroids in gastrointestinal graft-versus-host disease". *Blood* 2012. **119**: 5909-5917.
- Hagen LE, Schechter T, Luk Y, Brodovitch A, Gassas A, Doyle JJ: "High Alpha-1-Antitrypsin clearance predicts severity of gut graft-versus-host disease (GVHD) in children". *Pediatr Transplant*. 2011;15(6):659-663.
- Umar SB and DiBaise JK: "Protein-losing enteropathy: case illustrations and clinical review". *Am J Gastroenterol* 2010. **105**: 43-49.
- Braamskamp MJ, Dolman KM, and Tabbers MM, Clinical practice: "Protein-losing enteropathy in children". *Eur J Pediatr* 2010. **169**: 1179-1185.
- Biancone L, Onali S, Calabrese E, Petruzzello C, Zorzi F, Condino G, Sica GS, and Pallone F: "Non-invasive techniques for assessing postoperative recurrence in Crohn's disease". *Dig Liver Dis* 2008. **40 Suppl 2**: S265-S270.
- Faust D, Raschke K, Hormann S, Milovic V, Stein J: "Regulation of alpha-1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells". *Clin Exp Immunol*. 2002 May; **128**(2):279-84.
- Faust D, Spirchez Z, Armbruster FP, and Stein J: "Determination of alpha-1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line". *Z Gastroenterol* 2001. **39**: 769-774.
- Becker K, Frieling T, and Haussinger D: "Quantification of fecal Alpha-1-Antitrypsin excretion for assessment of inflammatory bowel diseases". *Eur J Med Res* 1998. **3**: 65-70.
- Choudhary S, Gibson PR, Deacon MC, Young GP: "Measurement of faecal Alpha-1-Antitrypsin: methodologies and clinical application". *J Gastroenterol Hepatol*. 1996 Apr;11(4):311-8.
- Arndt B, Schürmann G, Betzler M, Herfarth C, Schmidt-Gayk H: "Assessment of Crohn's disease activity and alpha 1-antitrypsin in faeces". *Lancet*. 1992 Oct 24, **340**(8826):1037.
- Boege F and Fischbach W: "Structural heterogeneity of faecal Alpha-1-Antitrypsin shown by immunoblot

analysis in patients with Crohn's disease". *Gut* 1991. **32**: 496-499.

- Strygler B, Nicar MJ, Santangelo WC, Porter JL, Fordtran JS: "Alpha-1-Antitrypsin excretion in stool in normal subjects and in patients with gastrointestinal disorders". *Gastroenterology*. 1990 Nov; **99**(5):1380-7.
- Karbach U, Ewe K: "Enteric protein loss in various gastrointestinal diseases determined by intestinal Alpha-1-Antitrypsin clearance". *Z Gastroenterol*. 1989 Jul; **27**(7):362-5.
- Bjarnason I. Gastrointestinal safety of NSAIDs and over-the-counter analgesics. *Int J Clin Pract Suppl*. 2013 (178):37-42.
- Bjarnason I, Zanelli G, Prouse P, Smethurst P, Smith T, Levi S, et al. Blood and protein loss via small-intestinal inflammation induced by non-steroidal anti-inflammatory drugs. *Lancet*. 1987;2(8561):711-4.

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
RÉFÉRENCES


- Rodriguez-Otero P, Porcher R, Peffault de LR, Contreras M, Bouhnik Y, Xhaard A, Andreoli A, Ribaud P, Kapel N, Janin A, Socie G and Robin M: "Fecal calprotectin and Alpha-1-Antitrypsin predict severity and response to corticosteroids in gastrointestinal graft-versus-host disease". *Blood* 2012. **119**: 5909-5917.
- Hagen LE, Schechter T, Luk Y, Brodovitch A, Gassas A, Doyle JJ: "High Alpha-1-Antitrypsin clearance predicts severity of gut graft-versus-host disease (GVHD) in children". *Pediatr Transplant*. 2011;15(6):659-663.
- Umar SB and DiBaise JK: "Protein-losing enteropathy: case illustrations and clinical review". *Am J Gastroenterol* 2010. **105**: 43-49.
- Braamskamp MJ, Dolman KM, and Tabbers MM, Clinical practice: "Protein-losing enteropathy in children". *Eur J Pediatr* 2010. **169**: 1179-1185.
- Biancone L, Onali S, Calabrese E, Petruzzello C, Zorzi F, Condino G, Sica GS, and Pallone F: "Non-invasive techniques for assessing postoperative recurrence in Crohn's disease". *Dig Liver Dis* 2008. **40 Suppl 2**: S265-S270.
- Faust D, Raschke K, Hormann S, Milovic V, Stein J: "Regulation of alpha-1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells". *Clin Exp Immunol*. 2002 May; **128**(2):279-84.
- Faust D, Spirchez Z, Armbruster FP, and Stein J: "Determination of alpha-1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line". *Z Gastroenterol* 2001. **39**: 769-774.
- Becker K, Frieling T, and Haussinger D: "Quantification of fecal Alpha-1-Antitrypsin excretion for assessment of inflammatory bowel diseases". *Eur J Med Res* 1998. **3**: 65-70.
- Choudhary S, Gibson PR, Deacon MC, Young GP: "Measurement of faecal Alpha-1-Antitrypsin: methodologies and clinical application". *J Gastroenterol Hepatol*. 1996 Apr; **11**(4):311-8.
- Arndt B, Schürmann G, Betzler M, Herfarth C, Schmidt-Gayk H: "Assessment of Crohn's disease activity and alpha 1-antitrypsin in faeces". *Lancet*. 1992 Oct 24, **340**(8826):1037.
- Boege F and Fischbach W: "Structural heterogeneity of faecal Alpha-1-Antitrypsin shown by immunoblot analysis in patients with Crohn's disease". *Gut* 1991. **32**: 496-499.
- Strygler B, Nicar MJ, Santangelo WC, Porter JL, Fordtran JS: "Alpha-1-Antitrypsin excretion in stool in normal

- subjects and in patients with gastrointestinal disorders". *Gastroenterology*. 1990 Nov; **99**(5):1380-7.
- Karbach U, Ewe K: "Enteric protein loss in various gastrointestinal diseases determined by intestinal Alpha-1-Antitrypsin clearance". *Z Gastroenterol*. 1989 Jul; **27**(7):362-5.
- Bjarnason I. Gastrointestinal safety of NSAIDs and over-the-counter analgesics. *Int J Clin Pract Suppl*. 2013 (178):37-42.
- Bjarnason I, Zanelli G, Prouse P, Smethurst P, Smith T, Levi S, et al. Blood and protein loss via small-intestinal inflammation induced by non-steroidal anti-inflammatory drugs. *Lancet*. 1987;2(8561):711-4.

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




- Karbach U, Ewe K: "Enteric protein loss in various gastrointestinal diseases determined by intestinal Alpha-1-Antitrypsin clearance". *Z Gastroenterol*. 1989 Jul; **27**(7):362-5.
- Bjarnason I. Gastrointestinal safety of NSAIDs and over-the-counter analgesics. *Int J Clin Pract Suppl*. 2013 (178):37-42.
- Bjarnason I, Zanelli G, Prouse P, Smethurst P, Smith T, Levi S, et al. Blood and protein loss via small-intestinal inflammation induced by non-steroidal anti-inflammatory drugs. *Lancet*. 1987;2(8561):711-4.

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
- Rodriguez-Otero P, Porcher R, Peffault de LR, Contreras M, Bouhnik Y, Xhaard A, Andreoli A, Ribaud P, Kapel N, Janin A, Socie G and Robin M: "Fecal calprotectin and Alpha-1-Antitrypsin predict severity and response to corticosteroids in gastrointestinal graft-versus-host disease". *Blood* 2012. **119**: 5909-5917.
- Hagen LE, Schechter T, Luk Y, Brodovitch A, Gassas A, Doyle JJ: "High Alpha-1-Antitrypsin clearance predicts severity of gut graft-versus-host disease (GVHD) in children". *Pediatr Transplant*. 2011;15(6):659-663.
- Umar SB and DiBaise JK: "Protein-losing enteropathy: case illustrations and clinical review". *Am J Gastroenterol* 2010. **105**: 43-49.
- Braamskamp MJ, Dolman KM, and Tabbers MM, Clinical practice: "Protein-losing enteropathy in children". *Eur J Pediatr* 2010. **169**: 1179-1185.
- Biancone L, Onali S, Calabrese E, Petruzzello C, Zorzi F, Condino G, Sica GS, and Pallone F: "Non-invasive techniques for assessing postoperative recurrence in Crohn's disease". *Dig Liver Dis* 2008. **40 Suppl 2**: S265-S270.
- Faust D, Raschke K, Hormann S, Milovic V, Stein J: "Regulation of alpha-1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells". *Clin Exp Immunol*. 2002 May; **128**(2):279-84.
- Faust D, Spirchez Z, Armbruster FP, and Stein J: "Determination of alpha-1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line". *Z.Gastroenterol* 2001. **39**: 769-774.
- Becker K, Frieling T, and Haussinger D: "Quantification of fecal Alpha-1-Antitrypsin excretion for assessment of inflammatory bowel diseases". *Eur J Med Res* 1998. **3**: 65-70.
- Choudhary S, Gibson PR, Deacon MC, Young GP: "Measurement of faecal Alpha-1-Antitrypsin: methodologies


and clinical application". *J Gastroenterol Hepatol*. 1996 Apr;11(4):311-8.

- Arndt B, Schürmann G, Betzler M, Herfarth C, Schmidt-Gayk H: "Assessment of Crohn's disease activity and alpha 1-antitrypsin in faeces". *Lancet*. 1992 Oct 24, **340**(8826):1037.
- Boege F and Fischbach W: "Structural heterogeneity of faecal Alpha-1-Antitrypsin shown by immunoblot analysis in patients with Crohn's disease". *Gut* 1991. **32**: 496-499.
- Strygler B, Nicar MJ, Santangelo WC, Porter JL, Fordtran JS: "Alpha-1-Antitrypsin excretion in stool in normal subjects and in patients with gastrointestinal disorders". *Gastroenterology*. 1990 Nov; **99**(5):1380-7.
- Karbach U, Ewe K: "Enteric protein loss in various gastrointestinal diseases determined by intestinal Alpha-1-Antitrypsin clearance". *Z Gastroenterol*. 1989 Jul; **27**(7):362-5.
- Bjarnason I. Gastrointestinal safety of NSAIDs and over-the-counter analgesics. *Int J Clin Pract Suppl*. 2013 (178):37-42.
- Bjarnason I, Zanelli G, Prouse P, Smethurst P, Smith T, Levi S, et al. Blood and protein loss via small-intestinal inflammation induced by non-steroidal anti-inflammatory drugs. *Lancet*. 1987;2(8561):711-4.

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




1. Rodriguez-Otero P, Porcher R, Peffault de LR, Contreras M, Bouhnik Y, Xhaard A, Andreoli A, Ribaud P, Kapel N, Janin A, Socie G and Robin M: "Fecal calprotectin and Alpha-1-Antitrypsin predict severity and response to corticosteroids in gastrointestinal graft-versus-host disease". *Blood* 2012. **119**: 5909-5917.
2. Hagen LE, Schechter T, Luk Y, Brodovitch A, Gassas A, Doyle JJ: "High Alpha-1-Antitrypsin clearance predicts severity of gut graft-versus-host disease (GVHD) in children". *Pediatr Transplant*. 2011;15(6):659-663.
3. Umar SB and DiBaise JK: "Protein-losing enteropathy: case illustrations and clinical review". *Am J Gastroenterol* 2010. **105**: 43-49.
4. Braamskamp MJ, Dolman KM, and Tabbers MM, Clinical practice: "Protein-losing enteropathy in children". *Eur J Pediatr* 2010. **169**: 1179-1185.
5. Biancone L, Onali S, Calabrese E, Petruzzello C, Zorzi F, Condino G, Sica GS, and Pallone F: "Non-invasive techniques for assessing postoperative recurrence in Crohn's disease". *Dig Liver Dis* 2008. **40 Suppl 2**: S265-S270.
6. Faust D, Raschke K, Hormann S, Milovic V, Stein J: "Regulation of alpha-1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells". *Clin Exp Immunol*. 2002 May; 128(2):279-84.
7. Faust D, Spirchez Z, Armbruster FP, and Stein J: "Determination of alpha-1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line". *Z Gastroenterol* 2001. **39**: 769-774.
8. Becker K, Frieling T, and Haussinger D: "Quantification of fecal Alpha-1-Antitrypsin excretion for assessment of inflammatory bowel diseases". *Eur J Med Res* 1998. **3**: 65-70.
9. Choudhary S, Gibson PR, Deacon MC, Young GP: "Measurement of faecal Alpha-1-Antitrypsin: methodologies and clinical application". *J Gastroenterol Hepatol*. 1996 Apr;11(4):311-8.
10. Arndt B, Schürmann G, Betzler M, Herfarth C, Schmidt-Gayk H: "Assessment of Crohn's disease activity and alpha 1-antitrypsin in faeces". *Lancet*. 1992 Oct 24, 340(8826):1037.
11. Boege F and Fischbach W: "Structural heterogeneity of faecal Alpha-1-Antitrypsin shown by immunoblot analysis in patients with Crohn's disease". *Gut* 1991. **32**: 496-499.
12. Strygler B, Nicard MJ, Santangelo WC, Porter JL, Fordtran JS: "Alpha-1-Antitrypsin excretion in stool in normal subjects and in patients with gastrointestinal disorders". *Gastroenterology*. 1990 Nov; 99(5):1380-7.
13. Karbach U, Ewe K: "Enteric protein loss in various gastrointestinal diseases determined by intestinal Alpha-1-Antitrypsin clearance". *Z Gastroenterol*. 1989 Jul; 27(7):362-5.
14. Bjarnason I. Gastrointestinal safety of NSAIDs and over-the-counter analgesics. *Int J Clin Pract Suppl*. 2013 (178):37-42.
15. Bjarnason I, Zanelli G, Prouse P, Smethurst P, Smith T, Levi S, et al. Blood and protein loss via small-intestinal inflammation induced by non-steroidal anti-inflammatory drugs. *Lancet*. 1987;2(8561):711-4.

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REFERENCE

1. Rodriguez-Otero P, Porcher R, Peffault de LR, Contreras M, Bouhnik Y, Xhaard A, Andreoli A, Ribaud P, Kapel N, Janin A, Socie G and Robin M: "Fecal calprotectin and Alpha-1-Antitrypsin predict severity and response to corticosteroids in gastrointestinal graft-versus-host disease". *Blood* 2012. **119**: 5909-5917.
2. Hagen LE, Schechter T, Luk Y, Brodovitch A, Gassas A, Doyle JJ: "High Alpha-1-Antitrypsin clearance predicts severity of gut graft-versus-host disease (GVHD) in children". *Pediatr Transplant*. 2011;15(6):659-663.
3. Umar SB and DiBaise JK: "Protein-losing enteropathy: case illustrations and clinical review". *Am J Gastroenterol* 2010. **105**: 43-49.
4. Braamskamp MJ, Dolman KM, and Tabbers MM, Clinical practice: "Protein-losing enteropathy in children". *Eur J Pediatr* 2010. **169**: 1179-1185.
5. Biancone L, Onali S, Calabrese E, Petruzzello C, Zorzi F, Condino G, Sica GS, and Pallone F: "Non-invasive techniques for assessing postoperative recurrence in Crohn's disease". *Dig.Liver Dis* 2008. **40 Suppl 2**: S265-S270.
6. Faust D, Raschke K, Hormann S, Milovic V, Stein J: "Regulation of alpha-1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells". *Clin Exp Immunol*. 2002 May; 128(2):279-84.
7. Faust D, Spirchez Z, Armbruster FP, and Stein J: "Determination of alpha-1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line". *Z.Gastroenterol* 2001. **39**: 769-774.
8. Becker K, Frieling T, and Haussinger D: "Quantification of fecal Alpha-1-Antitrypsin excretion for assessment of inflammatory bowel diseases". *Eur J Med Res* 1998. **3**: 65-70.
9. Choudhary S, Gibson PR, Deacon MC, Young GP: "Measurement of faecal Alpha-1-Antitrypsin: methodologies and clinical application". *J Gastroenterol Hepatol*. 1996 Apr;11(4):311-8.
10. Arndt B, Schürmann G, Betzler M, Herfarth C, Schmidt-Gayk H: "Assessment of Crohn's disease activity and alpha 1-antitrypsin in faeces". *Lancet*. 1992 Oct 24, 340(8826):1037.
11. Boege F and Fischbach W: "Structural heterogeneity of faecal Alpha-1-Antitrypsin shown by immunoblot analysis in patients with Crohn's disease". *Gut* 1991. **32**: 496-499.
12. Strygler B, Nicar MJ, Santangelo WC, Porter JL, Fordtran JS: "Alpha-1-Antitrypsin excretion in stool in normal subjects and in patients with gastrointestinal disorders". *Gastroenterology*. 1990 Nov; 99(5):1380-7.
13. Karbach U, Ewe K: "Enteric protein loss in various gastrointestinal diseases determined by intestinal Alpha-1-Antitrypsin clearance". *Z Gastroenterol*. 1989 Jul; 27(7):362-5.
14. Bjarnason I. Gastrointestinal safety of NSAIDs and over-the-counter analgesics. *Int J Clin Pract Suppl*. 2013 (178):37-42.
15. Bjarnason I, Zanelli G, Prouse P, Smethurst P, Smith T, Levi S, et al. Blood and protein loss via small-intestinal inflammation induced by non-steroidal anti-inflammatory drugs. *Lancet*. 1987;2(8561):711-4.