HbA1c Direct Calibrator 5 Level Series (assayed calibration hemolysate)

Lyophilized calibration hemolysate for the quantitative in vitro determination of Hemoglobin A1c in blood on photometric systems

REF	Cont.
Y04603	1x 1 mL (Level 1) + 4x 1 mL (Level 2 – 5) finished kit
Y04603SV	1x 1 mL (Level 1) + 4x 1 mL (Level 2 – 5) labelled vials

COMPOSITION

The HbA1c Calibrator Level 1 is a liquid buffer with preservative: 0.095% sodium azide. The HbA1c Calibrators Level 2 – 5 are hemolysates prepared from packed human erythrocytes. Stabilizers to maintain hemoglobin in the reduced state for the accurate calibration of the HbA1c procedure.

PREPARATION

Level 1 is ready for use. Level 2 – 5:

1. Open the vial carefully, avoiding any loss of the lyophilized material.
2. Add exactly 1 ml of deionized water (inaccurate reconstitution of the calibrators and error in assay technique can cause erroneous results).
3. Close the vial carefully and gently mix for 10 minutes, or until all material has dissolved, avoiding the formation of foam. Do not shake!
4. Reconstituted calibrators are now ready for use. Do not lyse calibrators.

STABILITY AND STORAGE

Storage:	at 2 – 8 °C
Stability:	up to the expiration date if sealed tightly
Stability after reconstitution:	30 days at 2 – 8 °C 3 months at -20 °C

CLOSE IMMEDIATELY AFTER USE. PROTECT FROM LIGHT AND HEAT. Do not use the product if there is visible evidence of microbial growth in the vial. Improper handling and/or storage of the calibrators can affect results.

SETPOINT VALUE

The setpoint of the calibrators were obtained by assaying representative samples of the entire lot against materials referenced to NGSP values using the DIALAB HbA1c Direct Reagent. For calibrator values according to IFCC, please use values reported for IFCC method.

WARNINGS AND PRECAUTIONS

1. For in Vitro Diagnostic use.
2. Although this product has been tested and found non-reactive for HBsAg, HCV and anti-HIV 1&2, no known test can offer assurance that products derived from human blood will not transmit disease. Therefore all human serum products and patient specimens should be handled in the same manner as an infectious agent.
3. Do not pipette by mouth. Avoid contact with skin and mucous membranes.

WASTE MANAGEMENT

Please refer to local legal requirements.

REFERENCES

1. Trivelli, L.A., Ranney, H.M., and Lai, H.T., New Eng. J. Med. 284,353 (1971).
2. Gonen, B., and Runenstein, A.H., Diabetologia 15, 1 (1978).
3. Gabbay, K.H., Hasty, K., Breslow, J.L., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J.Clin. Endocrinol. Metab. 44, 859 (1977).
4. Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978)

LOT SPECIFIC ASSAY DATA

Hemoglobin A1c calibrator set with zero calibrator and four levels of haemoglobin A1c. Values and expiry date are lot specific.

LOT: 06205		EXP: 2023-06-30
Assigned Values according to NGSP [%]		
Parameter	Level	Assigned Value
HbA1c	Level 1	0.0 %
HbA1c	Level 2	4.2 %
HbA1c	Level 3	6.8 %
HbA1c	Level 4	9.5 %
HbA1c	Level 5	13.2 %
Assigned Values according to IFCC [mmol/mol Hb] *		
Parameter	Level	Assigned Value
HbA1c	Level 1	0 mmol/mol
HbA1c	Level 2	22 mmol/mol
HbA1c	Level 3	51 mmol/mol
HbA1c	Level 4	80 mmol/mol
HbA1c	Level 5	121 mmol/mol
* HbA1c in mmol/mol according to IFCC = (HbA1c in % NGSP – 2.15) / 0.0915		

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