HbA1c Direct Control Set

(en) English

REF

Content

Y04604 4x 1 mL:
 2x 1 mL Level 1
 2x 1 mL Level 2

For professional in vitro diagnostic use only.

INTENDED USE

HbA1c Direct Control Set is a set of lyophilized control hemolysates used for tests for the quantitative in vitro determination of Haemoglobin A1c in blood on photometric systems. For professional laboratory use.

DIAGNOSTIC SIGNIFICANCE

Please refer to the corresponding section of the used HbA1c Direct reagent.

TEST PRINCIPLE

Please refer to the corresponding section of the used Dialab HbA1c Direct reagent.

REAGENT COMPOSITION

HbA1c Direct Control Set is a set of lyophilized Hemolysates prepared from packed human erythrocytes with stabilizers. The controls provide two levels of HbA1c, one level in the normal range and the other level in the elevated range.

MATERIAL REQUIRED BUT NOT PROVIDED

- Deionized water for the preparation of controls and haemolysis of patient samples.
- HbA1c Direct Reagent, e.g.:

REF	Name	Content
Y04602	HbA1c Direct	4x 22.5 mL R1 + 1x 30 mL R2
Y04606	HbA1c Direct	4x 7.5 mL R1 + 1x 10 mL R2
Y04601	HbA1c Direct	4x 45 mL R1 + 1x 60 mL R2

- Calibrator:

REF	Name	Content
Y04603	HbA1c Direct Calibrator 5 Level Series	5 x 1 mL

- Clinical Laboratory Analyser or Photometer
- General laboratory equipment

REAGENT PREPARATION

- Open the vial carefully, avoiding any loss of the lyophilized material.
- Add exactly 1 mL of deionized water (inaccurate reconstitution of the controls and error in assay technique can cause erroneous results).
- Close the vial carefully and gently mix for 10 minutes, or until the material is completely dissolved, avoiding the formation of foam. Do not shake!
- Reconstituted controls are now ready for use. Do not lyse controls.

STORAGE AND STABILITY

Storage:	at 2 – 8 °C
Stability:	up to the expiration date
Stability after reconstitution:	30 days at 2 – 8 °C if sealed tightly 3 months at - 20 °C

CLOSE IMMEDIATELY AFTER USE. PROTECT FROM LIGHT AND HEAT.

Do not use the product in case of suspect of microbial growth in the vial. Improper handling and/or storage of the calibrators and controls can affect results.

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Although this product has been tested and found non-reactive for HBsAg, HCV and anti-HIV 1&2, no known test can offer assurance that products derived from human blood will not transmit disease. Therefore, all human serum products and patient specimens should be handled in the same manner as an infectious agent. Do not pipette by mouth. Avoid contact with skin and mucous membranes.
- Please refer to the safety data sheets (SDS) and take the necessary precautions for the use of laboratory reagents.
- In case of product malfunction or altered appearance that could affect the performance, contact the manufacturer.
- In the event of an incident related to the device, report it to the manufacturer and your competent authority as required.

SPECIMEN COLLECTION AND STORAGE

Not applicable

TEST PROCEDURE

- Please refer to the steps in the section for reagent preparation.
- Measure controls in the normal and elevated range with HbA1c reagent in the same manner, like haemolysed human samples in accordance with the instructions of the reagent and the instrument manual.
- Compare received results with specified ranges defined for each control.

INTERPRETATION OF RESULTS

The assay results with the controls should be within their expectation. The limits in the lot-specific table shown below are to be used as a guideline in determining the accuracy of the assay procedure. Target values according to NGSP in % and IFCC in mmol/mol Hb as well as suggested ranges can be found in the corresponding table below. Assay values may vary slightly with different reagents and methodologies used.

QUALITY CONTROL AND CALIBRATION

Please refer to the instructions for use of the used HbA1c reagent, calibrator and instrument. It is recommended to use quality controls after each calibration, at new measuring sessions, after each change of reagents and every time where a potential increased risk of malfunction can be assumed or when results indicate a potential malfunction and to repeat the calibration where required. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Please refer to the corresponding section of HbA1c Direct reagent package insert.

TRACEABILITY

Values are based on reference material from NGSP/DCCT.

EXPECTED VALUES

For the reagent, please refer to the corresponding section of HbA1c Direct reagent package insert. Lot specific assigned values of the controls are shown in the table below.

LIMITATIONS

The product should not be used after the expiry date.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Trivelli, L.A., Ranney, H.M., and Lai, H.T., New Eng. J. Med. 284,353 (1971).
- Gonen, B., and Runenstein, A.H., Diabetologia 15, 1 (1978).
- Gabbay, K.H., Hasty, K., Breslow, J.L., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J.Clin. Endocrinol. Metab. 44, 859 (1977).
- Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978).

USED SYMBOLS

Symbol	Description
	Content

LOT-SPECIFIC VALUES

Haemoglobin A1c control set with two levels of haemoglobin A1c. Values and expiry date are lot specific. Method: direct immunoturbidimetric latex test.

LOT	07297	2025-04-30
Assigned Values according to NGSP		
Parameter	Assigned Value [%]	
HbA1c Control Level 1, Lot. 07297	5.7 (4.8 – 6.6)	
HbA1c Control Level 2, Lot. 07298	9.8 (8.3 – 11.3)	
Assigned Values according to IFCC *		
Parameter	Assigned Value [mmol/mol Hb]	
HbA1c Control Level 1, Lot. 07297	39 (29 - 49)	
HbA1c Control Level 2, Lot. 07298	84 (67 - 100)	
* HbA1c in mmol/mol according to IFCC = (HbA1c in % NGSP – 2.15) / 0.0915		

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- Kalibrator:

REF	Name	Inhalt
Y04603	HbA1c Direct Calibrator 5 Level Series	5 x 1 mL

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