Lyophilized Calibration Hemolysate

HEMOGLOBIN CALIBRATOR

(Assayed Calibration Hemolysate)

Lyophilized calibration blood for the use in tests for the quantitative and qualitative in vitro determination of total hemoglobin in blood with the cyanmethemoglobin hemoglobin procedure.

REF

Cont.

Y04705SV

1 x 2 ml single vial

CALIBRATOR COMPOSITION

The lyophilized haemoglobin standard is hemolysate prepared from human erythrocytes. The standard provides a level of haemoglobin concentration within the linearity range of the cyanmethemoglobin hemoglobin procedure. Stabilizers are added to maintain haemoglobin in the reduced state.

Reactive ingredients: Hemoglobin prepared from human erythrocytes.

Preservatives added.

CALIBRATOR PREPARATION

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Add exactly 2 ml of deionized water (inaccurate reconstitution of the control and error in assay technique can cause erroneous results).
3. Close the vial carefully and allow to stand for 15 minutes with occasional gentle mixing.

CALIBRATOR STABILITY AND STORAGE

Storage: at 2 – 8°C

Stability: until date of expiration date if sealed tightly.

Stability after reconstitution: at 2-8°C 30 days

CLOSE IMMEDIATELY AFTER USE.

PROTECT FROM LIGHT AND HEAT.

Do not use the product if there is visible evidence of microbial growth in the vial. The presence of obvious turbidity after reaction with cyanmethemoglobin reagent would indicate deterioration, and the product should not be used.

Improper handling and/or storage of the calibrator can affect results.

WARNINGS AND PRECAUTIONS

1. This product is for in vitro diagnostic use only.
2. Handle with the same precautions used for human blood samples. The human blood used in this product has been tested and found to be non-reactive for hepatitis B surface antigen, hepatitis C antigen, and HIV antibody.
3. Do not pipette by mouth. Avoid contact with skin and mucous membranes.

NOTE:

This calibrator should be used to calibrate the determination of unknown specimens using the cyanmethemoglobin assay.

Controls should be included to validate the calibration of the assay.

Limitations: among factors that might cause inaccurate results are improper pipetting, inadequate mixing.

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