G6PDH CONTROL SET

3 levels

Lyophilized control blood for the use in tests for the quantitative and qualitative in vitro determination of Glucose-6-phosphate dehydrogenase (G6PDH) activity in human blood.

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COMPOSITION

Lyophilized human blood preparations with specified activities of Glucose-6-phosphate dehydrogenase.

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CONTROL PREPARATION

- Carefully remove the cap and rubber plug avoiding any loss of contents.
- Add exactly 0.5 mL of deionized water (inaccurate reconstitution of the control can cause erroneous results).
- Close the vial carefully and allow the control to stand for 10 minutes at room temperature
- Swirl the vial carefully and mix thoroughly making sure that all traces of dry material have dissolved. Avoid foaming. Do not shake!
- Again let stand for 5 minutes at room temperature.
- The controls should be treated and assayed in the same manner as patient samples.

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TEST PROCEDURE

Please refer to the reagent package insert for instructions for use.

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MATERIALS REQUIRED BUT NOT PROVIDED

Lyse Reagent, G6PDH reagent.
General laboratory equipment.

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STABILITY AND STORAGE

Storage: at 2 – 6 °C
 Stability: until date of expiration date if sealed tightly.
 Stability after reconstitution: at 2 – 6 °C 10 days
 Close immediately after use. Avoid contamination!

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WARNINGS AND PRECAUTIONS

- For professional use only.
- For in vitro diagnostic use only.
- Handle products with care and precision and with the same precautions used for human blood samples.
- Each individual donation intended for use in manufacture of the G6PDH controls was tested for Hepatitis B surface antigen (HBsAg), Hepatitis C virus (HCV) and anti-HIV by FDA required tests. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human blood should be considered potentially infectious.
- Refer to safety data for additional risk information

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WASTE MANAGEMENT

Please refer to local legal requirements.

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