## HbA1c Direct

Diagnostic reagent for quantitative in vitro determination of HbA1c in human blood by turbidimetric assay

REF	Content
Y04602	4x 22.5 mL R1 + 1x 30 mL R2
Y04606	4x 7.5 mL R1 + 1x 10 mL R2
Y04601	4x 45 mL R1 + 1x 60 mL R2
YA1004	1x 30 mL R1 + 1x 10 mL R2
YT1204	1x 30 mL R1 + 1x 10 mL R2
Y06911	1x 30 mL R1 + 1x 10 mL R2
Y06917	1x 30 mL R1 + 1x 10 mL R2
YE1804	2x 30 mL R1 + 1x 20 mL R2

Additionally offers

Y04605	5x 100 mL	Hemolysis Reagent
Y04603	5x 1 mL	HbA1c Direct Calibrator 5 Level Series
Y04604	4x 1 mL	HbA1c Direct Control Set

### GENERAL INFORMATION

Method	Immunoturbidimetric
Wavelength	660 nm
Assay Temperature	37 °C
Sample	Whole blood with EDTA
Measuring range	approx. 0 – 14 %

### INTENDED USE

The determination of HbA1c is most commonly performed for the evaluation of glycemic control in diabetes mellitus. HbA1c values provide an indication of glucose levels over the preceding 4-8 weeks. A higher HbA1c value indicates poorer glycemic control.

### REAGENT COMPOSITION

<b>Reagent 1</b>	
Latex	
Sodium azide	0.95 g/L
<b>Reagent 2</b>	
Anti-human HbA1c mouse monoclonal antibody	
Stabilizers	

### REAGENT PREPARATION

Reagent 1 and 2 are ready to use. Mix gently before use.

### REAGENT STABILITY AND STORAGE

Conditions:	Protect from light. Close immediately after use. Do not freeze the reagents!
Storage:	at 2 – 8 °C
Stability:	up to the expiration date
R1 and R2 are stable for 30 days after opening if stored at 2 – 8 °C	

### SAMPLE COLLECTION AND PREPARATION

Use fresh EDTA whole blood.  
 To determine HbA1c, a hemolysate must be prepared for each sample:

1. Dispense 1 mL Hemolysis Reagent into a test tube.
2. Add 10 µL of well mixed EDTA whole blood into the tube and mix.
3. Allow to stand for 5 minutes or until complete lysis is evident.

Hemolysate stability: at 2 – 8 °C 72 hours

Discard contaminated specimens.

### GENERAL TEST PROCEDURE

Solve calibrators and controls, do not lyse them. Perform lysis with samples only.

Use the HbA1c Direct Calibrator 5 Level Series to generate a calibration curve.

Pipette into test tubes:	Calibrators	Samples/Controls
HbA1c Reagent 1	180 µL	180 µL
Calibrators/Controls/Samples	5 µL	5 µL
Mix. Incubate for 2 minutes at assay temperature. Then add:		
HbA1c Reagent 2	60 µL	60 µL
Mix and read A1 of calibrators and samples/controls at 660 nm. Incubate for 5 minutes, then read A2 at 660 nm. Calculate: $\Delta A = (A2 - A1)$		

### CALCULATION (light path 1 cm)

HbA1c results according to NGSP for the samples and controls are determined using the prepared calibration curve.

For calculating results according to IFCC, use IFCC calibrator values (see insert calibrator series) or use the following equation:

$$\text{IFCC (mmol/mol)} = (\text{HbA1c in \% NGSP} - 2.15) / 0.0915$$

### REFERENCE RANGE

According to NGSP: < 6 % for a non-diabetic  
 < 7 % for glycemic control of a person with diabetes

According to IFCC: < 42,0 mmol/mol Hb for a non-diabetic  
 < 53,0 mmol/mol Hb for glycemic control of a person with diabetes

It is recommended that each laboratory establishes its own normal range. In using Hemoglobin A1c to monitor diabetic patients, results should be interpreted individually (the patient should be monitored against him or herself). There is a 3-4 week time lag before Hemoglobin A1c reflects changes in blood glucose level.

### TEST PRINCIPLE

This method utilizes the interaction of antigen and antibody to determine directly the HbA1c in whole EDTA blood. HbA1c in test samples is adsorbed onto the surface of latex particles, which react with Anti-HbA1c (antigen-antibody reaction) and gives agglutination. The amount of agglutination is measured as absorbance. The HbA1c value is obtained from a calibration curve.

### DIAGNOSTIC IMPLICATION

Throughout the circulatory life of the red cell, HbA1c is formed continuously by the addition of glucose to the N-terminal end of the haemoglobin beta chain. The process, which is non-enzymatic, reflects the average exposure of haemoglobin to glucose over an extended period. In a classical study, Trivelli et al<sup>1</sup> showed HbA1c in diabetic subjects to be elevated 2-3 folds over the levels found in normal individuals. Several investigators have recommended that HbA1c serves as an indicator of metabolic control of the diabetic, since HbA1c levels approach normal values for diabetics in metabolic control.<sup>2,3,4</sup>  
 HbA1c has been defined operationally as the "fast fraction" haemoglobins (HbA1a, A1b, A1c) that elute first during column chromatography with cation-exchange resins. The non-glycosylated haemoglobin, which consists of the bulk of the haemoglobin has been designated HbA0.

### PERFORMANCE CHARACTERISTICS

#### LINEARITY

The Hemoglobin A1c assay range is 0 % - 14 % (NGSP).

#### PRECISION [%CV]

	Low	Medium	High
Intra-Assay	1.14	0.72	1.22
Inter-Assay	1.85	2.05	2.97

#### ACCURACY [%]

Control	Assigned Value	Measured Value
DIALAB L1	5.8 (4.9 – 6.7)	6.0
DIALAB L2	10.1 (8.6 – 11.6)	10.1
BIORAD L1	5.48 (4.39 – 6.58)	5.81
BIORAD L2	9.42 (7.54 – 11.3)	9.08

#### METHOD COMPARISON

A comparison study of this Hemoglobin A1c procedure and Roche gave the following results:  $y = 0.9286x + 0.26 / r = 0.9656$

#### SENSITIVITY

The change of absorbance units per concentration unit was calculated for the up going part of the nonlinear curve. Result was measured: 0.0458744 Absorbance/Concentration units.

#### INTERFERING SUBSTANCES

No interference of:  
 High Lipids: 3000 formazin turbidity units  
 Free Bilirubin: 30 mg/dL  
 Conjugated Bilirubin: 30 mg/dL  
 The reagent cannot detect HbF therefore elevated levels of HbF can lead to underestimation of HbA1c results. Patients suffering from Thalassemia might show lower results due to short lifespan of HbA1c in samples.

#### QUALITY CONTROL

All control material with % HbA1c values determined by this method can be used. We recommend the DIALAB HbA1c Direct Control Set.

#### CALIBRATION

The assay requires the use of a HbA1c calibrator. We recommend the DIALAB HbA1c Direct Calibrator 5 Level Series.

#### AUTOMATION

Applications for automated systems are available upon request.

#### WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Sodium azide has been reported to form lead or copper azide in laboratory plumbing, which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.
3. Each donor unit used in the preparation of the calibrators and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.
4. This assay should not be used for the diagnosis of diabetes mellitus.

#### WASTE MANAGEMENT

Please refer to local legal requirements.

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