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REF

Content

Z02500CE	- 30 tests individually packed (30 x Ref. No: Z02500B). - 1 package insert
Z02500B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	competitive immunochromatographic assay
Antigen / Antibodies	monoclonal anti-methamphetamine-antibody-coupled particles and methamphetamine-protein conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Cut-off	1000 ng/ml

INTENDED USE

The "DIAQUICK" MET Dipstick (urine) is a lateral flow, one step immunoassay for the qualitative detection of methamphetamine in human urine at a cut-off concentration of 1000 ng/ml. This product is used to obtain a visual, qualitative result and is intended for professional use. This assay provides only a preliminary analytical test result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) has been established as the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Methamphetamine is a potent sympathomimetic agent with therapeutic applications. The drug can be taken orally, injected, or inhaled. Acute higher doses lead to enhanced stimulation of the central nervous system and induce euphoria, alertness, reduced appetite, and a sense of increased energy and power. Cardiovascular responses to Methamphetamine include increased blood pressure and cardiac arrhythmias. More acute responses include anxiety, paranoia, hallucinations, psychotic behavior, and eventually, depression and exhaustion. The effects of Methamphetamine generally last 2-4 hours, and the drug has a half-life of 9-24 hours in the body. Methamphetamine is excreted in the urine primarily as amphetamine and oxidized and deaminated derivatives. However, 10-20% of Methamphetamine is excreted unchanged. Thus, the presence of the parent compound in the urine indicates Methamphetamine use. Methamphetamine is generally detectable in the urine for 3-5 days, depending on urine pH level.

TEST PRINCIPLE

The "DIAQUICK" MET Dipstick (urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Methamphetamine, if present in the urine specimen below 1,000 ng/ml, will not saturate the binding sites of antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized methamphetamine conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the methamphetamine level exceeds 1,000 ng/ml because it will saturate all the binding sites of anti-methamphetamine antibodies.

To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS PROVIDED

- test strips
- package insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded in a proper biohazard container after testing.

STORAGE

The kit can be stored at room temperature or refrigerated (2-30°C). The test strip is stable through the expiration date printed on the sealed pouch. The test strip must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE

Allow the test strip, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Remove the test strip from the sealed pouch and use it as soon as possible.
- With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. The result should be read at 5 minutes.
- Do not interpret the result after 10 minutes.


INTERPRETATION OF RESULTS

NEGATIVE:* Two lines appear. One line should be in the control region (C), and another line should be in the test region (T). This negative result indicates that the methamphetamine concentration is below the detectable level of 1000 ng/ml.

***NOTE:** The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the methamphetamine concentration is above the detectable level of 1000 ng/ml.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

1. This product is designed for use with urine only.
2. Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
3. The test is a qualitative urine assay and is not for determining quantitative concentration levels or the level of intoxication.
4. Adulterants such as bleaching or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" MET Dipstick (urine) and a leading commercially available Methamphetamine rapid test. Testing was performed on 300 clinical specimens. Ten percent of the specimens employed were either at -25% or +25% level of the cut-off concentration of 1,000 ng/ml methamphetamine. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method		Other MET Rapid Test		Total Results
DIAQUICK MET Dipstick	Results	Positive	Negative	
	Positive	145	0	145
	Negative	3	152	155
Total Results		148	152	300
% Agreement		98%	100%	99%

When compared to GC/MS at the cut-off of 1000 ng/ml, the following results were tabulated:

Method		GC/MS		Total Results
DIAQUICK MET Dipstick	Results	Positive	Negative	
	Positive	135	10	145
	Negative	1	154	155
Total Results		136	164	300
% Agreement with GC/MS Analysis		99%	94%	96%

Analytical Sensitivity

A drug-free urine pool was spiked with methamphetamine at the following concentrations: 0 ng/ml, 500 ng/ml, 750 ng/ml, 1,000 ng/ml, 1,250 ng/ml and 1,500 ng/ml. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

MET Concentration (ng/ml)	Percent of cut-off	n	Visual Result	
			Negative	Positive
0	0%	30	30	0
500	-50%	30	30	0
750	-25%	30	30	0
1000	Cut-off	30	18	12
1250	+25%	30	1	29
1500	+50%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" MET Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/ml)
p-Hydroxymethamphetamine	30,000
D-Methamphetamine	1,000
L-Methamphetamine	8,000
(±)-3,4-Methylenedioxyamphetamine	2,000
Mephentermine	50,000

Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing no methamphetamine, 25% methamphetamine above and below the cut-off and 50% methamphetamine above and below the 1000 ng/ml cut-off was provided to each site. The following results were tabulated:

MET conc. (ng/ml)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
500	15	15	0	14	1	13	2
750	15	11	4	10	5	10	5
1250	15	8	7	4	11	6	9
1500	15	1	14	1	14	0	15

Effect of Urinary Specific Gravity

Fifteen urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 500 ng/ml and 1500 ng/ml of methamphetamine respectively. The "DIAQUICK" MET Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with methamphetamine to 500 ng/ml and 1500 ng/ml. The spiked, pH-adjusted urine was tested with the "DIAQUICK" MET Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or methamphetamine positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" MET Dipstick (urine) at a concentration of 100 µg/ml.

Non Cross-Reacting Compounds

Acetaminophen	Estrone-3-sulfate	Oxymetazoline
4-Acetamidophenol	Erythromycin	Papaverine
Acetophenetidin	β-Estradiol	Penicillin-G
N-Acetylprocainamide	Estrone-3-sulfate	Pentobarbital
Acetylsalicylic acid	Ethyl-p-aminobenzoate	Perphenazine
Aminopyrine	Fenfluramine	Phencyclidine
Amitypyline	Fenpropfen	Phenelzine
Amobarbital	Furosemide	Phenobarbital
Amoxicillin	Gentisic acid	Phentermine
Ampicillin	Hemoglobin	L-Phenylephrine
L-Ascorbic acid	Hydralazine	β-Phenylethylamine
D-Amphetamine	Hydrochlorothiazide	Phenylpropanolamine
D,L-Amphetamine	Hydrocodone	Prednisolone
L-Amphetamine	Hydrocortisone	Prednisone
Apomorphine	p-Hydroxyamphetamine	Procaine
Aspartame	O-Hydroxyhippuric acid	Promazine
Atropine	3-Hydroxytyramine	Promethazine
Benzilic acid	Ibuprofen	D, L-Propranolol
Benzolic acid	Imipramine	D-Propoxyphene
Benzoylcholine	Ipronioid	D-Pseudoephedrine
Benzphetamine	(±) - Isoproterenol	Quinacrine
Bilirubin	Isoxsuprine	Quinidine
(±) - Brompheniramine	Meprobamate	Thiamine
Caffeine	Methadone	Thioridazine
Cannabidiol	Methoxyphenamine	D, L-Tyrosine
Chloralhydrate	(+) 3,4-Methylenedioxy-amphetamine	Tolbutamide
Chloralidoxide	3,4-Methylenedioxyethyl-amphetamine	Trans-2-phenyl-cyclopropylamine
Chlorothiazide	Methyphenidate	Triamterene Trifluoperazine
(±) Chlorpheniramine	Morphine-3-β-D-glucuronide	Trimethoprim
Chlorpromazine	Nalidixic acid	Trimipramine
Chlorquine	Naltrexone	Tryptamine
Cholesterol Clomipramine	Niacinamide Nifedipine	D, L-Tryptophan
Clonidine	Norethindrone	Tyramine
Cocacethylene	D-Norpropoxyphene	Tetrahydrocortisone, 3
Cocainemethochloride	Codeine	Acetate
(-) Cotinine	Noscapine	Tetrahydrocortisone 3 (β-D glucuronide)
Creatinine	D,L-Octopamine	Tetrahydrozoline
Deoxycorticosterone	Oxalic acid	Uric acid
Dextromethorphan	Oxazepam	Verapamil
Digoxin Diphenhydramine	Oxolinic acid	Zomepirac
Doxylamine	Oxycodone	
Egogone hydrochloride	Quinine	
Egogone methylester	Ranitidine	
(1R,2S)-(-)-Ephedrine	Salicylic acid	
L-Epinephrine	Secobarbital	
Ketoprofen	Serotonin (5-Hydroxytyramine)	
Labetalol	Sulfamethazine	
Levorphanol	Sulindac	
Loperamide	Temazepam	
Maprotiline	Tetracycline	
Meperidine		

REFERENCES

- Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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Substanz	Konzentration (ng/ml)
p-Hydroxymethamphetamine	30,000
D-Methamphetamine	1,000
L-Methamphetamine	8,000
(±)-3,4-Methylenedioxy-methamphetamine	2,000
Mephentermin	50,000

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LITERATUR

- Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73,1986

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