"DIAQUICK" BZO Dipstick (Benzodiazepines)

for Urine Samples

REF

Content

- Z02501CE** - 30 tests individually packed
 (30 x Ref. No: Z02501B).
 - 1 package insert
- Z02501B** - 1 test individually packed
 - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	Competitive immunochromatographic assay
Antigen	BZO conjugate,
Antibodies	Anti-BZO monoclonal antibody colloidal gold conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Sensitivity	300 ng/ml

INTENDED USE

The "DIAQUICK" BZO Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Oxazepam (major metabolite) in urine at a cut-off concentration of 300 ng/mL. This test will detect other Benzodiazepines, please refer to Analytical Specificity table in this package insert. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Benzodiazepines are the most widely used anxiolytic drugs. They are used extensively as anti-anxiety agents, hypnotics, muscle relaxants and anti-convulsants. They are taken orally or sometimes by injection. Benzodiazepines are metabolized in the liver, some of their metabolites also exhibit pharmacological activities. Both, Benzodiazepines and their metabolites, are excreted in the urine. The use of Benzodiazepines can result in drowsiness and confusion, moreover they potentiate alcohol and other CNS depressants. Psychological and physical dependence on benzodiazepines can be developed if high doses of the drug are given over a period of time. Risk of physical dependence increases if Benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception. Only trace amounts (less than 1%) of most Benzodiazepines are excreted unaltered in the urine; most of the concentration in urine is conjugated drug. The typical detection period for the Benzodiazepines in the urine is 3-7 days. The "DIAQUICK" BZO Dipstick (urine) is an easy to use urine test that can be performed without the use of an instrument. The test utilizes the antibody to selectively detect elevated levels of Benzodiazepines in urine. The "DIAQUICK" BZO Dipstick (urine) yields a positive result when the Benzodiazepines in urine exceeds the cut-off level.

TEST PRINCIPLE

The "DIAQUICK" BZO Dipstick (urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Benzodiazepines, if present in the urine specimen below the cut-off level, will not saturate the binding sites of the antibody in the test strip. The antibody coated particles will then be captured by immobilized Benzodiazepines-protein conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Benzodiazepines level exceeds the cut-off level, because it will saturate all the binding sites of anti-Benzodiazepines antibody. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- Do not use kit beyond the expiration date.
- To avoid cross contaminations use a new sample container for each sample
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded according to federal, state and local regulations.

STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability.

SAMPLE COLLECTION AND PREPARATION

Collect a urine sample in clean, dry container, either plastic or glass, without any preservatives. Urine specimens may be refrigerated (2-8°C) and stored up to 2 days. For longer storage, freeze samples (-20°C or below). Bring frozen or refrigerated samples to room temperature before testing. Urine samples exhibiting visible precipitates should be filtered, centrifuged or allowed to settle. Use only clear aliquots for testing.

ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

- Open the foil pouch by tearing along the notch and remove the test strip.
- Dip the stick in the urine specimen for 10 to 15 seconds submerging only up to the max mark.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear.
- The result should be read after 5 minutes.

Important: Do not interpret the result after 10 minutes. Waiting longer than 10 minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result at five minutes.


INTERPRETATION OF RESULTS

NEGATIVE:* Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Benzodiazepines concentration is below the detectable cut-off level.

***NOTE:** The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Benzodiazepines concentration exceeds the detectable cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

- This product is designed for use with human urine only.
- Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
- The test is a qualitative assay and is not for determining quantitative concentration levels or the level of intoxication.
- Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.
- Test does not distinguish between drugs of abuse and certain medications

PERFORMANCE CHARACTERISTICS

ACCURACY

A side-by-side comparison was conducted using the "DIAQUICK" BZO Dipstick (urine) and a leading commercially available BZO rapid test. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Other Rapid Test		Total Results		
	Positive	Negative			
"DIAQUICK" BZO Dipstick	Results Positive	131	4	135	
	Results Negative	14	149	163	
Total Results			145	153	298
% Agreement with this Rapid Test Kit			90%	97%	94%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

Method	GC/MS		Total Results		
	Positive	Negative			
"DIAQUICK" BZO Dipstick	Results Positive	131	8	139	
	Results Negative	4	157	161	
Total Results			135	165	300
% Agreement with GC/MS Analysis			97%	95%	96%

ANALYTICAL SENSITIVITY

A drug-free urine pool was spiked with Oxazepam at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates 100% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Oxazepam Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	27	3
300	Cut-off	30	11	19
375	+25%	30	5	25
450	+50%	30	0	30

ANALYTICAL SPECIFICITY

The following table lists compounds that are positively detected in urine by "DIAQUICK" BZO Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
Alprazolam	196
a-hydroxyalprazolam	1,262
Bromazepam	1,562
Chlordiazepoxide	1,562
Clobazam	98
Clonazepam	781
Clorazepate dipotassium	195
Delorazepam	1,562
Desalkylflurazepam	390
Diazepam	195
Estazolam	2,500
Flunitrazepam	390
(±) Lorazepam	1,562
RS-Lorazepam glucuronide	156
Midazolam	12,500
Nitrazepam	98
Norchlordiazepoxide	195
Nordiazepam	390
Oxazepam	300
Temazepam	98
Triazolam	2,500

PRECISION

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Oxazepam, 25% Oxazepam above and below the cut-off and 50% Oxazepam above and below the 300ng/mL cut-off was provided to each site. The following results were tabulated:

Oxazepam conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
150	15	15	0	13	2	13	2
225	15	6	9	7	8	13	2
375	15	0	15	1	14	3	12
450	15	0	15	0	15	0	15

EFFECT OF URINARY SPECIFIC GRAVITY

Fifteen (15) urine specimens with specific gravity range from 1.001 to 1.032 were spiked with 150 ng/mL and 450 ng/mL of Oxazepam respectively. The "DIAQUICK" BZO Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity does not affect the test results.

EFFECT OF THE URINARY PH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Oxazepam to 150 ng/mL and 450 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" BZO Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH does not interfere with the performance of the test.

CROSS-REACTIVITY

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Oxazepam positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" BZO Dipstick (urine) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

Acetaminophen	Estro-3-sulfate	Oxymetazoline
Acetophenetidin	Ethyl-p-aminobenzoate	Papaverine
N-Acetylprocainamide	Fenoprofen	Penicillin-G
Acetylsalicylic acid	Furosemide	Pentazocine
Aminopyrine	Gentisic acid	Pentobarbital
Amitypyline	Hemoglobin	Perphenazine
Amobarbital	Hydralazine	Phencyclidine
Ampicillin	Hydrochlorothiazide	Phenazone
L-Ascorbic acid	Hydrocodone	Phenobarbital
DL-Amphetamine sulfate	Hydrocortisone	Phentermine
Apomorphine	O-Hydroxyhippuric acid	Trans-2-phenylcydo-propylamine hydrochloride
Aspartame	p-Hydroxyamphetamine	L-Phenylephrine
Atropine	p-Hydroxy-methamphetamine	β-Phenylethylamine
Benzic acid	methamphetamine	Phenylpropanolamine
Benzoic acid	3-Hydroxytyramine	Prednisolone
Benzoylcegonine	Ibuprofen	Prednisone
Benzphetamine	Imipramine	Procaine
Bilirubin	lproniazid	Promazine
(±) - Brompheniramine	(±) - Isoproterenol	Promethazine
Caffeine	Isoxsuprine	DL-Propranolol
Cannabidiol	Ketamine	D-Propoxyphene
Cannabinol	Ketoprofen	D-Pseudoephedrine
Chloralhydrate	Labelalol	Sulindac
Chloramphenicol	Loperamide	Tetracycline
Chlorothiazide	Maprotiline	Tetrahydrocortisone, 3-acetate
(±) - Chlorpheniramine	MDE	
Chlorpromazine	Meperidine	
Chlorquine	Meprobamate	
Cholesterol	Methadone	
Ciomiapramine	(L) Methamphetamine	
Clonidine	Methoxyphenamine	
Cocaeethylene	(±) - 3,4-Methylenedioxy-amphetamine hydrochloride	
Cocaine hydrochloride	(±) - 3,4-Methylenedioxy-meth-amphetamine hydrochloride	
Codeine	Morphine-3-β-D glucuronide	Tetrahydrocortisone 3-(β-D-glucuronide)
Cortisone	Morphine Sulfate	Tetrahydrozoline
(-) Cotinine	Nalidixic acid	Thiamine
Creatinine	Naloxone	Thioridazine
Deoxycorticosterone	Naltrexone	DL-Tyrosine
Dextromethorphan	Naproxen	Tolbutamide
Diclofenac	Niacinamide	Triamterene
Diflunisal	Nifedipine	Trifluoperazine
Digoxin	Norcodeine	Trimethoprim
Diphenhydramine	Norethindrone	Trimipramine
Doxylamine	D-Norpropoxyphene	Tryptamine
Ecgoinone hydrochloride	Noscapine	DL-Tryptophan
Ecgoinone methylester	DL-Octopamine	Tyramine
(-) - α - Ephedrine	Oxalic acid	Uric acid
[1R,2S] (-) Ephedrine	Oxolinic acid	Verapamil
(L) - Epinephrine	Oxycodone	Zomepirac
Erythromycin		
β-Estradiol		

REFERENCES

- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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44-Acetylaminophenol	5,5-Diphenylhydantoin	Lithium carbonate	Trans-2- phenyl cyclopropylamine
Acetone	Disopyramide	Loperamide	L-Phenylephrine
Acetophenetidin	Doxylamine	Meprobamate	B-Phenylethylamine
Acetylsalicylic acid	Ecgonine hydrochloride	Mephentermine	Phenylpropanolamine
N-Acetylprocainamide	Ecgonine methylester	Mephentermine	(D,L-norephedrine)
Albumin	EDDP	Methadone	(±) Phenylpropanolamine
Aminopyrine	Efavirenz (Sustiva)	D-Methamphetamine	Prednisolone
Amiripryline	EMDP	L-Methamphetamine	Prednisone
Amobarbital	Ephedrine	Methaqualone	5 beta- pregnane3alpha17alpha- 21triol 21
Amoxapine	(1r,2s)-(-)-Ephedrine	Methoxyphenamine	Procaine
Amoxicillin	(-)-ψ-Ephedrine	(-) 3,4-Methylenedioxy- amphetamine (MDA)	Promazine
L-Amphetamine	(±)-Epinephrine	Erythromycin	Promethazine
Ampicillin	Erythromycin	β-Estradiol	D,L-Propranolol
Apomorphine	β-Estradiol	(+) 3,4 Methylendioxy- methamphetamine	D-Propoxyphene
Aspartame	Estrone-3-sulfate	Ethanol (Ethyl alcohol)	Methaqualone
Atropine	Ethanol (Ethyl alcohol)	Ethyl-p-aminobenzoate	Metoprolol
Benzilic acid	Etodolac	Famprofazone	Morphine sulfate
Benzolic acid	Famprofazone	Fenfluramine	Morphine-
Benzoylcegonine	Fenfluramine	Fenpropfen	3-β-D-glucuronide
Benzphetamine	Bilirubin	Fentanyl	Nalidixic acid
Bilirubin	(±)-Brompheniramine	Buspirone	Naloxone
Buspirone	Buspirone	Caffeine	Naltrexone
Caffeine	Benzilic acid	Cannabidiol	D (+) Glucose
Benzilic acid	Benzilic acid	Cannabiol	Guaiacol Glyceryl Ether
Etodolac	Benzilic acid	Chloralhydrate	Guaiacol Glyceryl Ether
Famprofazone	Benzilic acid	Chloramphenicol	carbamate
Fenfluramine	Benzilic acid	Chlordiazepoxide	Chloroquine
Fenpropfen	Benzilic acid	Chloroquine	Hemoglobin
3-β-D-glucuronide	Benzilic acid	Chlorothiazide	Hydralazine
Ranitidine	Benzilic acid	(+) Chlorpheniramine	Hydrochlorothiazide
Temazepam	Benzilic acid	(±)-Chlorpheniramine	Hydrocodone
Riboflavin	Benzilic acid	Chlorpromazine	Hydrocortisone
Salicylic acid	Benzilic acid	Chlorprothixene	Hydrocortisone
Secobarbital	Benzilic acid	Cholesterol	Hydromorphone
Serotonin	Benzilic acid	Cimetidine	p-Hydroxyamphetamine
(5-hydroxytryramine)	Benzilic acid	Clomipramine	O-Hydroxyhippuric acid
Sodium chloride	Benzilic acid	Clonidine	Noscapine
Sulfamethazine	Benzilic acid	Cocaine	p-Hydroxymethamphetamine
Sulindac	Benzilic acid	Cocaine HCl	p-Hydroxymorephedrine
Ternazepam	Benzilic acid	Codeine	Hydroxyzine
Tetrazepam	Benzilic acid	Cortisone	Oxalic acid
Tetracycline	Benzilic acid	(-) Cotinine	Oxazepam
Tetrahydrocortisone,	Benzilic acid	Creatinine	Oxolinic acid
3-acetate	Benzilic acid	Cyclobarbitol	Oxycodone
Tetrahydrozoline	Benzilic acid	Cyclobenzaprine	Oxymetazoline
Thebaine	Benzilic acid	Deoxycorticosterone	Oxymorphone
Theophylline	Benzilic acid	(-) Deoxyephedrine	Papaverine
Thiamine	Benzilic acid	R (-) Deprenyl HCl	Penicilline
Thioridazine	Benzilic acid	Dextromethorphan	Penicillin-G
(chlorpromazine)	Benzilic acid	Diazepam	Pentazocine
L-Thyroxine	Benzilic acid	Diclofenac	Perphenazine
Tolbutamide	Benzilic acid	Dicyclomine	Phenacetyl
Cis-Tramadol	Benzilic acid	Diffunisal	L-Ephedrine
Trazodone	Benzilic acid	Digoxin	L-Epinephrine
Triamterene	Benzilic acid	4-Dimethylamino-antipyrine	Levopropriolol
Trifluoperazine	Benzilic acid	Diphenhydramine	Lidocaine
Trimethoprim	Benzilic acid		Lindane
Trimipramine	Benzilic acid		(hexachlorocyclohexane)
Tryptamine	Benzilic acid		Phentemine
D, L-Tryptophan	Benzilic acid		
Tyramine	Benzilic acid		
D, L-Tyrosine	Benzilic acid		
Uric acid	Benzilic acid		
Verapamil	Benzilic acid		
Zomepirac	Benzilic acid		

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