# "DIAQUICK" THC Dipstick (Tetrahydrocannabinol) for urine samples

## REF

### Content

<b>Z02502CE</b>	- 30 tests individually packed (30 x Ref. No: Z02502B). - 1 package insert
<b>Z02502B</b>	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.  
 For diagnosis and therapeutic monitoring only.  
 For use by medical professionals only.**

## GENERAL INFORMATION

<b>Method</b>	competitive immunochromatographic assay
<b>Antigen / Antibodies</b>	mouse monoclonal anti-Marihuana antibody-coupled particles Marihuana-protein conjugate
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	human urine
<b>Results</b>	within 5 minutes at room temperature
<b>Sensitivity</b>	50ng/ml

## INTENDED USE

The "DIAQUICK" THC Dipstick (urine) is a rapid chromatographic immunoassay for the detection of 11-nor-D9 -THC-9 COOH (THC metabolite) in human urine at a cut-off concentration of 50 ng/mL. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrophotometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

## CLINICAL SIGNIFICANCE

THC (D9-tetrahydrocannabinol) is the primary active ingredient in cannabinoids (Marihuana). When smoked or orally administered, it produces euphoric effects. Users have impaired short term memory and slowed learning. They may also experience transient episodes of confusion and anxiety. Long term relatively heavy use may be associated with behavioral disorders. The peak effect of smoking Marihuana occurs in 20-30 minutes and the duration is 90-120 minutes after one cigarette. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 3-10 days after smoking. The main metabolite excreted in the urine is 11-nor-D9-tetrahydrocannabinol-9-carboxylic acid (D9-THC-COOH). The "DIAQUICK" THC Dipstick (urine) is an easy to use urine test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Marihuana in urine. The "DIAQUICK" THC Dipstick (urine) yields a positive result when the concentration of Marihuana in urine exceeds 50 ng/mL. This is the suggested cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration.

## TEST PRINCIPLE

The "DIAQUICK" THC Dipstick (urine) is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Marihuana, if present in the urine specimen below 50 ng/mL, will not saturate the binding sites of the antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized Marihuana conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Marihuana level is above 50 ng/mL because it will saturate all the binding sites of anti-Marihuana antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

## WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after the expiration date.
- The dipstick should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used strip should be discarded in a proper biohazard container after testing.

## STORAGE

The kit can be stored at room temperature or refrigerated (2-30°C). The test strip is stable through the expiration date printed on the sealed pouch. The test strip must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

## SAMPLE COLLECTION AND PREPARATION

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

## SPECIMEN STORAGE

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

## ASSAY PROCEDURE

**Allow the test strip, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.**

- Remove the test strip from the sealed pouch and use it as soon as possible.
- With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.


## INTERPRETATION OF RESULTS

**NEGATIVE:**\* Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Marihuana concentration is below the detectable level of 50 ng/mL.

**\*NOTE:** The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

**POSITIVE:** One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Marihuana concentration is above the detectable level of 50 ng/mL.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

## QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability.

The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

## LIMITATIONS

1. This product is designed for use with urine only.
2. Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
3. The test is a qualitative urine assay and is not for determining quantitative concentration levels or the level of intoxication.
4. Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when the drug is present but below the cut-off level of the test.
6. The Test does not distinguish between drugs of abuse and certain medications.

## PERFORMANCE CHARACTERISTICS

### Accuracy

A three way side-by-side comparison was conducted using the "DIAQUICK" THC Dipstick (urine) and a leading commercially available THC rapid test. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Results	Other THC Rapid Test		Total Results
		Positive	Negative	
DIAQUICK THC Dipstick	Positive	140	0	140
	Negative	3	157	160
Total Results		143	157	300
% Agreement with this Rapid Test Kit		98%	100%	99%

When compared to GC/MS at 50 ng/mL, the following results were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
DIAQUICK THC Dipstick	Positive	118	22	140
	Negative	4	156	160
Total Results		122	178	300
% Agreement with GC/MS at 50 ng/mL		97%	88%	92%

When compared to GC/MS at 25 ng/mL, the following results were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
DIAQUICK THC Dipstick	Positive	135	5	140
	Negative	6	154	160
Total Results		141	159	300
% Agreement with GC/MS at 25 ng/mL		96%	97%	96%

### Analytical Sensitivity

A drug-free urine pool was spiked with 11-nor- $\Delta^9$ -Tetrahydrocannabinol-9-carboxylic acid at the following concentrations: 75 ng/mL, 62.5 ng/mL, 37.5 ng/mL, 25 ng/mL, and 0 ng/mL. The result demonstrates 100% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

THC Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
25	50%	30	30	0
37.5	75%	30	12	18
50	Cut-off	30	1	29
62.5	125%	30	1	29
75	150%	30	0	30

### Analytical Specificity

The following table lists compounds and their respective concentrations in urine that yield a positive result in the "DIAQUICK" THC Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
Cannabinol	20,000
11-nor- $\Delta^9$ -THC-9 COOH	30
11-nor- $\Delta^9$ -THC-9 COOH	50
$\Delta^8$ -THC	15,000
$\Delta^9$ -THC	15,000

### Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no THC, 25% THC above and below the cut-off, and 50% THC above and below the 50 ng/mL cut-off was provided to each site. For the specimens below the 25 % cut-off concentration, the 3 sites demonstrated 98% agreement with each other. For the -25% to +25% cut-off specimens, the 3 sites demonstrated 83% agreement with each other. For specimens above the 25% cut-off concentration, the 3 sites demonstrated 94% agreement with each other. For all results, the 3 sites were found to have a 91% agreement with each other.

### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Marijuana positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" THC Dipstick (urine) at a concentration of 100  $\mu$ g/mL.

### Non Cross-Reacting Compounds

4-Acetamidophenol	Fenoprofen	Pentazocine
Acetophenetidin	Furosemide	Pentobarbital
N-Acetylprocainamide	Genistic acid	Perphenazine
Acetylsalicylic acid	Hemoglobin	Phencyclidine
Aminopyrine	Hydralazine	Phenelzine
Amitypyline	Hydrochlorothiazide	Phenobarbital
Amobarbital	Hydrocodone	Phentermine
Amoxicillin	Hydrocortisone	L-Phenylephrine
Ampicillin	O-Hydroxyhippuric acid	$\beta$ -Phenylethylamine
L-Ascorbic acid	3-Hydroxytyramine	$\beta$ -Phenylthylamine
D,L-Amphetamine	Ibuprofen	Phenylpropanolamine
L-Amphetamine	Imipramine	Prednisolone
Apomorphine	Iproniazid	Prednisone
Aspartame	(±) - Isoproterenol	Procaine
Atropine	Isoxsuprine	Promazine
Benzilic acid	Ketamine	Propoxyphene
Benzoic acid	Ketoprofen	D,L-Propranolol
Benzoylcegonine	Labetalol	D-Propoxyphene
Benzphetamine	Levorphanol	D-Pseudoephedrine
Bilirubin	Loperamide	Quinidine
(±) - Brompheniramine	Maprotiline	Quinine
Caffeine	Meprobamate	Ranitidine
Cannabidiol	Methadone	Salicylic acid
Chloralhydrate	Methoxyphenamine	Secobarbital
Chloramphenicol	(+) 3,4-Methylenedioxy-amphetamine	Serotonin
Chlordiazepoxide		
Chlorothiazide	(+) 3,4-Methylenedioxy-methamphetamine	Sulfamethazine
(±) Chlorpheniramine		Sulindac
Chlorpromazine	Methyphenidate	Temazepam
Chlorquine	Methpyrion	Tetracycline
Cholesterol	Morphine-3- $\beta$ -D-glucuronide	Tetrahydrocortisone, Acetate
Clomipramine		Tetrahydrocortisone 3 ( $\beta$ -D glucuronide)
Clonidine	Nalidixic acid	Tetrahydrozoline
Cocaine hydrochloride	Nalorphine	Thebaine
Codeine	Naloxone	Thiamine
Cortisone	Naltrexone	Thioridazine
(-) Cotinine	Naproxen	D, L-Tyroxine
Creatinine	Niacinamide	Tolbutamide
Deoxycorticosterone	Nifedipine	Triamterene
Dextromethorphan	Norcodein	Trifluoperazine
Diazepam	Norethindrone	Trimethoprim
Diclofenac	D-Norpropoxyphene	Trimipramine
Diffunisal	Noscaphine	Tryptamine
Digoxin	D,L-Octopamine	D, L-Tryptophan
Diphenhydramine	Oxalic acid	Tyramine
Doxylamine	Oxazepam	PRD, L-Tyrosine
Ecgonine hydrochloride	Oxolinic acid	Uric acid
Ecgonine methylester	Oxycodone	Verapamil
(-) Y Ephedrine	Oxymetazoline	Zomepirac
Erythromycin	p-Hydroxy-methamphetamine	
$\beta$ -Estradiol	Papaverine	
Estrone-3-sulfate	Penicillin-G	
Ethyl-p-aminobenzoate		

### REFERENCES

- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

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