# "DIAQUICK" AMP Dipstick (Amphetamine) for Urine Samples

REF	Content
Z02504CE	- 30 tests individually packed (30 x Ref. No: Z02504B). - 1 package insert
Z02504B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.  
 For diagnosis and therapeutic monitoring only.  
 For use by medical professionals only.**

## GENERAL INFORMATION

<b>Method</b>	competitive immunochromatographic assay
<b>Antigen</b>	AMP conjugate,
<b>Antibodies</b>	Anti-AMP monoclonal antibody colloidal gold conjugate
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	human urine
<b>Results</b>	within 5 minutes at room temperature
<b>Sensitivity</b>	1,000 ng/ml

## INTENDED USE

The "DIAQUICK" AMP Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Amphetamines in human urine. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

## CLINICAL SIGNIFICANCE

Amphetamines (Amphetamine, Methamphetamine, and the structurally related „designer“ drugs are sympathomimetic amines whose biological effects include potent central nervous system (CNS) stimulation, anorectic, hyperthermic and cardiovascular properties. They are usually taken orally, intravenously or by smoking. Amphetamines are readily absorbed from the gastrointestinal tract and are then either deactivated by the liver or excreted unchanged in the urine. Amphetamine is metabolized to deaminated (hippuric and benzoic acids) and hydroxylated metabolites. Methamphetamine is partially metabolized to amphetamine, its major active metabolite. Amphetamines increase the heart rate and blood pressure and suppress the appetite. Some studies indicate that heavy abuse may result in permanent damage to certain essential nerve structures in the brain.

The "DIAQUICK" AMP Dipstick (urine) is an easy to use urine test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Amphetamine in urine. The "DIAQUICK" AMP Dipstick (urine) yields a positive result when Amphetamines in urine exceed 1,000 ng/mL.

## TEST PRINCIPLE

The "DIAQUICK" AMP Dipstick (urine) is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Amphetamine, if present in the urine specimen below 1,000 ng/mL, will not saturate the binding sites of the antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized Amphetamine conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Amphetamine level exceeds 1,000 ng/mL because it will saturate all the binding sites of anti-Amphetamine antibodies. A drug-positive urine specimen will not generate a colored line in the test line region, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

## WARNINGS AND PRECAUTIONS

- Do not use kit beyond the expiration date.
- To avoid cross contaminations use a new sample container for each sample
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded according to federal, state and local regulations.

## STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability.

## SAMPLE COLLECTION AND PREPARATION

Collect a urine sample in clean, dry container, either plastic or glass,

without any preservatives. Urine specimens may be refrigerated (2-8°C) and stored up to 2 days. For longer storage, freeze samples (-20°C or below).

Bring frozen or refrigerated samples to room temperature before testing. Urine samples exhibiting visible precipitates should be filtered, centrifuged or allowed to settle. Use only clear aliquots for testing.

## ASSAY PROCEDURE

**Bring an unopened foil pouch and the sample specimen to room temperature.**

- Open the foil pouch by tearing along the notch and remove the test strip.
- Dip the stick in the urine specimen for 10 to 15 seconds submerging only up to the max mark.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear.
- The result should be read after 5 minutes.

**Important:** Do not interpret the result after 10 minutes. Waiting longer than 10 minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result at five minutes.


## INTERPRETATION OF RESULTS

**NEGATIVE:** Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Amphetamine concentration is below the detectable level (1,000 ng/mL).

**NOTE:** The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

**POSITIVE:** One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Amphetamine concentration exceeds the detectable level (1,000 ng/mL).

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test strip. If the problem persists, discontinue using this batch immediately and contact your local distributor.

## QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability.

The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

## LIMITATIONS

- This product is designed for use with human urine only.
- Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
- The test is a qualitative urine assay and is not intended for determining quantitative concentration levels or the level of intoxication.
- Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.
- Test does not distinguish between drugs of abuse and certain medications

## PERFORMANCE CHARACTERISTICS

### Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" AMP Dipstick (urine) and a leading commercially available AMP rapid test. Testing was performed on 300 clinical specimens. Ten percent of the specimens employed were either at -25% or +25% level of the cut-off concentration of 1,000 ng/mL Amphetamine. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Other AMP Rapid Test		Total Results
	Positive	Negative	
"DIAQUICK" AMP Dipstick	141	0	141
	5	154	159
<b>Total Results</b>		146	300
<b>% Agreement with this Rapid Test Kit</b>		97%	100%

When compared at 1,000 ng/mL cut-off with GC/MS, the following results were tabulated:

Method	GC/MS		Total Results
	Positive	Negative	
"DIAQUICK" AMP Dipstick	Results Positive	132	141
	Negative	4	159
<b>Total Results</b>	<b>136</b>	<b>164</b>	<b>300</b>
<b>% Agreement with GC/MS Analysis</b>	<b>97%</b>	<b>95%</b>	<b>96%</b>

#### Analytical Sensitivity

A drug-free urine pool was spiked with Amphetamine at the following concentrations: 0 ng/mL, 500 ng/mL, 750 ng/mL, 1,000 ng/mL, 1,250 ng/mL and 1,500 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

AMP Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
500	-50%	30	30	0
750	-25%	30	22	8
1,000	Cut-off	30	12	18
1,250	+25%	30	2	28
1,500	+50%	30	0	30

#### Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" AMP Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
D-Amphetamine	1,000
D,L-Amphetamine sulfate	3,000
L-Amphetamine	50,000
(±) 3,4-MethylenedioxyAmphetamine	2,000
Phentermine	3,000

#### Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Amphetamine, 25% Amphetamine above and below the cut-off, and 50% Amphetamine above and below the 1,000 ng/mL cut-off was provided to each site. The results are given below:

Amphetamine concentration (ng/mL)	n	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
500	15	15	0	15	0	14	1
750	15	13	2	11	4	11	4
1250	15	6	9	4	11	4	11
1500	15	2	13	1	14	1	14

#### Effect of Urinary Specific Gravity

Fifteen (15) urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 500 ng/mL and 1,500 ng/mL of Amphetamine respectively. The "DIAQUICK" AMP Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

#### Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Amphetamine to 500 ng/mL and 1,500 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" AMP Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

#### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Amphetamine positive urine. The following compounds were found not to cross-react when tested at a 100 µg/ml concentration in urine.

4-Acetamidophenol	Estrone-3-sulfate	Oxolinic acid
Acetophenetidin	Ethyl-p-aminobenzoate	Oxycodone
N-Acetylprocainamide	Fenfluramine	Oxymetazoline
Acetylsalicylic acid	Fenpropfen	Papaverine
Aminopyrine	Furosemide	Penicillin-G
Amitypyline	Genitistic acid	Pentazocine
Amobarbital	Hemoglobin	Pentobarbital
Amoxicillin	Hydralazine	Perphenazine
Ampicillin	Hydrochlorothiazide	Phencyclidine
L-Ascorbic acid	Hydrocodone	Phenelzine
Apomorphine	Hydrocortisone	Phenobarbital
Aspartame	p-HydroxyAmphetamine	L-Phenylephrine
Atropine	O-Hydroxyhippuric acid	β-Phenylethylamine
Benzilic acid	p-Hydroxy-methamphetamine	Phenylpropanolamine
Benzoic acid	3-Hydroxytyramine	Prednisolone
Benzoylcegonine	Ibuprofen	Prednisone
Benzphetamine	Imipramine	Procaine
Billirubin	(±) - Brompheniramine	Promazine
(±) - Brompheniramine	(±) -Isoproterenol	Promethazine
Caffeine	Isoxsuprine	D,L-Propranolol
Cannabidiol	Ketamine	D-Propoxyphene
Cannabinol	Ketoprofen	D-Pseudoephedrine
Chloralhydrate	Labetalol	Quinidine
Chloramphenicol	Levorphanol	Quinine
Chlordiazepoxide	Loperamide	Ranitidine
Chlorothiazide	Maprotiline	Salicylic acid
(±) Chlorpheniramine	Meperidine	Secobarbital
Chlorpromazine	Meprobamate	Serotonin
Cholesterol	D-methAmphetamine	Sulfamethazine
Clomipramine	(L)-methAmphetamine	Sulindac
Clonidine	Methoxyphenamine	Temazepam
Cocaine hydrochloride	3,4-Methylenedioxyethyl-Amphetamine	Tetracycline
Codeine	(+) 3,4-Methylenedioxy-methAmphetamine	Tetrahydrocortisone, 3 Acetate
Cortisone	(-) Cotinine	Tetrahydrocortisone
Deoxycorticosterone	Methylphenidate	Tetrahydrozoline
Dextromethorphan	Morphine-3-β-D-glucuronide	Thebaine
Diazepam	Nalidixic acid	Thiamine
Diclofenac	Naloxone	Thioridazine
Diffunisal	Naltrexone	Tolbutamine
Digoxin	Naproxen	Triamterene
Diphenhydramine	Niacinamide	Trifluoperazine
Doxylamine	Nifedipine	Trimethoprim
Ecgonine hydrochloride	Norcodein	Trimipramine
Ecgonine methyl ester	Norethindrone	D, L-Tryptophan
(IR,2S)-(-)-Ephedrine	D-Norpropoxyphene	Tyramine
L-Ephedrine	Noscopamine	D, L-Tyrosine
(-) Y Ephedrine	D,L-Octopamine	Uric acid
Erythromycin	Oxalic acid	Verapamil
β-Estradiol	Oxazepam	Methadone

#### REFERENCES

- Baselt, R.C., Disposition of Toxic Drugs and Chemicals in Man, 2nd Ed., Biomedical Publ., Davis, CA, p 488 (1982).
- Cody, J.T., and Schwarzhoff, R., J. Anal. Toxicol., 17:2630 (1993).
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
- Federal Register, Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Test Programs, 69, 11970-11979, 1988, Federal Register 59, 110, 29908-29931, 1994.
- Dasgupta, A., Saldana, S., Kinnaman, G., Smith, M., and Johansen, K., Clin. Chem., 39 (1):104-108 (1993).
- FDA Guidance for Labeling Urine Drugs of Abuse Screening Testing, Kshitij Mohan, 7/21/1987.
- Froelich, P.E., Gross, G.: Separation and Detection of Phencyclidine in Urine by Gas Chromatography. J. Chromatograph 1977;137:135-143.
- Gupta R.C. et al: Determination of Phencyclidine in Urine and Illicit Street Drug Samples. Clin. Toxicol 1975; 8:611-621.
- Oellerich, M.: Enzyme Immunoassays In Clinical Chemistry: Present Status and Trends. J. Clin. Chem. Biochem., 1980;18:197-208.

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