"DIAQUICK" MOP Dipstick (Morphine) for urine samples

REF	Content
Z02505CE	- 30 tests individually packed (30 x Ref. No: Z02505B). - 1 package insert
Z02505B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	competitive immunochromatographic assay
Antigen / Antibodies	monoclonal anti-morphine-antibody-coupled particles and protein conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Cut-off	300 ng/ml

INTENDED USE

The "DIAQUICK" MOP Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of morphine in human urine at a cut-off concentration of 300 ng/ml. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Opioid analgesics comprise a large group of substances which control pain by depressing the central nervous system. Morphine is the prototype compound of this group. Morphine is excreted unmetabolized, and is also the major metabolic product of codeine and heroin. Morphine is detectable in the urine for several days after an Opiate dose. The test utilizes a monoclonal antibody to selectively detect elevated levels of Morphine in urine. The "DIAQUICK" MOP Dipstick (urine) yields a positive result the Morphine in urine reaches 300 ng/mL. This is the suggested cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration.

TEST PRINCIPLE

The "DIAQUICK" MOP Dipstick (urine) is an immunoassay based on the principle of competitive binding. Drugs, which may be present in the urine specimen, compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Morphine, if present in the urine specimen below 300 ng/mL, will not saturate the binding sites of the antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized Morphine conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Morphine level is at or above 300 ng/mL because it will saturate all the binding sites of anti-Morphine antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS REQUI PROVIDED

- test strips
- package insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded in a proper biohazard container after testing.

STORAGE

The kit can be stored at room temperature or refrigerated (2-30°C). The test strip is stable through the expiration date printed on the sealed pouch. The test strip must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE

Allow the test strip, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Remove the test strip from the sealed pouch and use it as soon as possible.
- With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. The result should be read at 5 minutes.
 Do not interpret the result after 10 minutes.


INTERPRETATION OF RESULTS

NEGATIVE:* Two lines appear. One line should be in the control region (C), and another line should be in the test region (T). This negative result indicates that the morphine concentration is below the detectable level of 300 ng/ml.

***NOTE:** The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Methadone concentration is above the detectable level of 300 ng/mL.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

1. This product is designed for use with urine only.
2. Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
3. The test is a qualitative urine assay and is not for determining quantitative concentration levels or the level of intoxication.
4. Adulterants such as bleaching or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" MOP Dipstick (urine) and a leading commercially available MOP rapid test. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method		Other MOP Rapid Test		Total Results
DIAQUICK MOP Dipstick	Results	Positive	Negative	
	Positive	150	0	150
	Negative	0	150	150
Total Results		150	150	300
% Agreement with this Rapid Test Kit		100%	100%	100%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

Method		GC/MS		Total Results
DIAQUICK MOP Dipstick	Results	Positive	Negative	
	Positive	141	9	150
	Negative	0	150	150
Total Results		141	159	300
% Agreement with GC/MS Analysis		100%	94%	97%

Analytical Sensitivity

A drug-free urine pool was spiked with morphine at the following concentrations: 0 ng/ml, 150 ng/mL, 225 ng/ml, 300 ng/ml, 375 ng/ml and 450 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

MOP Concentration (ng/mL)	Percent of cut-off	n	Visual Result	
			Negative	Positive
0	0%	30	30	0
150	-50%	30	30	0
225	-25%	30	25	5
300	Cut-off	30	17	13
375	+25%	30	1	29
450	+50%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" MOP Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
Codeine	300
Ethylmorphine	6,250
Hydrocodone	50,000
Hydromorphone	3,125
Levophanol	1,500
6-Monoacetylmorphine	400
Morphine	300
Morphine 3-β-D-glucuronid	1,000
Norcodeine	6,250
Normorphine	100,000
Oxycodone	30,000
Oxymorphone	100,000
Procaine	15,000
Thebaine	6,250

Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing no morphine, 25% morphine above and below the cut-off and 50% morphine above and below the 300 ng/mL cut-off was provided to each site. The following results were tabulated:

Morphine conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
150	15	13	2	13	2	15	0
225	15	3	12	7	8	10	5
375	15	0	15	1	14	0	15
450	15	0	15	0	15	0	15

Effect of Urinary Specific Gravity

Fifteen urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 150 ng/ml and 450 ng/ml of Morphine respectively. The "DIAQUICK" MOP Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Morphine to 150 ng/mL and 450 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" MOP Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Morphine-positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" MOP Dipstick (urine) at a concentration of 100 µg/ml.

Non Cross-Reacting Compounds

Acetaminophen	Estrone-3-sulfate	Oxymetazoline
Acetophenetidin	Ethyl-p-aminobenzoate	Papaverine
N-Acetylprocainamide	Fenpropfen	Penicillin-G
Acetylsalicylic acid	Furosemide	Pentazocine
Aminopyrine	Genistic acid	Pentobarbital
Amitypyline	Hemoglobin	Perphenazine
Amobarbital	Hydralazine	Phenylephrine
Amoxicillin	Hydrochlorothiazide	Phenelzine
Ampicillin	Hydrocodone	Phenobarbital
L-Ascorbic acid	Hydrocortisone	Phentermine
DL-Amphetamine sulfate	O-Hydroxyhippuric acid	Trans-2-phenylcyclopropylamine hydrochloride
Apomorphine	p-Hydroxyamphetamine	L-Phenylephrine
Aspartame	p-Hydroxy-methamphetamine	β-benyltetraamine
Atropine	3-Hydroxytyramine	Phenyl-propanolamine
Benzilic acid	Ibuprofen	Prednisolone
Benzoic acid	Imipramine	Prednisone
Benzoylcegonine	Iproniazid	Procaine
Benzphetamine	(±) - Isoproterenol	Promazine
Bilirubin	(±) - Brompheniramine	Promethazine
(±) - Brompheniramine	Caffeine	DL-Propanolol
Caffeine	Cannabidiol	D-Propoxyphene
Cannabidiol	Cannabitol	D-Pseudoephedrine
Cannabitol	Chloralhydrate	Quinacrine
Chloralhydrate	Chloramphenicol	Quinidine
Chloramphenicol	Chlorothiazide	Quinine
Chlorothiazide	(±) - Chlorpheniramine	Ranitidine
(±) - Chlorpheniramine	Chlorpromazine	Salicylic acid
Chlorpromazine	Chloroquine	Secobarbital
Chloroquine	Cholesterol	Serotonin
Cholesterol	Clomipramine	Sulfamethazine
Clomipramine	Clonidine	Sulindac
Clonidine	Cocaine hydrochloride	Temazepam
Cocaine hydrochloride	Codeine hydrochloride	Tetracycline
Codeine hydrochloride	Cotinine (-)	Tetrahydro-cortisone, 3-acetate
(-) Cotinine	Creatinine	Tetrahydro-cortisone 3-(β-D-glucuronide)
Creatinine	Deoxycorticosterone	Tetrahydrozoline
Deoxycorticosterone	Dextromethorphan	Thebaine
Dextromethorphan	Diazepam	Thiamine
Diazepam	Diclofenac	Thioridazine
Diclofenac	Diffunisal	DL-Tyrosine
Diffunisal	Digoxin	Tolbutamide
Digoxin	Diphenhydramine	Triamterene
Diphenhydramine	EDDP	Trifluoperazine
EDDP	EMDP	Trimethoprim
EMDP	Egonine hydrochloride	Trimipramine
Egonine hydrochloride	Egonine methylester	Tryptamine
Egonine methylester	(-) - ψ - Ephedrine	DL-Tryptophan
(-) - ψ - Ephedrine	[1R,2S] (-) Ephedrine	Tyramine
[1R,2S] (-) Ephedrine	(L) - Epinephrine	Uric acid
(L) - Epinephrine	Erythromycin	Verapamil
Erythromycin	β-Estradiol	Zomepirac
β-Estradiol		

REFERENCES

1. Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216
2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
3. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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Acetaminophen	Estrone-3-sulfate	Oxymetazoline
Acetophenetidin	Ethyl-p-aminobenzoate	Papaverine
N-Acetylprocainamide	Fenoprofen	Penicillin-G
Acetylsalicylic acid	Furosemide	Pentazocine
Aminopyrine	Genisic acid	Pericardial
Amitypyline	Hemoglobin	Perphenazine
Amobarbital	Hydralazine	Phencyclidine
Amoxicillin	Hydrochlorothiazide	Phenelzine
Ampicillin	Hydrocodone	Phenobarbital
L-Ascorbic acid	Hydrocortisone	Phentermine
DL-Amphetamine sulfate	O-Hydroxyhippuric acid	Trans-2-phenylcyclopropylamine hydrochloride
Aprorphine	p-Hydroxyamphetamine	L-Phenylephrine
Asparfame	p-Hydroxy-methamphetamine	β-henylethylamine
Atropine	3-Hydroxytyramine	Phenyl-propanolamine
Benzilic acid	Ibuprofen	Prednisolone
Benzoic acid	Imipramine	Prednisone
Benzoylecgonine	Iproniazid	Procaine
Benzphetamine	(±) - Isoproterenol	Promazine
Bilirubin	(±) - Brompheniramine	Promethazine
(±) - Brompheniramine	Caffeine	DL-Propriolol
Cannabidiol	Cannabidiol	D-Propoxyphene
Cannabiol	Chloralhydrate	D-Pseudoephedrine
Chloralhydrate	Chloramphenicol	Quinacrine
Chloramphenicol	Chlorothiazide	Quinidine
Chlorothiazide	(±) - Chlorpheniramine	Quinine
(±) - Chlorpheniramine	Chlorpromazine	Ranitidine
Chlorpromazine	Chlorquine	Salicylic acid
Chlorquine	Cholesterol	Salicylic acid
Cholesterol	Clomipramine	Secobarbital
Clomipramine	Clonidine	Serotonin
Clonidine	Cocacethylene	Sulfamethazine
Cocacethylene	Cocaine hydrochloride	Sulfindac
Cocaine hydrochloride	Codeine	Temazepam
Codeine	Cortisone	Tetracycline
Cortisone	(-) Cotinine	Tetrahydro-cortisone, 3-acetate
(-) Cotinine	Creatinine	Tetrahydro-cortisone 3-(β-D-glucuronide)
Creatinine	Deoxycorticosterone	Tetrahydrozoline
Deoxycorticosterone	Dextromethorphan	Thebaine
Dextromethorphan	Diazepam	Thiamine
Diazepam	Diclofenac	Thioridazine
Diclofenac	Diffunisal	DL-Tyrosine
Diffunisal	Digoxin	Tolbutamide
Digoxin	Diphenhydramine	Triamterene
Diphenhydramine	EDDP	Trifluoperazine
EDDP	EMDP	Trimethoprim
EMDP	Ecgonine hydrochloride	Trimipramine
Ecgonine hydrochloride	Ecgonine methylester	Trypamine
Ecgonine methylester	(-) -Ψ -Ephedrine	DL-Tryptophan
(-) -Ψ -Ephedrine	[1R,2S] (-) -Ephedrine	Tyramine
[1R,2S] (-) -Ephedrine	(L) - Epinephrine	Uric acid
(L) - Epinephrine	Erythromycin	Verapamil
Erythromycin	β-Estradiol	Zomepirac

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