"DIAQUICK" BAR Dipstick (Barbiturates)

for urine samples

REF	Content
Z02506CE	- 30 tests individually packed (30 x Ref. No: Z02506B). - 1 package insert
Z02506B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	Competitive immunochromatographic assay
Antigen	BAR conjugate,
Antibodies	Anti-BAR monoclonal antibody colloidal gold conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Sensitivity	300 ng/ml

INTENDED USE

The "DIAQUICK" BAR Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Barbiturates in urine at a cut-off concentration of 300 ng/mL of Secobarbital. This test will detect other Barbiturates, please refer to Analytical Specificity table in this package insert. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of Barbiturates leads to tolerance and physical dependence. Short acting Barbiturates taken at 400 mg/day for 2-3 months produces a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death. Only a small amount (less than 5%) of most Barbiturates are excreted unaltered in the urine. The detection period for the Barbiturates in the urine is 4-7 days. The "DIAQUICK" BAR Dipstick (urine) is an easy to use urine test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Barbiturates in urine. The "DIAQUICK" BAR Dipstick (urine) yields a positive result when the Barbiturates in urine exceeds the cut-off level.

TEST PRINCIPLE

The "DIAQUICK" BAR Dipstick (urine) is an immunoassay based on the principle of competitive binding. Drugs that may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Barbiturates, if present in the urine specimen below the cut-off level, will not saturate the binding sites of the antibody in the Dipstick. The antibody coated particles will then be captured by immobilized Barbiturates-protein conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Barbiturates level exceeds the cut-off level because it will saturate all the binding sites of anti-Barbiturates antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- Do not use kit beyond the expiration date.
- Samples may be infectious; properly handle and dispose of all used reaction devices in a biohazard container.

STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability.

SAMPLE COLLECTION AND PREPARATION

Collect a urine sample in clean, dry container, either plastic or glass, without any preservatives. Urine specimens may be refrigerated (2-8°C) and stored up to 2 days. For longer storage, freeze samples (-20°C or below). Bring frozen or refrigerated samples to room temperature before testing. Urine samples exhibiting visible precipitates should be filtered,

centrifuged or allowed to settle. Use only clear aliquots for testing.

ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

- Open the foil pouch by tearing along the notch and remove the Dipstick.
- Dip the stick in the urine specimen submerging only up to the max mark.
- Read the test result at 5 minutes.


Important: Do not interpret a test result after more than 10 minutes. Waiting longer than five minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result at five minutes.

INTERPRETATION OF RESULTS

- Positive Result:**
 A colored band is visible in the Control Zone "C" but not in the Test Zone "T". A positive result indicates that the BAR concentration is at or above the sensitivity of the test.
- Negative Result:**
 A colored band is visible in the Control Zone "C", and in the Test Zone "T". A negative result indicates that the BAR concentration is below the sensitivity of the test.
- Invalid Result:**
 There is no distinct color band visible neither in the Test Zone "T" nor in the Control Zone "C". The test should be considered invalid and it is recommended that the sample be retested.

Note: A very faint line on the test region indicates that the barbiturate in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is NEGATIVE.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

- The "DIAQUICK" BAR Dipstick (urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/ mass spectrometry (GC/MS) is the preferred confirmatory method.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- The test does not distinguish between drugs of abuse and certain medications.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using "DIAQUICK" BAR Dipstick (urine) and a commercially available BAR rapid test. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Other BAR Rapid Test		Total Results
	Positive	Negative	
"DIAQUICK" BAR Dipstick	Results		
	Positive	126	1
	Negative	0	165
Total Results		126	166
% Agreement with this Rapid Test Kit		>99%	>99%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

Method	GC/MS		Total Results
	Positive	Negative	
"DIAQUICK" BAR Dipstick	Results		
	Positive	122	4
	Negative	10	156
Total Results		132	160
% Agreement with GC/MS Analysis		92%	98%

Analytical Sensitivity

A drug-free urine pool was spiked with Secobarbital at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates 96% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Secobarbital Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	27	3
300	Cut-off	30	22	8
375	+25%	30	7	23
450	+50%	30	2	28
600	100%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" BAR Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
Secobarbital	300
Amobarbital	300
Alphenol	150
Aprobarbital	200
Butabarbital	75
Butalbital	2,500
Butethal	100
Cyclopentobarbital	600
Pentobarbital	300
Phenobarbital	100

Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens, containing no Secobarbital, 25% Secobarbital above and below the cut-off, and 50% Secobarbital above and below the 300 ng/mL cut-off was provided to each site. The following results were tabulated:

Secobarbital conc. (ng/mL)	n	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
150	15	13	2	15	0	15	0
225	15	5	10	7	8	10	5
375	15	2	13	5	10	5	10
450	15	0	15	1	14	1	14

Effect of Urinary Specific Gravity

Fifteen (15) urine samples with specific gravity ranging from 1.001 to 1.032 were spiked with 150 ng/mL and 450 ng/mL of Secobarbital respectively. The "DIAQUICK" BAR Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens.

The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Secobarbital to 150 ng/mL and 450 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" BAR Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Barbiturates positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" BAR Dipstick (urine) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

Acetaminophen	Estrone-3-sulfate	Oxolinic acid
Acetophenetidin	Ethyl-p-aminobenzoate	Oxycodone
N-Acetylprocainamide	Fenpropfen	Oxymetazoline
Acetylsalicylic acid	Furosemide	Papaverine
Aminopyrine	Gentisic acid	Penicillin-G
Amitypyline	Hemoglobin	Pentazocine hydrochloride
Amoxicillin	Hydralazine	Perphenazine
Ampicillin	Hydrochlorothiazide	Phencyclidine
L-Ascorbic acid	Hydrocodone	Phenelzine
DL-Amphetamine sulfate	Hydrocortisone	Phentermine
Apomorphine	O-Hydroxyhippuric acid	Trans-2-phenylcyclo-propylamine
Aspartame	p-Hydroxyamphetamine	hydrochloride
Atropine	p-Hydroxy-methamphetamine	L-Phenylephrine
Benzilic acid	3-Hydroxytyramine	β-Phenylethylamine
Benzoic acid	Ibuprofen	Phenylpropanolamine
Benzoyllecgonine	Imipramine	Prednisolone
Benzphetamine	Iproniazid	Prednisone
Bilirubin	(±) - Brompheniramine	Procaine
(±) - Brompheniramine	Caffeine	Promazine
Caffeine	Cannabidiol	Promethazine
Isoxsuprine	Cannabinol	DL-Propranolol
Ketamine	Chloralhydrate	D-Propoxyphene
Ketoprofen	Chlorazepate	D-Pseudoephedrine
DL-Betaxolol	Chlorazepate	Quinacrine
Hydrochlorothiazide	Chlorazepate	Quinidine
Loperamide	(±) - Chlorpheniramine	Quinine
Maprotiline	Chlorpromazine	Rantidine
MDE	Chlorquine	Salicylic acid
Meperidine	Cholesterol	Serotonin
Meprobamate	Chlormipramine	Sulfamethazine
Methadone	Clonidine	Sulindac
(L) Methamphetamine	Cocaehtylene	Methoxyphenamine
Temazepam	Cocaine hydrochloride	(±) - 3,4-Methylenedioxy-amphetamine hydrochloride
Tetracycline	Codeine	(±) - 3,4-Methylenedioxy-methamphetamine hydrochloride
Tetrahydrocortisone, 3-acetate	Cortisone	Morphine-3-β-D glucuronide
Tetrahydrocortisone 3-(β-D-glucuronide)	(-) Cotinine	Morphine Sulfate
(β-D-glucuronide)	Creatinine	Naldixic acid
Tetrahydrozoline	Deoxy corticosterone	Naloxone
Thiamine	Dextromethorphan	Naltrexone
Thioridazine	Diazepam	Naproxen
DL-Tyrosine	Diclofenac	Niacinamide
Tolbutamide	Diffunisal	Nifedipine
Triamterene	Digoxin	Norcodein
Trifluoperazine	Diphenhydramine	Norethindrone
Trimethoprim	Doxylamine	D-Norpropoxyphene
Nitrimpramine	Eggonine hydrochloride	(-) - Ψ-Ephedrine
Tryptamine	Eggonine methylester	[1R,2S] (-) Ephedrine
DL-Tryptophan	(-) - Ψ-Ephedrine	(L) - Epinephrine
Tyramine	[1R,2S] (-) Ephedrine	Ethronycin
Uric acid	(L) - Epinephrine	p-Estradiol
Verapamil	Ethronycin	
Zomepirac	p-Estradiol	

REFERENCES

1. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
3. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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