"DIAQUICK" PCP Dipstick (Phencyclidine) for urine samples

REF

Content

- Z02560CE** - 30 tests individually packed
 (30 x Ref. No: Z02560B).
 - 1 package insert
- Z02560B** - 1 test individually packed,
 - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	competitive immunochromatographic assay
Antigen / Antibodies	PCP conjugate, Anti-PCP monoclonal antibody conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Sensitivity	25 ng/ml

INTENDED USE

The "DIAQUICK" PCP Dipstick (urine) is a rapid qualitative, competitive binding immunoassay for determination of phencyclidine hydrochloride in urine. The test provides only preliminary data which should be confirmed by other methods (such as GC/MS). The test is intended for in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Phencyclidine, commonly known as PCP or "angel dust", is used primarily as a recreational drug for its hallucinogenic effects. It is generally self-administered by intravenous injection or by inhalation and concentrates fastest in fatty tissues and the brain. Excretion patterns vary widely with the individual. Methods historically used for detecting phencyclidine in biological fluids include thin-layer chromatography, gas chromatography, ultraviolet spectroscopy, enzyme immunoassay and radioimmunoassay. While confirmation techniques other than GC/MS may be adequate for some drugs of abuse, GC/MS is generally accepted as a vigorous confirmation technique for all drugs since it provides the best level of confidence in the result. The "DIAQUICK" PCP Dipstick (urine) is an easy, fast and visually read method without the need for instrumentation. The test system employs unique monoclonal and polyclonal antibodies to selectively identify phencyclidine in urine samples with a high degree of sensitivity.

TEST PRINCIPLE

The "DIAQUICK" PCP Dipstick (urine) is a one-step immunoassay in which a chemically labeled drug (drug conjugate) competes with drug which may be present in urine for limited antibody binding sites. In the absence of drug in the urine, a colored antibody-colloid gold conjugate moves along with the sample solution to the immobilized drug conjugate zone on the test band region and is attached to the drug conjugate forming a visible line. Therefore, the formation of a visible precipitant on the test zone occurs when the test urine is negative for the drug. When drug is present in the urine, the drug/metabolite antigen competes with drug conjugate on the test band region for the limited antibody sites on the anti-PCP monoclonal antibody-colloid gold conjugate. When an adequate amount of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored antibody-colloid gold conjugate to the drug conjugate zone on the test band region. Therefore, absence of the color band on the test region indicates a positive result. A control or reference band with a different antigen/antibody reaction is added to the immunochromatographic membrane. It serves as internal control and shows that the test is performed properly.

WARNINGS AND PRECAUTIONS

- Do not use kit beyond the expiration date.
- Samples may be infectious; properly handle and dispose of all used reaction devices in a biohazard container.

STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability.

SAMPLE COLLECTION AND PREPARATION

Collect a urine sample in clean, dry container, either plastic or glass, without any preservatives. Urine specimens may be refrigerated (2-8°C) and stored up to 2 days. For longer storage, freeze samples (-20°C or below). Bring frozen or refrigerated samples to room temperature before testing. Urine samples exhibiting visible precipitates should be filtered, centrifuged or allowed to settle. Use only clear aliquots for testing.

ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

- Open the foil pouch by tearing along the notch and remove the test strip.
- Dip the stick in the urine specimen for 10 to 15 seconds submerging only up to the max mark.
- Place the strip on a non absorbent surface and start the timer.
- Read the test result after 5 minutes.

Important: Do not interpret a test result after more than 10 minutes. Waiting longer than five minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result at five minutes.

INTERPRETATION OF RESULTS

Positive Result:

A colored band is visible in the control line region "C" but not in the test line region "T". A positive result indicates that the PCP concentration is at or above the sensitivity of the test.


Negative Result:

A colored band is visible in the control line region "C", and in the test line region "T". A negative result indicates that the PCP concentration is below the sensitivity of the test.

Invalid Result:

There is no distinct color band visible neither in the test line region "T" nor in the control line region "C". The test should be considered invalid and it is recommended that the sample be retested.

Note: A very faint line on the test region indicates that the phencyclidine in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is **NEGATIVE**.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

- This product is designed for use with human urine only.
- Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
- The test is a qualitative urine assay and is not for determining quantitative concentration levels or the level of intoxication.
- Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" PCP Dipstick (Urine) and a leading commercially available PCP rapid test. Testing was performed on 212 clinical specimens. Ten percent of the specimens employed were either at -25% or +25% level of the cut-off concentration of 25 ng/mL Phencyclidine. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Other PCP Rapid Test		Total Results
	Positive	Negative	
DIAQUICK PCP Dipstick	56	0	56
	1	155	156
Total Results	57	155	212
% Agreement with this Rapid Test Kit	98%	100%	99%

When compared at 25 ng/mL cut-off with GC/MS, the following results

were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
		DIAQUICK PCP Dipstick	Positive	
	Negative	0	157	157
Total Results		50	162	212
% Agreement with GC/MS Analysis		100%	97%	98%

Analytical Sensitivity

A drug-free urine pool was spiked with Phencyclidine at the following concentrations: 0 ng/mL, 12.5 ng/mL, 18.75 ng/mL, 25 ng/mL, 31.25 ng/mL and 37.5 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

PCP Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
12.5	-50%	30	30	0
18.75	-25%	30	19	11
25	Cut-off	30	16	14
31.25	+25%	30	6	24
37.5	+50%	30	0	30

Analytical Specificity

The following table lists the compound that is positively detected in urine by the "DIAQUICK" PCP Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
4-HydroxyPhencyclidine	12,500
Phencyclidine	25

Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Phencyclidine, 25% Phencyclidine above and below the cut-off, and 50% Phencyclidine above and below the 25 ng/mL cut-off was provided to each site. The following results were tabulated:

Phencyclidine conc. (ng/mL)	n	Site A		Site B		Site C	
		-	+	-	+	-	+
		0	15	15	0	15	0
12.5	15	15	0	14	1	14	1
18.75	15	11	4	13	2	10	5
31.25	15	8	7	5	10	1	14
37.5	15	4	11	0	15	0	15

Effect of Urinary Specific Gravity

Fifteen (15) urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 12.5 ng/mL and 37.5 ng/mL of Phencyclidine respectively. The "DIAQUICK" PCP Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Phencyclidine to 12.5 ng/mL and 37.5ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" PCP Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.


Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Phencyclidine positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" PCP Dipstick (urine) at a concentration of 100 µg/mL.

Acetophenetidin	I-Cotinine	Ketamine	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labelalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diffunisal	Methoxyphenamine	Sulfamethazine
I-Ascorbic acid	Digoxin	Methyphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone,
Atropine	β-Estradiol	Niacinamide	3-Acetate
Benzoic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Benzoic acid	Erythromycin	Norethindrone	Tetrahydrozoline
Bilirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Gentisic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloralhydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine	Isoxsuprine	d,l-Propranolol	Zomepirac
Cortisone			

REFERENCES

- Baselt, R.C., Disposition of Toxic Drugs and Chemicals in Man, 2nd Ed., Biomedical Publ., Davis, CA, p 488 (1982).
- Cody, J.T., and Schwarzhoff, R., J. Anal. Toxicol., 17:2630 (1993).
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA) Research Monograph 73. 1986.
- Federal Register, Department of Health and Human Services. Mandatory Guidelines for Federal Workplace Drug Test Programs, 69, 11970-11979, 1988, Federal Register 59, 110, 29908-29931, 1994.
- Dasgupta, A., Saldana, S., Kinnaman, G., Smith, M., and Johansen, K., Clin. Chem., 39 (1):104-108 (1993).
- FDA Guidance for Labeling Urine Drugs of Abuse Screening Testing, Kshitij Mohan, 7/21/1987.
- Froelich, P.E., Gross, G.: Separation and Detection of Phencyclidine in Urine by Gas Chromatography. J. Chromatograph 1977;137:135-143.
- Gupta R.C. et al: Determination of Phencyclidine in Urine and Illicit Street Drug Samples. Clin. Toxicol 1975; 8:611-621.
- Oellerich, M.: Enzyme Immunoassays In Clinical Chemistry: Present Status and Trends. J. Clin. Chem. Biochem., 1980;18:197-208.
- Anilino, Pitts, F.N.: Phencyclidine (PCP): A Review and Perspectives. CRC Crit. Rev. Toxicol 1982;10:145-177.

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