## "DIAQUICK" MDMA Dipstick (Methylenedioxy-methamphetamine) for urine samples

REF	Content
Z05610CE	- 30 tests individually packed (30 x Ref. No: Z05610 B). - 1 package insert
Z05610B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.  
 For diagnosis and therapeutic monitoring only.  
 For use by medical professionals only.**

### GENERAL INFORMATION

<b>Method</b>	competitive immunochromatographic assay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	human urine
<b>Results</b>	within 5 to 10 minutes at room temperature
<b>Sensitivity</b>	500 ng/ml

### INTENDED USE

The "DIAQUICK" MDMA Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Methylenedioxy-methamphetamine (primary ingredient of Ecstasy) in human urine. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

### CLINICAL SIGNIFICANCE

Methylenedioxy-methamphetamine (Ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlander, 1990). The most pervasive effect of MDMA, occurring in virtually all people who have taken a reasonable dose of the drug, is to produce a clenching of the jaws. The "DIAQUICK" MDMA Dipstick (urine) yields a positive result when Methylenedioxy-methamphetamine in urine exceeds 500 ng/mL.

### TEST PRINCIPLE

The "DIAQUICK" MDMA Dipstick (urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Methylenedioxy-methamphetamine, if present in the urine specimen below 500 ng/mL, will not saturate the binding sites of antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized Methylenedioxy-methamphetamine conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Methylenedioxy-methamphetamine level exceeds 500 ng/mL because it will saturate all the binding sites of anti-Methylenedioxy-methamphetamine antibodies. A drug-positive urine specimen will not generate a colored line in the test line region, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

### WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded according to federal, state and local regulations.

### STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability.

### SAMPLE COLLECTION AND PREPARATION

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

### SPECIMEN STORAGE

Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For long-term storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

### ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

1. Open the foil pouch by tearing along the notch and remove the test strip.
2. With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing the strip.
3. Read the result at 5 minutes. **Do not read the results after more than 10 minutes!**


**Important:** Do not interpret a test result after more than 10 minutes. Waiting longer than five minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result at five minutes.

### INTERPRETATION OF RESULTS

1. **Positive Result:**  
 A colored band is visible in the control line region "C" but not in the test line region "T". A positive result indicates that the MDMA concentration is at or above the sensitivity of the test.
2. **Negative Result:**  
 A colored band is visible in the control line region "C", and in test line region "T". A negative result indicates that the MDMA concentration is below the sensitivity of the test.
3. **Invalid Result:**  
 There is no distinct color band visible neither in the test line region "T" nor in the control line region "C". The test should be considered invalid and it is recommended that the sample be retested.

**Note:** A very faint line on the test region indicates that the phencyclidine in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is NEGATIVE.

### QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability.

The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

### LIMITATIONS

1. This product is designed for use with human urine only.
2. Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
3. The test is a qualitative urine assay and is not for determining quantitative concentration levels or the level of intoxication.
4. Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.

### PERFORMANCE CHARACTERISTICS

#### Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" MDMA Dipstick (urine) and a leading commercially available Methylenedioxy-methamphetamine rapid test. Testing was performed on 240 clinical specimens. Ten percent of the specimens employed were either at -25% or +25% level of the cut-off concentration of 500 ng/mL Methylenedioxy-methamphetamine. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method		Other MDMA Rapid Test		Total Results
DIAQUICK MDMA Dipstick	Results	Positive	Negative	
		Positive	90	1
	Negative	0	149	149
Total Results		90	150	240
% Agreement with this Rapid Test Kit		100%	99%	99%

When compared at 500 ng/mL cut-off with GC/MS, the following results were tabulated:

Method		GC/MS				% agreement with GC/MS
		Negative	-25% of Cut-off	+25% of Cut-off	> +25% of Cut-off	
DIAQUICK MDMA Dipstick	Positive	0	3	6	82	97
	Negative	147	2	0	0	100

#### Analytical Sensitivity

A drug-free urine pool was spiked with Methylenedioxyamphetamine at the following concentrations: 0 ng/mL, 250 ng/mL, 375 ng/mL, 500 ng/mL, 625 ng/mL and 750 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

MDMA Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0%	30	30	0
250	-50%	30	30	0
375	-25%	30	26	4
500	Cut-off	30	17	13
625	+25%	30	4	26
750	+50%	30	0	30

#### Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" MDMA Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
(±) 3,4-Methylenedioxyamphetamine HCl	500
(±) 3,4-Methylenedioxyamphetamine HCl	3,000
3,4-Methylenedioxyethyl-amphetamine	300

#### Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing no MDMA, 25% MDMA above and below the cut-off and 50% MDMA above and below the 500 ng/mL cut-off were provided to each site. The results are given below:

MDMA concentration (ng/mL)	n	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
250	15	15	0	15	0	15	0
375	15	15	0	15	0	15	0
625	15	6	9	4	11	7	8
750	15	0	15	0	15	0	15

#### Effect of Urinary Specific Gravity

Fifteen (15) urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 250 ng/mL and 750 ng/mL of Methylenedioxyamphetamine. The "DIAQUICK" MDMA Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

#### Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Methylenedioxyamphetamine to 250 ng/mL and 750 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" MDMA Dipstick (urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

#### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Methylenedioxyamphetamine positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" MDMA Dipstick (urine) at a concentration of 100 µg/mL.

#### Non Cross-Reacting Compounds

4-Acetamidophenol	β-Estradiol	Pentobarbital
Acetophenetidin	Estrone-3-sulfate	Perphenazine
N-Acetylprocainamide	Ethyl-p-aminobenzoate	Phencyclidine
Acetylsalicylic acid	Fenpropfen	Phenazine
Aminopyrine	Furosemide	Phenobarbital
Amitypyline	Gentisic acid	Phentermine
Amobarbital	Hemoglobin	Trans-2-phenylcyclopropylamine hydrochloride
Amoxicillin	Hydralazine	
Ampicillin	Hydrochlorothiazide	
L-Ascorbic acid	Hydrocodone	L-Phenylephrine
D-Amphetamine	Hydrocortisone	β-Phenylethylamine
DL-Amphetamine sulfate	O-Hydroxyhippuric acid	Phenylpropanolamine
L-Amphetamine	p-Hydroxyamphetamine	Prednisolone
Apomorphine	p-Hydroxy-methamphetamine-mine	Prednisone
Aspartame		Procaine
Atropine	3-Hydroxytyramine	Promazine
Benzilic acid	Imipramine	Proprietary
Benzoic acid	Iproniazid	DL-Propranolol
Benzoylcegonine	(±) - Isoproterenol	D-Propoxyphene
Benzphetamine	Isoxsuprine	D-Pseudoephedrine
Bilirubin	Ketamine	Quinacrine
(±) - Brompheniramine	Ketoprofen	Quinidine
Bupropion	Levorphanol	Quinine
Caffeine	Loperamide	Ranitidine
Cannabidiol	Loperamide	Salicicylic acid
Cannabinol	Maprotiline	Secobarbital
Chloralhydrate	Meperidine	
Chloramphenicol	Mephentermine	Serotonin (5-Hydroxytryptamine)
Chlordiazepoxide	Meprobamate	Sulfamethazine
Chlorothiazide	Methamphetamine	Sulindac
(±) - Chlorpheniramine	Methadone	Sustiva
Chlorpromazine	Methoxyphenamine	Tamoxifen
Chlorquine	Methylphenidate	Tetracycline
Cholesterol		
Clomipramine	Morphine-3-β-D-glucuronide	Tetrahydrocortisone, 3- acetate
Clonidine	Morphine sulfate	Tetrahydrocortisone 3-(β-D glucuronide)
Cocaine	Nalidixic acid	Tetrahydrozoline
Cocaine hydrochloride	Naloxone	Thebaine
Codine	Naltrexone	Theophylline
Cortisone	Naproxen	Thiamine
(-) Cotinine	Niacinamide	Trans-2-phenylcyclopropylamine
Creatinine	Nifedipine	Thioridazine
Deoxycorticosterone	Nimesulidate	Tolbutamide
Dextromethorphan	Norcodein	Trazodone
Diclofenac	Norethindrone	DL-Tyrosine
Diazepam	D-Norpropoxyphene	Triamterene
Diffunisal	Noscapine	Trifluoperazine
Digoxin	D,L-Octopamine	Trimethoprim
Dicyclomine	Oxalic acid	Trimipramine
Diphenhydramine	Oxazepam	Tryptamine
5,5 - Diphenylhydantoin	Oxolinic acid	D L-Tryptophan
Doxylamine	Oxycodone	Tyramine
Egonine hydrochloride	Oxymetazoline	Uric acid
Egonine methylester	Papaverine	Verapamil
(-) - ψ-Ephedrine	Penicillin-G	
[1R,2S](-) Ephedrine	Pentazocine-hydrochloride	
(L) - Epinephrine	Erythromycin	Zomepirac

#### REFERENCES

- 1) A Handbook of Drug and Alcohol Abuse, Gail Winger, Third Edition, Oxford Press, 1992, page 146.
- 2) Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- 3) Hawks RL, Chiang CN. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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