“DIAQUICK” CRP Cassette

for human whole blood, serum or plasma samples

REF

Z06021

Content

- 20 tests individually packed + disposable pipette (20x REF Z06021B)
- 20 extraction tubes with buffer
- 1 work station
- 1 package insert

For professional in vitro diagnostic use only.

INTENDED USE

The „DIAQUICK” CRP Cassette is a rapid chromatographic immunoassay for the semi-quantitative detection of C-reactive protein in human whole blood, serum or plasma.

CLINICAL SIGNIFICANCE

The acute-phase protein CRP (C-reactive protein) is a classic marker for the diagnosis of an inflammation in the body. The synthesis of CRP happens in the liver and is also triggered by pre-inflammatory cytokines like Interleukine-6. CRP reaches the site of inflammation via plasma and actively contributes to the innate immune response.

Systemic inflammations are normally accompanied by a drastic increase of the CRP concentration in blood. Mainly bacterial infections, sterile tissue damage, e.g. after a surgery, traumas or extreme physical stress, tumors mainly in metastatic stage, some auto-immune diseases and malign systemic diseases (Hodgkin and non-Hodgkin lymphomas) can be the reason for a CRP increase.

In diagnostics CRP is often used to distinguish between viral and bacterial infections, as the amount of CRP increase is notably higher in bacterial infections. Often CRP is also determined to check for the success of a therapy, e.g. an antibiotic therapy.

PRINCIPLE

The “DIAQUICK” CRP Cassette is an immunochromatographic test, based on two specific antibodies against human CRP. The concentration-dependent formation of test lines allows a rapid semi-quantitative determination of CRP in whole blood, serum or plasma samples. The test sample is diluted with buffer solution and applied to the test cassettes. The sample now moves along the test strip. If the sample contains CRP, it binds to the first anti-CRP antibody, which is bound to colloidal gold of red colour. The red coloured CRP-antibody-gold-complex now migrates together with the sample fluid along the membrane, on which the second anti-CRP antibody is fixed in lines in different concentrations. The CRP-antibody-gold-complex is caught by the immobilized antibodies on the membrane, which results in the formation of red lines. The number of lines depends on the CRP concentration in the sample material. The more CRP is in the sample, the more red lines will appear.

As a functional control, a red control line will appear at the end of the membrane, which states a correct test procedure. The formation of the control line is independent of the CRP concentration of the sample. A failure of the red control line means that the result is invalid. The test should be repeated in this case.

REAGENTS

The test cassette contains CRP antibody conjugated colloid gold and CRP antibody coated on the membrane and buffer with 0.03 % Proclin 300.

MATERIALS PROVIDED

- test cassettes
- package insert
- dropper
- buffer
- extraction tubes
- dropper tips
- workstation

MATERIAL REQUIRED BUT NOT PROVIDED

- specimen collection container
- timer
- centrifuge (for plasma only)

For fingerstick whole blood:

- lancets
- heparinized glass capillary tubes (10 µL) and dispensing bulb

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. Do not use kit beyond the expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations
- Humidity and temperature can adversely affect results.

STORAGE

Store as packaged in the sealed pouch either at room temperature or refrigerated (2 – 30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE**. Do not use after the expiration date.

SAMPLE COLLECTION AND PREPARATION


- The “DIAQUICK” CRP Cassette (whole blood/serum/plasma) can be performed using whole blood (from venipuncture or fingerstick), serum, or plasma.
- To collect **Venipuncture Whole Blood specimens**:
 Collect anti-coagulated blood specimen (sodium or lithium heparin, potassium or sodium EDTA, sodium oxalate, sodium citrate) following standard laboratory procedures.
- To collect **Fingerstick Whole Blood specimens**:
 - Wash the patient’s hand with soap and warm water or clean with an alcohol swab. Allow to dry.

- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube**:
 - Touch the end of the capillary tube to the blood until filled to approx. 10 µL. Avoid air bubbles.
 - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen well (S) of the test cassette.
- Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2 – 8 °C for up to 3 days. For long term storage, specimens should be kept below -20 °C. Whole blood collected by venipuncture should be stored at 2 – 8 °C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

ASSAY PROCEDURE

Allow the test, specimen, buffer, and/or controls to reach room temperature (15 – 30 °C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
2. Place the test cassette on a clean and level surface.
 For **Serum or Plasma or Venipuncture Whole Blood specimens**: Use a disposable pipette to transfer 5 µL of the serum or plasma specimens to the extraction tube with buffer.
 For **Venipuncture Whole Blood specimens**: Use a disposable pipette to transfer 1 drop of whole blood (approximately 10 µL) to the extraction tube with buffer.
 For **Fingerstick Whole Blood specimens**: Fill the capillary tube and transfer approximately 10 µL of fingerstick whole blood specimen to the extraction tube with buffer.
3. Invert the tube several times to mix the specimens well.
4. **Add 3 drops of diluted specimen** to the specimen well and start the timer.
5. Wait for the coloured line(s) to appear. **Read results at 5 minutes**. Do not interpret the result after 8 minutes

INTERPRETATION OF RESULTS

POSITIVE: One coloured line appears in the control line region (C). One to three coloured lines appear in the test line regions T1 – T3.

Test line (T)	Semi-quantitative result
Only T1 appears	CRP level is ≥ 10 mg/L
T1 and T2 appear	CRP level is ≥ 40 mg/L
T1, T2 and T3 appear	CRP level is > 80 mg/L

NEGATIVE: One coloured line appears in the control line region (C). No line appears in the test line regions (T1 – T3). This result means that there was no CRP in the sample or the concentration was below the detectable range of 10 mg/L.

Test line (T)	Semi-quantitative result
No line appears	CRP level is < 10 mg/L

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

EXPECTED VALUES

CRP is a non-specific marker for inflammation and a cardiac risk marker. For ruling out cardiac risks, its expected value is less than 10 mg/L as per AHA. A CRP level above 80 mg/L signifies some other source of inflammation and/or infection.

QUALITY CONTROL

A procedural control is included in the test. A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The "DIAQUICK" CRP Cassette is for in vitro diagnostic use only. This test should be used for the detection of C-reactive protein in whole blood, serum or plasma specimen.
2. The "DIAQUICK" CRP Cassette will only indicate the semi-quantitative level of CRP in the specimen and should not be used as the sole criteria for evaluating cardiac risks or inflammatory conditions.
3. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
4. The presence of rheumatoid factors in the sample will probably interfere the test result.
5. There is a slight possibility that some whole blood specimens with a very high viscosity or which stored for more than 2 days may not run properly on the test device. Repeat the test with a serum or plasma specimen from the same patient using a new test device.
6. The elevated results of CRP in oral contraceptive (OC) users should be reported with caution as The American Physiological Society has recommended further studies on impact of OC use on CRP and inflammatory parameters.
7. High concentrations of CRP may produce a dose hook effect, resulting in false negative results. High dose hook effect has not been observed with this test up to 2000 mg/L of CRP.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The „DIAQUICK“ CRP Cassette has been tested in comparison with a leading commercial CRP EIA test using clinical specimens.

Method	Result	EIA				Total
		Negative		Positive		
		0-10 mg/L	10-40 mg/L	40-80 mg/L	≥ 80 mg/L	
CRP Cassette	0-10 mg/L	99	1	0	0	100
	10-40 mg/L	3	52	1	0	56
	40-80 mg/L	0	3	25	5	33
	≥ 80 mg/L	0	0	3	79	82
Total		102	56	29	84	271
% Relative Accuracy		97.1 %	92.9 %	86.2 %	94.0 %	98.5 %
			99.4 %			

Relative Sensitivity: 99.4 % (96.7 % – 100 %)*

Relative Specificity: 97.1 % (93.9 % - 98.1 %)*

Accuracy: 98.5 % (91.6 % - 99.4 %)*

* 95 % Confidence Interval

Precision

Intra-Assay

Assays were carried out to determine assay reproducibility using replicates of 10 tests in three different runs for each of three lots using CRP specimen levels at 0 mg/L, 10 mg/L, 40 mg/L, 80 mg/L and 200 mg/L. The specimens were correctly identified > 99.9 % of the time.

Inter-Assay


Between-run precision has been determined by using the five CRP specimen levels at 0 mg/L, 10 mg/L, 40 mg/L, 80 mg/L and 200 mg/L of CRP in 3 independent assays. Three different lots of the DIAQUICK™ CRP Cassette have been tested using these specimens. The specimens were correctly identified > 99.9 % of the time.

Interfering Substances

The following substances do not interfere with the test results at the indicated concentrations: Human Albumin at 20 mg/mL, Bilirubin at 6 mg/mL, Hemachrome at 10 mg/mL, Cholesterol at 5 mg/mL and Triglyceride at 15 mg/mL.

BIBLIOGRAPHY

1. Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417.
2. Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet: 980-983.
3. Macy EM, Hayes TE and Tracy RP (1997) Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. Clin. Chem. 43, 52-58.

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LITERATUR

- Morley JJ, Kushner (1982) Serum C-reactive protein levels in disease. In: Kushner I, Volanakis JE, Gewurz H, eds. C-reactive protein and the plasma protein response to tissue injury. Ann. NY Acad. Sci. 389: 406-417.
- Peltola HO (1982) C-reactive protein for rapid monitoring of infections of the central nervous system. Lancet: 980-983.
- Macy EM, Hayes TE and Tracy RP (1997) Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. Clin. Chem. 43, 52-58.