"DIAQUICK"

Rota/Adeno Combi Cassette

A rapid test for the qualitative detection of *Rota-* and *Adenovirus* antigen in human stool samples.

Cat.No.	Content
Z09600CE	- 5 tests individually packed (5x Ref.No. Z09600B) - 5 buffer tubes - 1 package insert
Z09601CE	- 30 tests individually packed (30x Ref.No. Z09600B) - 30 buffer tubes - 1 package insert
Z09600B	- 1 test individually packed - 1 buffer tube - 1 package insert

For professional in vitro diagnostic use only.

GENERAL INFORMATION

Method	sandwich type immunochromatographic assay
Shelf life	18 months from date of production
Storage	2-30°C
Sample	human stool samples
Results	after 5 minutes

SUMMARY

Rotavirus is the primary causative agent of paediatric gastroenteritis and diarrhoea worldwide. The improvement of food, water and hygiene has done nothing to decrease the incidence of rotavirus disease. Almost every child on the planet may get infected by the age of 5. According to studies 900,000 young children around the world die each year from rotaviruses. Many of them get the disease between December and April in the temperate climates of the northern hemisphere. Most of these deaths occur in developing countries. The infection usually begins with a fever. Soon the patient begins to vomit and has a nasty stomach-ache. The vomiting goes away, followed by watery diarrhoea that lasts from 3 to 9 days. Most of the time kids recover with little difficulty. Sometimes, severe dehydration can occur. The extreme dehydration that can be caused by rotaviruses is second only to the dehydration caused by cholera. The infection starts suddenly and lasts for an average of 4 to 6 days. Rotaviruses are extremely contagious. Only very few particles are needed to transmit the viruses. They originate from the stool, but are found throughout the environment wherever young children spend much time, especially during the winter months. They are resistant to disinfectants used to clean surfaces and to anti-bacterial hand-washing agents. Rotavirus particles remain active on human hands for at least 4 hours, on hard dry surfaces for 10 days, and on wet areas for weeks. Untreated, rotavirus infection may result in severe illness with dehydration and disturbances of the body's normal electrolyte balance, especially in babies and preschool children. Rotavirus is the cause of up to 50% of the hospitalized cases of diarrheal illness in infants and young children. Rotavirus induced dehydration is a major cause of infant morbidity in both developed and developing countries and a major cause of infant mortality in developing countries. The highest prevalence of the disease is experienced in temperate climates during the cooler months of the year. In tropical climates, rotavirus infection can occur all year round. The age groups most susceptible to the disease are infants and children.

Adenovirus is the second most common cause of viral gastro-enteritis in Children (10-15%). This virus may also cause respiratory diseases and, depending on the serotype, also diarrhoea, conjunctivitis, cystitis, etc. At least 47 serotypes of adenovirus have been described, all sharing a common hexon antigen. Serotypes 40 and 41 are the ones associated with gastro-enteritis. The main syndrome is diarrhoea that may last between 9 and 12 days associated with fever and vomits.

The DIAQUICK Rota/Adeno Combi Cassette is an in vitro qualitative immunochromatographic assay for the rapid detection of adenovirus and rotavirus antigens in human stool specimen. The test results are intended to aid in the diagnosis of rota- and adenovirus infection and to monitor the effectiveness of therapeutic treatment.

TEST PRINCIPLE

The DIAQUICK Rota/Adeno Combi Cassette is a sandwich solid phase immunochromatographic assay. To perform the test, an aliquot of diluted stool sample is added to the sample well of the test cassette. The sample flows through a pad containing antibodies against adenovirus and rotavirus coupled to red-colored colloidal gold. If the sample contains adenovirus or rotavirus antigens, the antigen will bind to the antibody coated on the colloidal gold particles to form antigen-antibody-gold complexes. These complexes move on the nitrocellulose membrane by capillary action toward the test line region on which adenovirus and rotavirus specific antibodies are immobilized separately. As the complexes reach the test line, they will bind to the antibody corresponding to the virus on the membrane to form of a line. A red control line will always appear in the result window to indicate that the test has been correctly performed and the test device functions properly. If virus is not present or lower than the detection limit of the test, only the control line will be visible. If the control line dose not developed, the test is invalid.

MATERIAL PROVIDED

- Individually pouched test cassettes
- Specimen collection tubes with 1,5 mL dilution buffer each for sample collection and dilution
- Package insert

MATERIAL REQUIRED BUT NOT PROVIDED

- Absorbent tissue paper to prevent solution from splashing.
- stool specimen collection units
- Timer

STORAGE AND STABILITY

Store kit as packaged either at room temperature or refrigerated (2-30°C). Under these conditions the test and the buffer is stable through the expiration date. The test must remain in the sealed pouch containing a desiccant until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

PRECAUTIONS

- For in vitro diagnostic use.
- Wear protective gloves while handling kit components and test specimens.
- Patient specimens and inactivated Positive Control may contain infectious agents and should be handled and disposed of as potential biohazards.
- Do not use kit components beyond expiration date.
- Dispose all used materials in an appropriate container. Treat as potential biohazard.

SPECIMEN COLLECTION, PREPARATION AND STORAGE

Specimen Collection and Storage

Stool samples must be taken as soon as the symptoms appear. Viral particles decrease in number after one week. Stool specimens should be collected in containers that do not contain media, preservatives, animal serum or detergents as any of these additives may interfere with the DIAQUICK Rota/Adeno Combi Cassette.

Specimens may be stored at 2-8°C for 2 days without interfering with the assay performance. For long-term storage of specimens, -20°C or colder is recommended. Repeated freezing and thawing of specimens is not recommended and may cause erroneous results. Do not store specimens in self-defrosting freezers.

Specimen Preparation

- Unscrew the sample bottle, use the attached applicator stick attached on the cap to transfer small piece of stool (4-6 mm in diameter; approx. 50-200 mg) into the sample bottle containing specimen preparation buffer.
- For liquid or semi-solid stools, add 100 µL of stool to the vial with an appropriate pipette.
- Replace the stick in the bottle and tighten securely. Mix stool sample with the buffer thoroughly by shaking the bottle for a few seconds.

ASSAY PROCEDURE

- 1) Bring all reagents, including test device, to room temperature before use.
- 2) Remove the test card from the sealed foil pouch.
- 3) Hold the sample bottle upright with the tip pointing away from the test performer, snap off the tip.
- 4) Hold the bottle in a vertical position over the sample well of the test card, deliver 3 drops (120 -150 µL) of diluted stool sample to the sample well.
- 5) Read the result between 5 - 10 minutes. A strong positive sample may show result earlier.

Note: Results after 10 minutes may not be accurate.

QUALITY CONTROL

- 1) The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
- 2) Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control material, which is not provided with this test kit, may be commercially available.

INTERPRETATION OF TEST RESULTS

For reading the test result, the colored lines that have appeared in the test result window are interpreted.

Positive result for Rotavirus, Adenovirus or both viruses

Rotavirus Positive: A colored line appears in the control line region (C) and another colored line appears in test result line region for Rotavirus (R).


Adenovirus Positive: A colored line appears in the control line region (C) and another colored line appears in test result line region for Adenovirus (A)

Rotavirus and Adenovirus Positive: A colored line appears in the control line region (C) and two other colored lines appear in the test result line regions for Rotavirus (R) and for Adenovirus (A), respectively.


Note: The intensity of the color in test line regions may vary depending on the concentration of aimed substances present in the specimen. Therefore, any shade of color in the test region should be considered positive. Besides, the substances level can not be determined by this qualitative test.

Negative Result

Negative result: One colored line appears in the control line region (C). No line appears in the test line regions Rota or Adeno.

Invalid result

Invalid result: The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

LIMITATIONS

- 1) The test is for qualitative detection of rota- and adenovirus antigen in stool sample and does not indicate the quantity of the antigens.
- 2) The test is for in vitro diagnostic use only.
- 3) The test result should be used only to evaluate with patients with signs and symptoms of the disease. A definitive clinical diagnosis should only be made by the physician after all clinical and laboratory findings have been evaluated.

EXPECTED VALUES

The DIAQUICK Rota/Adeno Combi Cassette detects the presence of adenovirus and rotavirus antigens in stool specimens. Expected values for any given population should be determined for each laboratory. The positivity rate of any given laboratory may vary depending on geographic location, season, and living environment.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity

The DIAQUICK Rota/Adeno Combi Cassette showed concordance with other commercial tests when tested with commercial performance panels.

Reproducibility

Reproducibility of the DIAQUICK Rota/Adeno Combi Cassette was determined using negative, low positive and high positive controls. These samples were tested in replicates of 10 in a blind study by 3 operators working independently in the same laboratory. The agreement of the expected result was 100%.

Cross Reactivity

The DIAQUICK Rota/Adeno Combi Cassette may cross with the rotavirus antigen from monkey and porcine.

LITERATURE

1. Al-Yousif, Y., J. Anderson, C. Chard-Bergstrom, A. Bustamante, M. Muenzenberger, K. Austin, and S. Kapil. 2001. Evaluation of a latex agglutination kit (Virogen Rotatest) for detection of bovine rotavirus in fecal samples. *Clin. Diagn. Lab. Immunol.* 8:496-498
2. Bellinzoni, R. C., J. Blackhall, H. R. Terzolo, A. R. Moreira, N. Auza, N. Mattion, G. L. Micheo, J. L. La Torre, and E. A. Scodeller. 1990. Microbiology of diarrhea in young beef and dairy calves in Argentina. *Rev. Argent. Microbiol.* 22:130-136.
3. Bendali, F., H. Bichet, F. Schelcher, and M. Sanaa. 1999. Pattern of diarrhea in newborn beef calves in southwest France. *Vet. Res.* 30:61-74.
4. Benfield, D. A., I. J. Stotz, E. A. Nelson, and K. S. Groom. 1984. Comparison of a commercial enzyme-linked immunosorbent assay with electron microscopy, fluorescent antibody, and virus isolation for the detection of bovine and porcine rotavirus. *Am. J. Vet. Res.* 45:1998-2002.
5. Chinsangaram, J., G. Y. Akita, A. E. Castro, and B. I. Osburn. 1993. PCR detection of group A bovine rotaviruses in feces. *J. Vet. Diagn. Investig.* 5:516-521.
6. Wood S.R. et al. 1997. Rapid detection and serotyping of adenovirus by direct immunofluorescence. *J. Med. Virol.* 51: 198-201

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BIBLIOGRAPHIE

1. Al-Yousif, Y., J. Anderson, C. Chard-Bergstrom, A. Bustamante, M. Muenzenberger, K. Austin, and S. Kapil. 2001. Evaluation of a latex agglutination kit (Virogen Rotatest) for detection of bovine rotavirus in fecal samples. Clin. Diagn. Lab. Immunol. 8:496-498
2. Bellinzoni, R. C., J. Blackhall, H. R. Terzolo, A. R. Moreira, N. Auza, N. Mattion, G. L. Micheo, J. L. La Torre, and E. A. Scodeller. 1990. Microbiology of diarrhea in young beef and dairy calves in Argentina. Rev. Argent. Microbiol. 22:130-136.
3. Bendali, F., H. Bichet, F. Schelcher, and M. Sanaa. 1999. Pattern of diarrhea in newborn beef calves in southwest France. Vet. Res. 30:61-74.
4. Benfield, D. A., I. J. Stoltz, E. A. Nelson, and K. S. Groon. 1984. Comparison of a commercial enzymelinked immunosorbent assay with electron microscopy, fluorescent antibody, and virus isolation for the detection of bovine and porcine rotavirus. Am. J. Vet. Res. 45:1998-2002.
5. Chinsangaram, J., G. Y. Akita, A. E. Castro, and B. I. Osburn. 1993. PCR detection of group A bovine rotaviruses in feces. J. Vet. Diagn. Investig. 5:516-521.
6. Wood S.R. et al. 1997. Rapid detection and serotyping of adenovirus by direct immunofluorescence. J. Med. Virol. 51: 198-201