DIAQUICK

Rota/Adeno Combi Cassette

An immunochromatographic lateral flow device for the qualitative detection of Rota- and Adenovirus antigens in human faecal samples

REF

Content

- Z09600CE** - 5 Tests individually packed in foil pouches (5x REF Z09600B)
 - 5 Collection Tubes, filled with 1.5 mL Buffer each
 - 1 Package Insert
- Z09601CE** - 30 Tests individually packed in foil pouches (30x REF Z09600B)
 - 30 Collection Tubes, filled with 1.5 mL Buffer each
 - 1 Package Insert

For professional in-vitro diagnostic use only

GENERAL INFORMATION

Method	immunochromatographic antigen detection
Shelf life	18 months from date of production
Storage	2-30 °C
Sample	fresh human faecal samples
Results	after 10 minutes at room temperature

TEST PRINCIPLE

The DIAQUICK Rota/Adeno Combi Cassette is a qualitative, membrane based immunochromatographic lateral flow device for the detection of Rotavirus and Adenovirus antigen in human faecal samples. The membrane is pre-coated with monoclonal antibodies against Rota- and Adenovirus in the test line region. After sample application, the sample reacts with colloidal gold particles coated with anti-Rota- and Adenovirus antibodies stored in the conjugate pad. The mixture migrates upward on the membrane by capillary action. In case of a positive result the specific antibodies present on the membrane will react with the conjugate mixture and generate visibly coloured lines. To serve as a procedural control, a visibly coloured line will always appear in the control line region, indicating that the proper specimen volume was added and correct membrane wicking has occurred.

REAGENT COMPOSITION

Rota/Adeno Combi Cassette

Strip, housed in a plastic cassette, with the following components:

- Nitro-cellulose membrane strip with coated antibodies against Rotavirus and Adenovirus in the test line region and anti-rabbit IgG in the control line region
- Conjugate pad with anti-Rota- and Adenovirus antibodies, as well as rabbit IgG conjugated to colloidal gold particles
- sample and absorbent pad

Rota/Adeno Combi Buffer

- Tris buffer, pH 7.6
- BSA
- Sodium azide < 0.1 %

ADDITIONALLY REQUIRED MATERIAL

- Specimen collection equipment
- Pipettes with disposable tips for 100 µL (for liquid samples)
- General laboratory equipment

REAGENT STABILITY AND STORAGE

Conditions: The cassette must remain in the sealed aluminium pouch until use.

Storage: at 2-30 °C; Do not freeze!

Stability: up to the expiration date, when stored in the sealed pouch

Do not use the test, if the pouch is damaged.


SAMPLE COLLECTION AND PREPARATION

Sample Collection

- Faecal samples must be taken as soon as the symptoms appear. Viral particles decrease in number after one week.
- Collect faeces in a clean, dry specimen collection container. Use a separate collection container for each sample. Frozen samples must be totally thawed and brought to room temperature before testing.
- Faecal specimens should be collected in containers that do not contain media, preservatives, animal serum or detergents as any of these additives may interfere with the test.

Sample Preparation

- Unscrew the cap of the specimen collection tube, then introduce the collection stick into the faecal specimen to pick up a sufficient quantity (approx. 4-6 mm in diameter or 50-200 mg). Do not scoop the faecal sample.
- For liquid stool samples, aspirate the faecal specimen with a pipette and add approx. 100 µL into the faeces collection tube.
- Screw on and tighten the cap on the faeces collection tube.
- Vigorously shake the collection tube to mix the sample with the dilution buffer.

SAMPLE STABILITY AND STORAGE

The DIAQUICK Rota/Adeno Combi Cassette must be used with human faecal samples.

Stability:	at 2-8 °C	undiluted	diluted in buffer
	at -20 °C	2 days	2 days
		2 months	--

The sample storage at room temperature cannot be recommended for this test. Avoid multiple freeze-thaw cycles, as this can cause erroneous results.

TEST PROCEDURE

- Bring reagents and samples to room temperature (18-25 °C). The sensitivity of the test may be reduced at low temperatures.
- Remove the DIAQUICK Rota/Adeno Combi Cassette from the sealed pouch.
- Vigorously shake the collection tube to assure good sample dispersion
- Take off the protective cap and break the tip of the tube with the tip pointed away from the user.
- Hold the bottle vertically and add 3 drops of solution to the specimen well (S). Avoid adding solid particles with the liquid. Start the timer.
- Read the result at 10 minutes after dispensing the sample. False positive results can appear if the test is read after >10 minutes.

Note: if the test does not start to run due to an excess of solid particles in the solution, try to stir the sample in the specimen well (S) with the collection stick.

QUALITY CONTROL

A procedural control is included in the test. A visible coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and a correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

READING AND INTERPRETATION

POSITIVE RESULTS

Rotavirus positive: A red coloured line is visible in the Rotavirus test line region (R) and another one in the control line region (C). This indicates the presence of Rotavirus antigen in the tested sample.

Adenovirus positive: A red coloured line is visible in the Adenovirus test line region (A) and another one in the control line region (C). This indicates the presence of Adenovirus antigen in the tested sample.

Rota & Adeno positive: Red coloured lines are visible in the Rota- and Adenovirus test line regions (A & R) and another line in the control line region (C). This indicates the presence of both Rota- and Adenovirus antigen in the tested sample.

NEGATIVE RESULTS – Only the control line appears.

Only the coloured line in the control line region (C) appears. This indicates that neither Rota- nor Adenovirus antigen could be detected in the tested sample or are present below the test cut-off.

INVALID RESULTS – The control line does not appear.

If the coloured control line does not appear, the result is invalid. Review the package insert and repeat the test with a new device. If the problem persists, discontinue using the test kit and contact your local distributor.

Note: The colour intensity of the test lines will vary depending on the antigen concentration of in the sample. However, neither the quantitative value nor the rate of increase in Rota- or Adenovirus antigen can be determined by this qualitative test.

DIAGNOSTIC IMPLICATIONS

Rotaviruses and adenoviruses are the main causes of severe infectious diarrhoea, especially in children under 5 years of age.

Rotavirus is the most important universal cause of dehydrating infant and childhood gastroenteritis. Annually, it causes approx. 450,000 deaths of children younger than the age of 5. Infection is universal in the first years of life and may be symptomatic or asymptomatic. Severe illness is most common between 6 months and 2 years of age. As gastroenteritis due to rotavirus is not readily distinguished from that caused by other agents on clinical grounds alone, it is important to perform diagnostic tests to decide on the further course of action in patient treatment.

Adenovirus is, after Rotavirus, the second most common cause of gastroenteritis in young children. The virus also causes febrile illness and respiratory tract infections. Adenoviruses are also ubiquitous, with most humans having serological evidence of prior infection by the age of 10. Most infections are mild and self-limited, however diagnosis is useful for the attending physician to decide on the further patient care, especially in combination with the Rotavirus test results.

Besides aiding in the diagnosis of rota- and adenovirus infections, the test results are also intended to monitor the effectiveness of therapeutic treatment.^{1,2}

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens or tests are handled.
- Handle all specimens as if they contained infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for the proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Used tests should be discarded according to local regulations in a proper biohazard container.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

Adenovirus

Adenovirus type 40 1x 10⁵ TCID₅₀/mL
 Adenovirus type 41 2x 10⁴ TCID₅₀/mL

Rotavirus

Rotavirus Wa 1x 10⁴ TCID₅₀/mL

Diagnostic Accuracy

The DIAQUICK Rota/Adeno Combi Cassette was compared to another commercially available rapid test. The following results were achieved:

Rotavirus

Method	Results	Other Rapid Test		Total Results
		Positive	Negative	
Rota/Adeno Combi Cassette	Positive	214	6	220
	Negative	6	824	830
Total Results		220	830	1050

Relative Sensitivity: 97.3 % Relative Specificity: 99.3 %
 Accuracy: 98.9 %

Adenovirus

Method	Results	Other Rapid Test		Total Results
		Positive	Negative	
Rota/Adeno Combi Cassette	Positive	197	6	203
	Negative	7	840	847
Total Results		204	846	1050

Relative Sensitivity: 96.6 % Relative Specificity: 99.3 %
 Accuracy: 98.8 %

Cross Reactivity

The DIAQUICK Rota/Adeno Combi Cassette may cross-react with simian and porcine rotavirus antigen.

Analytical Specificity

The following organisms were tested and did not show any interference with the with the DIAQUICK Rota/Adeno Combi Cassette:

- *Campylobacter jejuni* - *Klebsiella pneumonia* - *Staphylococcus aureus*
- *Candida albicans* - *Listeria monocytogenes* - *Shigella flexneri*
- *Clostridium perfringens A* - *Moracella catarrhalis* - *Shigella sonnei*
- *Citrobacter freundii* - *Neisseria gonorrhoeae* - *Streptococcus dysgalactiae*
- *Enterococcus faecalis* - *Pseudomonas aeruginosa* - *Streptococcus agalactiae*
- *Escherichia coli* - *Staphylococcus epidermidis* - *Streptococcus pyogenes*

EXPECTED VALUES

Expected values for any given population should be determined for each laboratory. The positivity rate of any given laboratory may vary depending on geographic location, season, and living environment.

WASTE MANAGEMENT

Please refer to local legal regulations.

LIMITATIONS

- The test will only qualitatively indicate the presence of rotavirus and adenovirus antigen in human faecal samples. Neither the quantitative value nor the rate of increase in antigen concentration can be determined by this test.
- This test delivers preliminary results only. All results must be confirmed with alternative clinical tests and information available to the attending physician. A definitive clinical diagnosis should only be made by the physician after all clinical and laboratory finding have been evaluated.
- The test devices are for one-time use only. Use a separate cassette for each sample.

LITERATURE

- Mandell et al. Principles and Practice of Infectious Diseases. 8th Ed. Elsevier Saunders, 2015
- Jayoung Kim MD et al. Evaluation of an Immunochromatographic Assay for the Rapid and Simultaneous Detection of Rotavirus and Adenovirus in Stool Samples. Ann Lab Med 2014;34:216-22
- Al-Yousif YJ et al 2001. Evaluation of a latex agglutination kit (Virogen Rotatest) for detection of bovine rotavirus in fecal samples. Clin. Diagn. Lab. Immunol. 8:496-498
- Bellinzoni RC et al 1990. Microbiology of diarrhea in young beef and dairy calves in Argentina. Rev. Argent. Microbiol. 22:130-136.
- Bendali FH 1999. Pattern of diarrhea in newborn beef calves in southwest France. Vet. Res. 30:61-74.
- Benfield DAJ et al 1984. Comparison of a commercial enzyme linked immunosorbent assay with electron microscopy, fluorescent antibody, and virus isolation for the detection of bovine and porcine rotavirus. Am. J. Vet. Res. 45:1998-2002.
- Chinsangaram JGY et al 1993. PCR detection of group A bovine rotaviruses in feces. J. Vet. Diagn. Investig. 5:516-521.
- Wood SR et al 1997. Rapid detection and serotyping of adenovirus by direct immunofluorescence. J. Med. Virol. 51: 198-201

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BIBLIOGRAPHIE

- Mandell et al. Principles and Practice of Infectious Diseases. 8th Ed. Elsevier Saunders, 2015
- Jayoung Kim MD et al. Evaluation of an Immunochromatographic Assay for the Rapid and Simultaneous Detection of Rotavirus and Adenovirus in Stool Samples. Ann Lab Med 2014;34:216-22
- Al-Yousif YJ et al 2001. Evaluation of a latex agglutination kit (Virogen Rotatest) for detection of bovine rotavirus in fecal samples. Clin. Diagn. Lab. Immunol. 8:496-498
- Bellinzoni RC et al 1990. Microbiology of diarrhea in young beef and dairy calves in Argentina. Rev. Argent. Microbiol. 22:130-136.
- Bendali FH 1999. Pattern of diarrhea in newborn beef calves in southwest France. Vet. Res. 30:61-74.
- Benfield DAJJ et al 1984. Comparison of a commercial enzyme linked immunosorbent assay with electron microscopy, fluorescent antibody, and virus isolation for the detection of bovine and porcine rotavirus. Am. J. Vet. Res. 45:1998-2002.
- Chinsangaram JGY et al 1993. PCR detection of group A bovine rotaviruses in feces. J. Vet. Diagn. Investig. 5:516-521.
- Wood SR et al 1997. Rapid detection and serotyping of adenovirus by direct immunofluorescence. J. Med. Virol. 51: 198-201