# "DIAQUICK" HIV 1&2 Ab Cassette

A rapid test for the qualitative detection of HIV-1 and HIV-2 antibodies in human whole blood, serum or plasma

Cat.No.	Content
Z09742CE	- 30 tests individually packed (30 x Ref. No: Z09742B) - Buffer tube; assay diluent sufficient for 30 tests - 1 package insert

For professional in vitro diagnostic use only.

## GENERAL INFORMATION

<b>Method</b>	sandwich type immunochromatographic assay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	human whole blood, serum or plasma
<b>Results</b>	after 10-20 min.

## SUMMARY

HIV (Human Immunodeficiency Virus) is recognized as the etiologic agent of Acquired Immune Deficiency Syndrome (AIDS). The virus is transmitted by sexual contact, exposure to infected blood, certain body fluids or tissues, and from mother to fetus or child during the perinatal period. HIV-1 has been isolated from patients with AIDS and AIDS related complex, and from healthy persons with high potential risk of developing AIDS. Patients with HIV-2 are found primarily in parts of West Africa. Its course is marked by increasing levels of viral replication and the emergence of more virulent viral strains. This process causes the destruction of the immune system. HIV infection is staged by CD4 cell counts and clinical symptoms. Not all people progress through all "stages" and the time frames may also vary greatly from person to person.

HIV-1 and HIV-2 are similar in their morphology, cell tropism, host interaction and generic structure. Serological studies showed that HIV-1 and HIV-2 have a lot of similar epitopes in core antigens but much less so in the envelope antigens.

These clinical diagnostic issues related to HIV are important in the detection of antibodies to HIV1/2 in human plasma or serum by immunoassays. Among the existing immunological methods, the ELISA format and the immunochromatographic format (rapid) to detect antibodies to HIV1/2 are available. Synthetic peptides and recombinant antigens defining the antigenic determinants of HIV proteins are widely used in ELISA procedures for detecting antibodies to HIV in serum. Recently, the use of recombinant antigens for both detector and capture materials markedly increased both the sensitivity and the specificity compared to a commercial rapid kit based on synthetic peptides.

## TEST PRINCIPLE

The DIAQUICK HIV 1&2 Ab Cassette is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 including subtype-O and HIV-2 simultaneously in human serum, plasma or whole blood.

The DIAQUICK HIV 1&2 Ab Cassette contains a membrane strip, which is precoated with recombinant HIV-1 capture antigen (gp41, p24) on test line 1 region and with recombinant HIV-2 capture antigen (gp36) on test line 2 region respectively. The recombinant HIV-1/2 antigen (gp41, p24 and gp36) colloidal gold conjugate and the sample move along the membrane chromatographically to the test region (T) and form a visible line as the antigen-antibody-antigen gold particle complex forms, with a high degree of sensitivity and specificity. The test line and control line in the result window have been clearly labelled "1" for test line 1, "2" for test line 2 and "C" for the control line. Both test lines and the control line in the result window are not visible before applying any sample. The control line is used for procedural control and should always appear if the test procedure is performed correctly.

## INTENDED USE

The DIAQUICK HIV 1&2 Ab Cassette is a rapid, qualitative test for the detection of antibodies to all isotypes (IgG, IgM, IgA) specific to HIV-1 and HIV-2 simultaneously in human serum, plasma or whole blood. The DIAQUICK HIV 1&2 Ab Cassette is intended for professional use and only as an initial screening test. Reactive samples should be confirmed by a supplemental assay such as ELISA or Western Blot test.

## MATERIAL PROVIDED

The DIAQUICK HIV 1&2 Ab Cassette contains the following items:

- 30 individually pouched test cassettes
- Buffer tube; assay diluent sufficient for 30 tests
- Package insert

Active ingredients of main components:

- 1 test strip includes:
  - gold conjugate: recombinant HIV-1 gp41, p24, HIV-2 gp36 antigen-gold colloid (1 ± 0.2 µg)
  - test line 1: recombinant HIV-1 antigen (gp41, p24) (0.625 ± 0.125 µg)
  - test line 2: recombinant HIV-2 antigen (gp36) (0.5 ± 0.1 µg)
  - control line: goat anti-HIV serum (0.75 ± 0.15 µg)
- Assay diluent: 50 mM Tris-HCl buffer, sodium azide (0.02/ w/v)

## MATERIAL REQUIRED BUT NOT PROVIDED

- 20 µL capillary pipette
- lancets
- Micropipette
- Alcohol swab

## STORAGE AND STABILITY

- The test device should be stored at room temperature (2-30°C). Do not freeze the kit or the components.
- After first opening the buffer tube, the assay diluent is stable until the expiration date, if kept at 2-30°C.
- The test device is sensitive to humidity as well as to heat.
- Perform the test immediately after removing the device from the foil pouch.
- Do not use the kit beyond the expiration date.
- The shelf-life of the kits is as indicated on the outer package.
- Check the desiccant for colour change of the humidity indicator and throw away the pouch if the colour indicates saturation.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not re-use the test device.

## WARNINGS AND PRECAUTIONS

### Precautions under the HIV regulations

- Negative results may not have detected recently acquired HIV infection.

### Safety precautions

- Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
- Do not pipette by mouth.
- Do not eat or smoke while handling specimens.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
- The DIAQUICK HIV 1&2 Ab assay diluent contains a proprietary anti-microbial agent which presents no hazard to the user if normal laboratory safety precautions are followed.

### Handling precautions

- Avoid splashing or aerosol formation.
- Do not mix and interchange different specimens.
- Serum or plasma specimens may be stored for 2 weeks at 2-8°C before testing.
- Anticoagulants such as heparin, EDTA and sodium citrate do not affect the test result.
- Use of hemolytic samples, rheumatoid factor-containing samples and lipemic or icteric samples can impair the test result.

### Warning

- For In vitro diagnostic use only. Do NOT RE-USE test device.
- The instructions for use must be followed exactly to get accurate results. Anyone performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
- Care should be taken to avoid contamination of the bottle tips when dropping assay diluent into the sample well.

## SPECIMEN COLLECTION AND STORAGE

### Whole blood

Collection by venipuncture:

- Collect whole blood in a collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture.
- If blood specimens are not immediately tested, they should be stored refrigerated at 2-8°C.
- When stored at 2-8°C, the blood specimens should be used within 3 days.
- Using blood specimens that have been stored for more than 3 days can cause non-specific reactions.

Collection using a lancet:

- Clean the area to be lanced with an alcohol swab.
- Squeeze the end of the fingertip and pierce with a sterile lancet.
- Take a 20µL-capillary pipette, immerse the open end in the blood drop and then release the pressure to draw blood into the capillary pipette to the black line.

### Plasma or Serum

- Plasma: Collect whole blood in a collection tube (containing anticoagulants such as heparin EDTA and sodium citrate) by venipuncture and then centrifuge the blood to get plasma specimens.
- Serum: Collect whole blood in a collection tube (NOT containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture, leave to settle for 30 min. for blood coagulation and then centrifuge the blood to get serum specimens of the supernatant.
- If plasma or serum specimens are not tested immediately, they should be stored refrigerated at 2-8°C. For storage periods longer than 2 weeks, freezing is recommended. The samples should be brought to room temperature prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

## TEST PROCEDURE

- Remove the test device from the foil pouch; place it on a flat, dry surface.
  - Using a capillary pipette: add 20 µL of blood specimen with a 20µL-capillary pipette into the sample well(s). **OR**
  - Using a micropipette: add 10 µL of plasma or serum specimen (20 µL of blood specimen) into the sample well(s).
- Add 4 drops of assay diluent into sample well(s).
- **Caution:** If you do not hold the bottle vertically, it can lead to inaccurate results. **Exactly 4 drops should be added.** Adding 5 or 6 drops (i.e. more than 4 drops) may cause reverse migration phenomenon and/or overall a little reddish unclear background due to limitation of absorbent pad capability.
- As the test begins to work, you will see purple colour moving across the result window in the centre of the test device.

- Interpret test results in 10–20 min. After adding the diluent, read the result after 10 min. but not more than 20 min.

**Caution:** If the test result is not legible after 10 min. due to background colour, read again 20 min. after adding the diluent. Do not read after more than 20 min.

#### INTERPRETATION OF TEST RESULTS (refer to figure)

- A colour band will appear in the left section of the result window to show that the test is working properly. This band is the control line (C).
- Colour bands will appear in the middle and right section of the result window. These bands are test line 2 and test line 1 (2, 1).

#### Negative Result

The presence of only the control line (C) within the result window indicates a negative result.

#### Positive Result

- Two lines appear within the result window, the control line (C) and the test line 1 (1). This indicates a positive result for HIV-1.
- Two lines appear within the result window, the control line (C) and test line 2 (2). This indicates a positive result for HIV-2.
- Three lines appear in the result window, the control line (C), the test line 1 (1) and the test line 2 (2). This indicates a positive result for HIV-1 and/or HIV-2.
  - If the colour intensity of test line 1 is darker than that of test line 2, you can interpret the result as HIV-1 positive.
  - If the colour intensity of test line 2 is darker than that of test line 1, you can interpret the result as HIV-2 positive.

**Caution:** Although a positive result for HIV-1 and HIV-2 in one patient is a rare case, it's possible as there is a homology in the amino acid sequence between HIV-1 and HIV-2. To determine the virus type or diagnose a co-infection accurately, you must perform a confirmatory test such as Western Blot etc.

#### Invalid Result

If the control line (C) is not present within the result window, this indicates an invalid result. The directions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen are re-tested.

#### LIMITATIONS

- Although a positive result may indicate infection with HIV-1 or HIV-2 virus, a diagnosis of AIDS can only be made on clinical grounds, if an individual meets the case definition for AIDS established by the CDC. For samples repeatedly tested as positive, more specific supplemental tests must be performed.
- A negative result does not eliminate the possibility of HIV-1 / HIV-2 infection. The specimen may contain low levels of antibodies to HIV-1 / HIV-2.
- Immunochromatographic testing alone cannot be used to diagnose AIDS even if the antibodies against HIV-1 and/or HIV-2 are present in a patient specimen.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by a physician after all clinical and laboratory findings have been evaluated.
- The whole blood application on seroconversion panels could not be evaluated due to unavailability of such panels made up of whole blood specimens

#### QUALITY CONTROL

The DIAQUICK HIV 1&2 Ab Cassette incorporates the characters "1", "2" and "C" on the device to identify the location of Test Line (HIV-1), Test Line (HIV-2) and Control Line in the result window. These lines are not visible before applying the sample. The Control Line is used for procedural control. A visible Control Line confirms that the diluent has been applied successfully and that the active ingredients on the strip are functional, but is no confirmation that the sample has been properly applied; it is not a positive sample control.

#### PERFORMANCE CHARACTERISTICS

##### Diagnostic sensitivity

HIV type	No. of samples tested	DIAQUICK HIV 1&2 Cassette	
		Negative	Positive
HIV-1	480	0	480
HIV-2	120	0	120
Total	600	0	600 (100%)

##### Diagnostic specificity

Category	No. of samples tested	DIAQUICK HIV 1&2 Cassette	
		Negative	Positive
Healthy donors	500	500	0
	250	250	0
	250	250	0
Possible interferences	500	498	2
Total	1500	1498 (99,87%)	2

##### Possible interferences

The table below shows the results of the DIAQUICK HIV 1&2 Ab Cassette tested on a variety of samples containing possible interfering substances.

Sample type	No. of samples tested	DIAQUICK HIV 1&2 Cassette	
		Negative	Positive
Hospitalized patients	200	198	2
Pregnant women	200	200	0
Potentially cross reacting samples	100	100	0

\* The results were negative for samples containing Anti-HBV(50), Anti-HCV(5), Anti-CMV(20), Anti-EBV(20), Rheumatoid factor(2), Anti-HTLV-I/HTLV-II(3).

##### Seroconversion panels

The sensitivity was evaluated on 30 seroconversion panels in comparison to another CE-marked rapid test. 14 HIV samples and 183 early HIV samples were evaluated. Both tests recognized 14/14 HIV samples and 43/183 early HIV samples.

#### Equivalence Studies

Different specimen types were tested with the DIAQUICK HIV 1&2 Ab Cassette. Used specimen types were serum, venous blood, fingerstick blood, heparin plasma, EDTA plasma and citrate plasma. For each sample type, 25 negative and 25 positive samples were tested and correctly identified by the test.

#### Precision

- **Within-Run:** the reproducibility was determined by testing 5 different replicates of 9 different specimens containing different concentration of antibody. All obtained results were 100% in agreement with the expected results.
- **Between-Run:** the reproducibility was determined by testing in replicates of 9 different specimens containing different concentration of antibody. All obtained results were 100% in agreement with the expected results.
- **Batch-to-Batch Performance:** the reproducibility was determined by testing in replicates of 9 different specimens containing different concentration of antibody with three different lots. All qualitative results obtained were 100% in agreement with the expected results.

#### LITERATURE

1. McCutchan, F.E. (2006). Global epidemiology of HIV. J. Med. Virol. 78, Suppl. 1, S7-S12.
2. Branson, B.M. (2003). Point-of Care Rapid Tests for HIV Antibodies. J. Lab. Med. 27, 288-295.
3. Ferreira Junior, O.C., Ferreira, C., Riedel, M., Widolin, M.G.V., and Barbosa-Junior, A. (2005). Evaluation of Rapid tests for anti-HIV detection in Brazil. AIDS 19 (suppl 4), S70-S75.
4. Constantine, N. (2006). HIV Antibody Assays. HIV InSite Knowledge Base Chapter.
5. Castilla, J., Pachon, I., Gonzalez, M.P., Amela, C., Munoz, L., Tello, O., Noguer, I., de Ory, F., Leon, P., Alonso, M., Gil, E., and Garcia-Saiz, A. (2000). Sero-prevalence of HIV and HTLV in a representative sample of the Spanish population. Epidemiol. Infect. 125, 159-162.
6. Manfredi, R., Calza, L., and Chiodo, F. (2001). HIV disease among immigrants coming to Italy from outside the European Union: a case control study of epidemiological and clinical features. Epidemiol. Infect. 127, 527-533.
7. Lot, F., Semaille, C., Cazein, F., Barin, F., Pinget, R., Pillonel, J., and Desenclos, J.C. (2004). Preliminary results from the new HIV surveillance system in France. Eurosurveillance 9, 10-12.
8. Soriano, V., Gomes, P., Heneine, W., Holguin, A., Doruana, M., Antunes, R., Mansinho, K., Switzer, W.M., Araujo, C., Shanmugam, V., Lourenco, H., Bonzalez-Lahoz, J., Antunes, F. (2000). Human immunodeficiency virus type 2 (HIV-2) in Portugal: clinical spectrum, circulating subtypes, virus isolation, and plasma viral load. J. Med. Virol. 61, 111-116.

#### SIGNS AND SYMBOLS

	In vitro diagnostic medical device		Manufacturer
	Do not re-use		Batch code
	Consult instructions for use		Catalogue number
	Temperature limitation		Contains sufficient for <n> tests
	Use By: Exp		Product Name

**PU**

For professional in vitro diagnostic use only.

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#### LITERATUR

1. McCutchan, F.E. (2006). Global epidemiology of HIV. J. Med. Virol. 78, Suppl. 1, S7-S12.
2. Branson, B.M. (2003). Point-of Care Rapid Tests for HIV Antibodies. J. Lab. Med. 27, 288-295.
3. Ferreira Junior, O.C., Ferreira, C., Riedel, M., Widolin, M.G.V., and Barbosa-Junior, A. (2005). Evaluation of Rapid tests for anti-HIV detection in Brazil. AIDS 19 (suppl 4), S70-S75.
4. Constantine, N. (2006). HIV Antibody Assays. HIV InSite Knowledge Base Chapter.
5. Castilla, J., Pachon, I., Gonzalez, M.P., Amela, C., Munoz, L., Tello, O., Noguez, I., de Ory, F., Leon, P., Alonso, M., Gil, E., and Garcia-Saiz, A. (2000). Seroprevalence of HIV and HTLV in a representative sample of the Spanish population. Epidemiol. Infect. 125, 159-162.
6. Manfredi, R., Calza, L., and Chiodo, F. (2001). HIV disease among immigrants coming to Italy from outside the European Union: a case control study of epidemiological and clinical features. Epidemiol. Infect. 127, 527-533.
7. Lot, F., Semaille, C., Cazein, F., Barin, F., Pinget, R., Pillonel, J., and Desenclos, J.C. (2004). Preliminary results from the new HIV surveillance system in France. Eurosurveillance 9, 10-12.
8. Soriano, V., Gomes, P., Heneine, W., Holguin, A., Doruana, M., Antunes, R., Mansinho, K., Switzer, W.M., Araujo, C., Shanmugam, V., Lourenco, H., Bonzalez-Lahoz, J., Antunes, F. (2000). Human immunodeficiency virus type 2 (HIV-2) in Portugal: clinical spectrum, circulating subtypes, virus isolation, and plasma viral load. J. Med. Virol. 61, 111-116.

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

- 1 Open the package and look for the following:
  - Test cassette, individually pouched
  - Assay Diluent
  - Package insert



- 2 Read carefully the instruction on how to use the DIAQUICK HIV 1&2 Ab Cassette.

- 3 Look at the expiry date at the back of the foil pouch. Use another kit, if expiry date has passed.


## I. Blood specimen (with a lancet)

Take a 20µL-capillary pipette, immerse the open end in the blood drop and then release the pressure to draw blood into the capillary pipette to the black line.

## II. Blood (by venipuncture), Plasma or Serum specimen

Add 10 µL of plasma or serum (20 µL of blood) into the sample well (S) using a micropipette.

Add 4 drops of assay diluent into the sample well (S). Exactly 4 drops should be added.

## Interpretation

Interpret test results in 10–20 min. After adding the diluent, read the result after 10 min. but not more than 20 min. If the test result is not legible after 10 min. due to background colour, read again 20 min. after adding the diluent. Do not read after more than 20 min.

### How to interpret test results:

#### Negative

#### Positive

##### HIV 1 Positive

When 2 lines appear (C, 1)

Strong

Medium

Weak

When 3 lines appear (C, 1, 2)

If line '1' is darker than line '2', it should be interpreted as HIV-1 positive only.

##### HIV 2 Positive

When 2 line appear (C, 2)

Strong

Medium

Weak

When 3 lines appear (C, 1, 2)

If line '2' is darker than line '1', it should be interpreted as HIV-2 positive only.

#### Remark (Both HIV 1 and HIV 2 Positive):

If the band thickness between line 1 and line 2 is so similar to each other, even very rare, it can be both positive for HIV-1 and HIV-2. In this case, confirm test using Western Blot is recommended to know exact virus type.

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