# “DIAQUICK” hCG Midstream (human chorionic gonadotropin) for urine samples

REF	Content
Z09800CE	30 x 1Test - 30 tests individually packed (30 x Ref. No: Z09800B). - 1 package insert
Z09800B	1 Test - 1 test individually packed - 1 package insert

For in vitro diagnostic use only

## GENERAL INFORMATION

<b>Method</b>	sandwich type immunochromatographic assay
<b>Shelf life</b>	18 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	urine
<b>Results</b>	within 2-5 minutes

## TEST PRINCIPLE

The urine sample is moved forward laterally on the test due to capillary forces. In the presence of hCG an immunologic reaction between hCG and anti-hCG antibodies takes place leading to a distinct red colored test line. Other antibodies form the control line. This reaction serves as a proof for the proper use and function of the device. This methodology is named immunochromatography.

## STORAGE

Store originally packed at room temperature (+2 to +30°C). Open the test immediately before use.

## INTENDED USE

In-vitro diagnostic medical device for the detection of hCG in urine.

## SAMPLE MATERIAL

Human urine sample; no pretreatment; no storage of sample;

## ACTIVE INGREDIENT

Monoclonal antibodies against hCG

## PERFORMANCE

Sensitivity: 20 mIU hCG/mL  
 Working range: 20 mIU hCG/mL– 24000 IU/mL  
 Diagnostic sensitivity: >99%  
 Diagnostic specificity: >99%

## ASSAY PROCEDURE


Allow the test device and urine specimen to equilibrate to room temperature (15-30°C) prior to testing.

1. Open the pouch immediately before performing the test
2. Remove the protection cap
3. Take the stick at the assigned handle and hold the tip into the urine stream for some seconds. Alternatively dip the tip for 5 seconds into a urine sample, which was collected before.
4. Put the protection cap over the tip, and lay down the test device on a plain surface.
5. On closer examination, you can see the moisture flowing through the membrane in the test window.
6. Interpret the result after 2-5 minutes. Please note, that the result should be interpreted after 10 minutes at the latest.  
 A test line, which is formed later the 10 minutes after performing the test, does not indicate a positive result.

## INTERPRETATION OF RESULTS

**Pregnant:** Two red colored lines are visible; one line in the round test window, and a second line in the square window

**Not Pregnant:** There is only one red colored line located in the square window

**Note:** If there is no colored line visible after 5 minutes, or if there is only a red colored line visible in the round window, the test is invalid. Please read further information at "reliability" mentioned below.

Both lines can be of different color intensity, which is normal. If you are not sure about the result after 2-5 minutes because of a very faint line, read the result again after 5 minutes. In many cases, the line intensity increases after this additional period, which makes interpretation easier.

## PRECAUTIONS

- Keep out of reach of children
- Do not eat the device
- Do not use the test in case of damaged pouch
- Do not use the device after the expiration date
- Used test can be put to normal household trash
- Additionally required materials: Clock and clean and dry cup for the collection of urine
- Some of the materials used (e.g. antibodies) are potentially infectious but are harmless in case of proper use of the device

## PERFORMANCE CHARACTERISTICS

### Accuracy

The DIAQUICK hCG Midstream Test has a proven reliability of more than 99%.

### Function

Already in the early stages of pregnancy, the body starts to produce the hormone hCG (human Chorionic Gonadotropine). This hormone is released to urine and can be detected with the DIAQUICK hCG Midstream Test. The DIAQUICK hCG Midstream Test has a sensitivity of 20 mIU/mL, and is capable of detecting pregnancy as early as 1 day after the first missed menses.

### Interfering Substances

Generally the consumption of drugs does not affect the test result. Exceptions are drugs containing hCG. These drugs can lead to false results.

## RELIABILITY

The DIAQUICK hCG Midstream Test is very high. Anyway, there may occur unclear results in few cases.

E.g. if the test line in the round window is very faint, the test should be repeated after 3 days. Therefore, use a new device in combination with morning urine. If you are pregnant, two distinctive lines will appear.

In case of a positive result, please contact your doctor!

If the test result is "not pregnant", and your period is not occurring anyway, we recommend to perform another test 3 days later. If this test says "not pregnant" again, but your period is still missing, please contact your doctor!

In case of positive results, please contact your doctor for helpful advices regarding an ideal progress of pregnancy.

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