"DIAQUICK" TRA Dipstick (Tramadol) for Urine Samples

REF	Content
Z10514CE	- 30 tests individually packed (30 x Ref. No: Z10514B) - 1 package insert
Z10514B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	competitive immunochromatographic assay
Antigen	TRA conjugate,
Antibodies	Anti-TRA monoclonal antibody colloidal gold conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Sensitivity	100 ng/ml

INTENDED USE

The "DIAQUICK" TRA Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Tramadol in human urine. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Liquid chromatography/mass spectrometry (LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Tramadol is a synthetic, centrally-acting analgesic, which is mainly used for treating moderate to moderately severe pain. Its mode of action of is not fully understood, yet. At least, two complementary mechanisms seem applicable: Tramadol possesses a low binding affinity to μ -opioid receptors and inhibits the reuptake of norepinephrine and serotonin.

The continuous use of tramadol in large doses may result in tolerance and physical dependency, ultimately resulting in its abuse. After oral administration, tramadol is extensively metabolized. Approximately 30% of the administered dose is excreted unchanged in urine, 60% are excreted as metabolites. The major pathways seem to be N- and O-demethylation, glucuronidation and sulfation in the liver. The half-life of tramadol in serum is 6-7 h in healthy individuals. Tramadol is detectable in urine 2-5 h after consumption.

The therapeutic concentration range in serum is normally 0.1-0.3 $\mu\text{g/ml}$ ⁵. Studies show that tramadol concentration in urine were considerably higher than in serum⁴. Therefore urine samples are very suitable for proof of tramadol intake in therapeutic drug monitoring or forensic toxicology.

TEST PRINCIPLE

The "DIAQUICK" TRA Dipstick (urine) is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Tramadol, if present in the urine specimen below 100 ng/mL, will not saturate the binding sites of the antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized Tramadol conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Tramadol level exceeds 100 ng/mL because it will saturate all the binding sites of anti-Tramadol antibodies. A drug-positive urine specimen will not generate a colored line in the test line region, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- Do not use kit beyond the expiration date.
- To avoid cross contaminations use a new sample container for each sample
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded according to local regulations.

REAGENTS

The test contains mouse monoclonal anti-Tramadol antibody-coupled particles and a Tramadol-protein conjugate. A goat antibody is employed in the control line system.

STORAGE

Store test kit at 2-30°C; **do not freeze**. Refer to the expiration date printed on the sealed pouch for stability. The test must remain in the sealed pouch until use. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION

The urine must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitations should be centrifuged, filtered or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

1. Open the foil pouch by tearing along the notch and remove the test strip. Use it as soon as possible.
2. Dip the stick vertically, with arrows pointing towards the sample in the urine specimen for 10 to 15 seconds submerging only up to the max mark.
3. Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear.
4. The result should be read after 5 minutes.

Important: Do not interpret the result after 10 minutes. Waiting longer than 10 minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result at five minutes.


INTERPRETATION OF RESULTS

NEGATIVE: Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Tramadol concentration is below the detectable level (100 ng/mL).

NOTE: The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Tramadol concentration exceeds the detectable level (100 ng/mL).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test strip. If the problem persists, discontinue using this batch immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The test is intended to be used with human urine samples, only.
2. The "DIAQUICK" TRA Dipstick (urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Liquid chromatography/mass spectrometry (LC/MS) is the preferred confirmatory method.^{2,3}
3. There is a possibility that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
4. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
5. A positive result does not indicate level or intoxication, administration route or concentration in urine.
6. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.

7. The test does not distinguish between drugs of abuse and certain medications.
8. A positive result may be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

Accuracy

Testing was performed on 300 clinical specimens. Presumptive positive results were confirmed by LC/MS. Negative urine specimens were screened initially by Predicate test, 10% negative specimens were confirmed by LC/MS. The following results were tabulated:

When compared at 100 ng/ml with LC/MS, the following results were tabulated:

Method	Results	LC/MS		Total Results
		Positive	Negative	
"DIAQUICK" TRA Dipstick	Positive	136	6	142
	Negative	1	157	158
	Total Results	137	163	300
% Agreement with LC/MS Analysis		99%	96%	97%

Analytical Sensitivity

A drug-free urine pool was spiked with Tramadol to the concentration at $\pm 50\%$ cut-off and $\pm 25\%$ cut-off. The data are summarized below:

TRA Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	90	90	0
50	-50%	90	90	0
75	-25%	90	90	0
100	Cut-off	90	61	29
125	+25%	90	21	69
150	+50%	90	2	88

Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" TRA Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
n-Desmethol-cis-tramadol	195
o-Desmethyl-cis-tramadol	6,250
Cis-tramadol	100
Phencyclidine	100,000
Procyclidine	100,000
d,l-O-Desmethyl venlafaxine	25,000

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Tramadol positive urine. The following compounds were found not to cross-react when tested at a 100 $\mu\text{g/ml}$ concentration in urine.

Non Cross-Reacting Compounds

Acetophenetidin	l-Cotinine	Ketamine	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labelalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diffunisal	Methoxyphenamine	Sulfamethazine
l-Ascorbic acid	Digoxin	Methylphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone,
Atropine	β -Estradiol	Niacinamide	3-acetate
Benzoic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Benzoic acid	Erythromycin	Norethindrone	Tetrahydrozoline
Billirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Gentisic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloral hydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine	Isoxsuprine	d,l-Propranolol	Zomepirac
Cortisone			

REFERENCES

1. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
2. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Davis, CA., 129, 2002.
3. (Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
4. Congress: Ethnic factors: Implications for Drug Therapy and Global Drug Development AGAH Annual Meeting '99, Heidelberg , ALLEMAGNE (07/02/1999) 1999, vol. 37, no 4, pp. 193-206 (27 ref.), pp. 175-183
5. Musshoff F. and Madea B. Fatality due to ingestion of tramadol alone, Forensic Science Int Vol 116:197-199 (2001)

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LITERATUR

- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Davis, CA., 129, 2002.
- (Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse*. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
- Congress: Ethnic factors: Implications for Drug Therapy and Global Drug Development AGAH Annual Meeting '99, Heidelberg, ALLEMAGNE (07/02/1999) 1999, vol. 37, no 4, pp. 193-206 (27 ref.), pp. 175-183
- Musshoff F. and Madea B. Fatality due to ingestion of tramadol alone, Forensic Science Int Vol 116:197-199 (2001)