# "DIAQUICK" LH Cassette (luteinizing hormone)

for urine samples

REF	Content
<b>Z98502CE</b>	- 30 tests individually packed, disposable pipette (30 x Ref. No: Z98502B). - 1 package insert
<b>Z07301CE</b>	- 10 tests individually packed, disposable pipette (10 x Ref. No Z98502B) - 1 package insert
<b>Z98502B</b>	- 1 test individually packed, disposable pipette - 1 package insert

### For in vitro diagnostic use only

#### GENERAL INFORMATION

<b>Method</b>	sandwich type immunochromatographic assay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	human urine
<b>Results</b>	within 5 minutes

#### INTENDED USE

The "DIAQUICK" LH Cassette is a rapid chromatographic immunoassay for the determination of the human Luteinizing Hormone (hLH) surge in the human urine, intended to predict the ovulation.

#### SUMMARY

Human Luteinizing Hormone (hLH) is a glycoprotein hormone produced and secreted by the basophilic cells of the anterior lobe of the pituitary gland. Secretion of LH is stimulated by hypothalamic luteinizing hormone releasing hormone (LRH). Gonadal steroids and inhibin have both positive and negative feedbacks on LH function. LH secretion follows a fluctuating pattern, more pronounced in women than in men. In women, LH levels vary during the menstrual cycle, with low base levels at the beginning (follicular phase) and end (luteal phase) of the cycle. Ovulation is preceded by a surge in LH (mid-cycle-phase). In postmenopausal women, LH levels rise because of lack of negative feedback from ovarian steroids. In men, LH stimulates production of testosterone in the Leydig cells of the testis, and is important for maintaining spermatogenesis. LH levels in children are low, and do not fluctuate before puberty. Determination of LH in women is useful in investigating disturbances in the menstrual cycle, primary and secondary amenorrhea, hirsutism and virilism. The LH surge is used as an indicator of ovulation. In men, assessment of LH is useful in evaluating infertility, hypogonadism, gynaecomastia and dysfunction of the hypothalamic and pituitary glands. While a woman's body always produces a small amount of Luteinizing hormone (LH), there is a sudden increase of this hormone during the middle of the menstrual cycle. This increase of the LH level called LH surge promotes the release of a mature egg from the ovary (ovulation). For most women, ovulation will occur within 24 - 36 hours after the maximum of the LH surge. Immediately after ovulation the egg is ready to be fertilized for a short time (appr. 12 hours). With this test, LH will be determined in urine during 6 days. In this time the cassette will detect if and when a LH surge occurs. The „most fertile days“ of the menstrual period starts after this LH surge.

#### TEST PRINCIPLE

The "DIAQUICK" LH Cassette is a qualitative, chromatographic sandwich immunoassay for the detection of human Luteinizing Hormone (hLH) in urine. When LH is present in the sample, a colored band with a specific antibody-LH-colored conjugate complex will form at the test region on the membrane. The intensity of the test band depends on the amount of LH in the urine sample. On the other hand a light colored band will always appear at the control region. The control line is forming a cut-off reference at a concentration of 35-40 mIU/ml hLH. If the intensity of the test line is equal or higher than the control line the test is positive and indicates a surge of hLH.

#### STORAGE

Store as packaged in the sealed pouch at room temperature (up to 30°C) or refrigerated (2 – 8°C). The kit is stable within the expiration date printed on the labeling. Do not use beyond the expiration date. Do not freeze!

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

#### SAMPLE COLLECTION AND PREPARATION

- Determine the time you will collect your urine. For best results, collect your urine at about the same time each day, for example, between 10:00AM and 8:00PM. Some women have found that their best specimen is after 12 noon. Do not collect your first urine after waking up.
- Reduce your liquid intake approximately 2 hours prior to urine collection.
- Be sure to write down the date, cycle day and time that you collected your urine in the "Test Results Chart" on the back of this instruction.


- You can store your urine for later testing in the day. Your urine can be stored at room temperature for up to 8 hours or in the refrigerator for up to 24 hours. Do not freeze it. For best results, test your urine on the same day that it is collected. If you refrigerate your urine, let it reach room temperature before testing (about 30 minutes). Do not shake the container. If a sediment forms at the bottom of the collection container, allow the sediment to settle. Use only urine from the top of the container.

#### ASSAY PROCEDURE

**Allow the test cassette, urine specimen and/or controls to equilibrate to room temperature (15-30°C) prior to testing.**

1. Determine the day you will begin testing. (See the below section: "WHEN TO START TESTING").
2. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
3. Place the test cassette on a clean and level surface. Hold the dropper vertically and transfer 3 full drops of urine (approx. 100µl) to the specimen well (S) of the test cassette, and then start the timer. Avoid trapping air bubbles in the specimen well (S).

**Wait for the red line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.**

#### INTERPRETATION OF RESULTS

**POSITIVE:** If two lines are visible, and the line in test region (T) is the same as or darker than the one in the control region (C). It indicates that you will probably ovulate in 24-36 hours.

**NEGATIVE:** If two lines are visible, but the line in the test region (T) is lighter than the one in the control region (C), or if there is no line in the test region (T). It indicates that no LH surge has been detected.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

#### QUALITY CONTROL

A positive result means your urine contains the level of LH that would indicate you are beginning to ovulate. You can assume that you will ovulate in 24-36 hours. A negative result means that no LH surge has been detected and you should continue daily testing. Not every woman ovulates mid-cycle, therefore, you may not see a positive result during the first 5 days of testing. Continue testing with the "DIAQUICK" LH Cassette (urine).

#### WHEN TO START TESTING

Determine the Length of Your Menstrual Cycle. our Menstrual Cycle Length is the number of days from the first day of your period (menstrual bleeding) to the last day before your next period starts. Think back over the last few months to decide what your usual cycle length has been.

Circle your usual cycle length on the WHEN TO START CHART below. Select the number directly underneath. Starting the first day of your last period, count ahead the selected number of days on your calendar. This is the day you should begin testing.

#### WHEN TO START CHART

##### Circle Your Usual Cycle Length

21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40
4	5	6	7	8	9	9	9	10	10	10	12	13	14	15	16	17	18	19	

Count ahead this many days STARTING WITH THE FIRST DAY of your LAST period. Note: If you are unsure about your cycle length, you may want to use your shortest cycle length when reading the chart. If you do this, you may need to test for more than 5 days.

Example: My usual cycle length is 28 days. My last period started on the third. The "When to Start Chart" shows that I will count ahead 9 days beginning with the third. When I count 9 days ahead on the calendar, I find that I will collect and test my urine starting on the 11<sup>th</sup>.

**Specimen Calendar**

Sun.	Mon.	Tue.	Wed.	Thu.	Fri.	Sat.
	1	2	(3)	4	5	6
7	8	9	10	<11>	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

( ) First day of your last period

< > Begin testing with the "DIAQUICK" LH Cassette (urine)

**PRECAUTIONS**

- For professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test cassette should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test cassette should be discarded in a proper biohazard container after testing.

**LIMITATIONS**

1. The test works only when the test procedures are precisely followed.
2. Do not reuse the test cassette.
3. For professional *in vitro* diagnostic use.
4. This test may not be used as a form of birth control.
5. The test results should not be affected by pain relievers, antibiotics and other common drugs. Medication containing hCG or LH may affect the test and should not be taken while using the "DIAQUICK" LH Cassette (urine). In addition, the test will not work properly if you are pregnant, in menopause, or are taking birth control pills.
6. Keep out of the reach of children.

**INTERFERING SUBSTANCES**

The "DIAQUICK" LH Cassette (urine) has been tested with commonly known drugs and hormones including FSH (700 mIU/mL), TSH (500 µIU/mL), and hCG (1000 mIU/mL). At the levels tested, none of these substances interfered with the expected test results.

**PERFORMANCE CHARACTERISTICS**

In laboratory studies show the sensitivity of the "DIAQUICK" LH Cassette (urine) is 40 mIU/mL and the accuracy is 98.7%.

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