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For in vitro diagnostic use only

INTENDED USE

The DIAQUICK LH Cassette is a rapid chromatographic immunoassay for the qualitative detection of human luteinising hormone (hLH) in urine, intended to predict ovulation.

SUMMARY

Human luteinising hormone (hLH) is a glycoprotein hormone produced and secreted by the basophilic cells of the anterior lobe of the pituitary gland. Secretion of LH is stimulated by hypothalamic luteinising hormone releasing hormone (LRH). Gonadal steroids and inhibin have both positive and negative feedbacks on LH function. LH secretion follows a fluctuating pattern, more pronounced in women than in men. In women, LH levels vary during the menstrual cycle, with low base levels at the beginning (follicular phase) and end (luteal phase) of the cycle. Ovulation is preceded by a surge in LH (mid-cycle-phase). In postmenopausal women, LH levels rise because of lack of negative feedback from ovarian steroids. In men, LH stimulates production of testosterone in the Leydig cells of the testis, and is important for maintaining spermatogenesis. LH levels in children are low, and do not fluctuate before puberty. Determination of LH in women is useful in investigating disturbances in the menstrual cycle, primary and secondary amenorrhea, hirsutism and virilism. The LH surge is used as an indicator of ovulation. In men, assessment of LH is useful in evaluating infertility, hypogonadism, gynaecomastia and dysfunction of the hypothalamic and pituitary glands. While a woman's body always produces a small amount of Luteinising hormone (LH), there is a sudden increase of this hormone during the middle of the menstrual cycle.

The DIAQUICK LH Cassette detects this increase of the LH level, called LH surge, which promotes the release of a mature egg from the ovary (ovulation). For most women, ovulation will occur within 24-36 hours after the maximum of the LH surge. Immediately after ovulation the egg is ready to be fertilized for a short time (approx. 12 hours). The most fertile days of the menstrual period start with this LH surge.

The test utilises a combination of antibodies including a monoclonal H antibody to selectively detect elevated levels of LH.

Important: The LH surge and ovulation may not occur in all menstrual cycles.

TEST PRINCIPLE

The DIAQUICK LH Cassette is a rapid chromatographic immunoassay for the qualitative detection of human luteinising hormone (hLH) in urine. When LH is present in the sample, a coloured band will form in the test line region on the membrane. On the test line, anti-LH antibodies are coated, which react with the LH in the sample and form a complex with the anti-LH-coated colloidal gold particles (conjugate).

The intensity of the test band depends on the amount of LH in the urine sample. A coloured band will always appear at the control region. The control line is forming a cut-off reference at a concentration of 40 mIU/mL hLH. If the intensity of the test line is equal or higher than the control line the test is positive and indicates a surge of hLH.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30 °C). The kit is stable within the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

MATERIALS PROVIDED

- test cassettes
- disposable pipettes
- package insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

REAGENTS

The test cassette contains anti-LH particles and anti-LH coated on the membrane.

SAMPLE COLLECTION AND PREPARATION

- Determine the optimal time of urine collection. For best results, collect urine at about the same time each day, for example, between 10:00AM and 8:00PM. Some women have found that their best specimen is after 12 noon. Do not collect the first urine after waking up.
- Reduce liquid intake approximately 2 hours prior to urine collection.
- Record the date, cycle day and time of urine collection. See the TEST RESULTS CHART at the end of this package insert.
- Urine can be stored at room temperature for up to 8 hours or at 2-8 °C for up to 24 hours. Do not freeze. For best results, test urine on the same day that it is collected. If refrigerated, let urine reach room temperature before testing. Do not shake the container. If sediment forms at the bottom of the collection container, allow the sediment to settle. Use only urine from the top of the container.

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Humidity and temperature can adversely affect results.
- The test cassette should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test cassette should be discarded according to local regulations.

WHEN TO START TESTING

Determine the length of the subject's menstrual cycle. The menstrual cycle length is the number of days from the first day of the period (menstrual bleeding) to the last day before the next period starts.

Circle the subject's usual cycle length on the chart below. Select the number directly underneath. Starting the first day of the last period, count ahead the selected number of days on a calendar. This is the day you should begin testing.

Menstrual Cycle Length

21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	
6	6	7	7	8	8	9	10	11	12	13	14	15	16	17	18	19	20	21

Days to Count Ahead

Finally, determine the day to start testing. Starting from and including the first day of the last period, count ahead the number of days indicated in the previous step. This is the day on which testing should begin. As a basic guideline, it is recommended to test once a day for five days.

Note: If uncertain about the length of the subject's menstrual cycle, use the shortest menstrual cycle length (21 days) when reading the chart. In this case, it may be necessary to test for more than 5 days.

Example: The subject's usual menstrual cycle length is 28 days. The chart indicates to count ahead 11 days from the subject's last period. The subject's last period started on the 3rd. Starting from and including the 3rd, count ahead 11 days to arrive at the 13th. Urine collection and testing should start on the 13th and proceed through the 17th. (See the Example Specimen Calendar below).

Example Specimen Calendar

Sun	Mon	Tue	Wed	Thu	Fri	Sat
	1	2	(3)	4	5	6
7	8	9	10	11	12	<13>
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			


() First day of the last period

< > Begin testing with the DIAQUICK LH Cassette

ASSAY PROCEDURE

Allow the test cassette, urine specimen and/or controls to equilibrate to room temperature (15-30 °C) prior to testing.

1. Determine the day you will begin testing. (See the above section "WHEN TO START TESTING").
2. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
3. Place the test cassette on a clean and level surface. Hold the dropper vertically and **transfer 3 full drops of urine** (approx. 120 µL) to the specimen well (S) of the test cassette, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
4. Wait for the coloured line(s) to appear. **Read the result at 5 minutes.** Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

POSITIVE: Two lines are visible, but the line in test region (T) is the same as or darker than the one in the control region (C). It indicates probable ovulation in 24-36 hours.

NEGATIVE: Two lines are visible, but the line in the test region (T) is lighter than the one in the control region (C), or there is no line in the test region (T). It indicates that no LH surge has been detected.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

A positive result indicates that the urine specimen contains the level of LH that would indicate the beginning of ovulation. Ovulation should begin in 24-36 hours.

A negative result indicates that no LH surge has been detected and daily testing should continue. Not every woman ovulates mid-cycle, therefore, a positive result may not appear during the selected days of testing. If you see no positive result, continue testing with the DIAQUICK LH Cassette for several more days. If negative results persist, test again during the next menstrual cycle.

LIMITATIONS

1. For professional in vitro diagnostic use only.
2. This test may not be used as a form of birth control.

3. The test results should not be affected by pain relievers, antibiotics and other common drugs. Medication containing hCG or LH may affect the test and should not be taken while using the DIAQUICK LH Cassette. In addition, the test will not work properly if you are pregnant, in menopause, or are taking birth control pills.
4. Keep out of the reach of children.

INTERFERING SUBSTANCES

The DIAQUICK LH Cassette has been tested with commonly known drugs and hormones including FSH (1 000 mIU/mL), TSH (1 000 µIU/mL), and hCG (100 mIU/mL). At the levels tested, none of these substances interfered with the expected test results.

PERFORMANCE CHARACTERISTICS

Laboratory studies show that the sensitivity of the DIAQUICK LH Cassette is 40 mIU/mL and the accuracy is 99.0 %.

QUESTIONS AND ANSWERS

1. Can the DIAQUICK LH Cassette be used to avoid pregnancy?

No, the test should not be used as a form of birth control.

2. Do alcohol or common medications affect the test?

No, but a physician should be consulted if taking any hormonal medication. Also, recent oral contraceptive use, breastfeeding or pregnancy could affect the results.

3. How long will the line remain visible?

The test should be read at 5 minutes for best results. A positive (surge) result will never disappear. The coloured line(s) may become darker and a tinted background may appear after several hours. Some negative results may later display a faint second colour line because of evaporation from the test line region, which prevents complete migration of the test chemicals. Therefore, do not read the result after 10 minutes and discard the test after reading the test.

4. Once a positive result is obtained, when is the best time to have intercourse?
 Ovulation is likely to occur within 24-36 hours. This is the most fertile time. Sexual intercourse within this time frame is advised.

5. Does this test replace the basal body temperature method (BBT)?

The shift in basal body temperature primarily indicates that ovulation has already occurred. The DIAQUICK LH Cassette indicates that ovulation is about to occur.

6. Does a positive result guarantee a pregnancy after intercourse within the fertile days?

There are many factors that can affect the ability to become pregnant. Often it may be necessary to use the test kit for 3-4 months before achieving pregnancy. A physician should be consulted if pregnancy is not achieved after 3-4 months.

7. How soon can pregnancy be confirmed?

Rapid pregnancy tests can provide results as early as the first missed period.

8. What time of the day should I perform the test? Do I need to use first morning urine?

First morning urine is not recommended, but any other time of day is suitable. For best results, try to collect urine at approximately the same time each day.

9. Will heavy intake of liquids affect the result?

Heavy intake of fluids prior to testing will dilute the hormone in urine. Limiting fluid intake for about two hours before collecting urine is suggested.

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