CORMAY IMMUNO-CONTROL I CORMAY IMMUNO-CONTROL II

CONTROL SERUM (LYOPHILISED)

CORMAY IMMUNO-CONTROL I 2 x 3 ml Cat. No 4-288
CORMAY IMMUNO-CONTROL II 2 x 3 ml Cat. No 4-290

CORMAY IMMUNO-CONTROL I (high level-H, low level L) and CORMAY IMMUNO-CONTROL II (high level-H, low level L) sera are lyophilised, pooled human sera which are as a multianalyte control sera when measuring various serum proteins by automatic analysers with turbidimetric capabilities. Assigned values for both the high and low controls vary from lot to lot and are given in values sheet.

RECONSTITUTION

Reconstitute the contents of one vial with 3 ml of distilled water. Allow to stand at room temperature for 20 minutes, gently swirling occasionally to facilitate dissolution. Once reconstituted treat as human serum specimen.

STABILITY AND STORAGE

1. Lyophilised serum remains stable when stored at 2-10°C until expiry date given on the product label. Protect from light.
2. The constituents in reconstituted control serum are stable for 2 weeks at 2-10°C.

NOTES

1. Product for in vitro diagnostic use only.
2. Once reconstituted control sera are ready for use.
3. Remove the cap carefully in order not to lose the powder in the bottle.
4. The controls are based on human serum that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
5. The reagents contain sodium azide (< 0,1%) as a preservative. Avoid contact with skin and mucous membranes.

WASTE MANAGEMENT

Please refer to local legal requirements.

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PZ CORMAY S.A.