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- Chartrand LJ, Seidman EG. Celiac disease is a lifelong disorder. Clin.Invest.Med., Vol. 19, 357-361, 1996
- Cornell HJ. Coeliac disease: A review of the causative agents and their possible mechanisms of action. Amino Acids, Vol. 10, 1-19, 1996
- Cronin CC, Feighery A, Ferriss JB, Liddy C, Shanahan F, Feighery C. High prevalence of celiac disease among patients with insulin-dependent (type I) diabetes mellitus. Am.J Gastroenterol., Vol. 92, 210-2212, 1997
- Jokinen J, Peters U, Maki M, Miettinen A, Collin P. Celiac sprue in patients with chronic oral ucosal symptoms. J Clin.Gastroenterol, Vol. 26, 23-26, 1998
- Taminiau JA. Celiac disease. Curr.Opin.Pediatr., Vol. 8, 483-486, 1996
- Williams CN. Celiac disease: past, present and future. Can.J Gastroenterol., Vol. 11, 647-649, 1997

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF IgA ANTIBODIES TO GLIADIN IN HUMAN BLOOD SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of IgA antibodies to gliadin in blood serum or plasma.

This kit is designed for measurement of antibodies to gliadin in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Celiac disease (CD) or gluten-sensitive enteropathy is a chronic disease characterized by impaired intestinal absorption due to mucosal lesions. The exact etiology of CD is unknown but it is clearly shown that gliadin – the alcohol soluble fraction of wheat gluten – is the toxic agent. High concentrations of antigliadin antibodies (AGA) found in blood, saliva and intestinal secretions are characteristic for untreated celiac patients. These antibodies gradually disappear after gluten exclusion from the patient's diet. AGA testing is a simple and inexpensive method to efficiently select candidates for mucosal biopsy of the duodenal-jejunal junction, the latter method being essential for the confirmation of the diagnosis of CD. Early detection of AGA in high risk populations would contribute to prevent the insidious consequences of chronic malabsorption. Individuals at risk include short-stature children, unexplained anemia, unexplained hypocalcemia or osteomalacia, delayed puberty, insulin-dependent diabetes mellitus, autoimmune thyroiditis and selective IgA deficiency.

It is well established that IgA-AGA are more specific than IgG-AGA. Nevertheless, combined IgG/IgA screening might be more effective since there is an unexplained but clear association between CD and selective IgA deficiency. In treated celiac patients without this deficiency IgA-AGA is the test of choice for monitoring diet compliance. Serum IgA-AGA level responds very quickly to the admission of gluten-free diet (levels drop below the cut-off level within two to six months) while IgG-AGA may take more than one year to become negative; breaking the diet causes more prompt elevation of IgA-AGA compared to IgG-AGA. Herpetiform dermatitis is a disease entity strongly associated to gluten-sensitive enteropathy, and AGA serology is not capable to distinguish between these two diseases.

Separately, there are well reported clinical conditions showing positive IgG-AGA and, rarely, IgA-AGA nonrelated to histologically proven CD, eg all kind of malabsorption syndromes, including Crohn's disease, ulcerative colitis, galactosidase deficiency, post-infection malabsorption etc. The patients with rheumatoid arthritis, Sjogren syndrome, systemic sclerosis and other connective tissue diseases show abnormally high prevalence of moderately elevated gliadin IgA and IgG. These findings may be considered as non-relevant to GI pathology; however, gliadin free diet may be implemented for their treatment.

According to data obtained by XEMA, first degree relatives of CD children patients (especially their parents) show very frequent positivity for both IgG-AGA and IgA-AGA while displaying virtually no symptoms of gastrointestinal diseases, herpetiform dermatitis, or other diseases supposed to cause this positivity.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards species specific Ig, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5. KIT COMPONENTS

5.1. Contents of the Kit

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1 SORB MTP	Gliadin IgG EIA strips, 8x12 wells	1	pcs		until exp. date
2 CAL 1-5	Calibrator set, 1.1 ml each. The set contains 5 calibrators: 0; 25; 50; 100; 200 U/ml	5	pcs	blue (C1 – colourless)	2 months
3 CONTROL	Control serum (1.1 ml)	1	pcs	colourless	2 months
4 CONJ HRP	Conjugate, 11 ml	1	pcs	purple	until exp. date
5 DIL SPE	EIA buffer, 50 ml	1	pcs	blue	until exp. date
6 SUBS TMB	Substrate solution, 11 ml	1	pcs	colourless	until exp. date
7 BUF WASH 21X	Washing solution concentrate 21x, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 1 month at +2...+8 °C or 5 days at RT
8 STOP	Stop solution, 11 ml	1	pcs	colourless	until exp. date
9 N003	Plate sealing tape	2	pcs		N/A
10 K181I	Instruction Gliadin IgA EIA	1	pcs		N/A
11 K181Q	QC data sheet Gliadin IgA EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0.

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute all samples using buffer DIL SPE (EIA buffer) 101 fold. (5 µl of sample + 500 µl of diluent) . Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5, control sample CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +18...+25 °C.
6	Prepare washing solution by 21x dilution of washing solution concentrate BUF WASH 21X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +18...+25 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells
11	Incubate 10-20 minutes at +18...+25 °C
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator
15	Apply point-by-point method for data reduction.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	5 µl of sample + 500 µl of diluent	100	1

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus antibodies to gliadin concentration.

9.3. Determine the corresponding concentration of antibodies to gliadin in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 U/ml	0.04
CAL 2	25 U/ml	0.44
CAL 3	50 U/ml	0.88
CAL 4	100 U/ml	1.58
CAL 5	200 U/ml	2.49

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for Gliadin IgA. Based on data obtained by XEMA, the following normal range is recommended (see below).

Patients with all forms of interstitial lung diseases (hypersensitivity pneumonitis, idiopathic fibrosing alveolitis, fibrotic stages of sarcoidosis, etc.) show very high frequency of Gliadin IgA (ca. 80% of cases) and/or Gliadin IgG (ca. 70% of cases). This phenomenon is not related to concomitant GI pathology and to the activity of lung disease itself and is resistant to gluten-free diet. The clinical role of these findings is unclear.

Sex, age	Units, U/ml	
	Lower limit	Upper limit
< 1 year*	-	35
1-6 yrs	-	46
6-12 yrs	-	40
>12 yrs	-	25

*antibodies of maternal origin

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical sensitivity

Sensitivity of the assay was assessed as being 5 U/ml.

11.2. Linearity

Linearity was checked by assaying dilution series of 5 samples with different antibodies to gliadin concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.3. Recovery

Recovery was estimated by assaying 5 mixed samples with known antibodies to gliadin concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Chartrand LJ, Seidman EG. Celiac disease is a lifelong disorder. Clin. Invest. Med., Vol. 19, 357-361, 1996
2. Cornell HJ. Coeliac disease: A review of the causative agents and their possible mechanisms of action. Amino Acids, Vol. 10, 1-19, 1996
3. Cronin CC, Feighery A, Ferriss JB, Liddy C, Shanahan F, Feighery C. High prevalence of celiac disease among patients with insulin-dependent (type I) diabetes mellitus. Am.J Gastroenterol., Vol. 92, 210-2212, 1997
4. Jokinen J, Peters U, Maki M, Miettinen A, Collin P. Celiac sprue in patients with chronic oral ucosal symptoms. J Clin.Gastroenterol, Vol. 26, 23-26, 1998
5. Taminiua JA. Celiac disease. Curr. Opin. Pediatr., Vol. 8, 483-486, 1996