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1. Schaelike M, Kossakiewicz M, Kossakiewicz A, Schild RL Examination of a first-trimester Down syndrome screening concept on a mix of 11,107 high- and low-risk patients at a private center for prenatal medicine in Germany. // *Ultrasound Obstet Gynecol.* 2009 May;33(5):518-23
2. Wortelboer EJ, Koster MP, Stoutenbeek P, Elvers LH, Loeber JG, Visser GH, Schielen PC. First-trimester Down syndrome screening performance in the Dutch population; how to achieve further improvement? // *Prenat Diagn.* 2009 Mar 17. [Epub ahead of print]
3. Linskens IH, Levitus M, Frans A, Schielen PC, van Vugt JM, Blankenstein MA, Dijkstra HM. Performance of free beta-human chorionic gonadotrophin (free beta-hCG) and pregnancy associated plasma protein-A (PAPP-A) analysis between Delfia Xpress and AutoDelfia systems in The Netherlands. // *Clin Chem Lab Med.* 2009;47(2):222-6
4. Schmidt P, Staboulidou I, Soergel P, Wüstemann M, Hillemanns P, Scharf A. Comparison of Nicolaides' risk evaluation for Down's syndrome with a novel software: an analysis of 1,463 cases. // *Arch Gynecol Obstet.* 2007 Jun;275(6):469-74. Epub 2007 Mar 1
5. Spencer K. Accuracy of Down's syndrome risks produced in a prenatal screening program. // *Ann Clin Biochem.* 1999 Jan;36 (Pt 1):101-3.

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Instruction for use

**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF PREGNANCY-ASSOCIATED PLASMA PROTEIN
A IN HUMAN BLOOD SERUM OR PLASMA**

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of pregnancy-associated plasma protein A in blood serum or plasma.

This kit is designed for measurement of pregnancy-associated plasma protein A in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

PAPP-A (pregnancy-associated plasma protein A) is a high molecular weight glycoprotein consisting of two subunits. In normal pregnancy, PAPP-A level in maternal blood increases during the first two trimesters. Functional significance of PAPP-A during pregnancy remains unclear.

Lowered levels of PAPP-A are observed in Down's syndrome (trisomy 21) during weeks 8-12; after week 14, PAPP-A levels become similar to those in normal pregnancies. Low PAPP-A levels are also found in other trisomies (18 and 13) and chromosomal abnormalities in the fetus and in complicated pregnancies.

Determination of PAPP-A level in the first trimester is used in the following combinations of tests:

1. PAPP-A + free beta-HCG
2. PAPP-A + free beta-HCG + USI of nuchal translucency

In men and non-pregnant women, PAPP-A level is extremely low – usually, it is below the sensitivity level of most immunoassays. Recently, some evidence has appeared to confirm a link between raised PAPP-A levels and increased risk of complications in patients with coronary disease.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human PAPP-A-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal to human PAPP-A, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5,0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP PAPP-A EIA strips, 8x12 wells	1	pcs		until exp. date
2	CAL 1 – 6 Calibrator set, 0.6 ml each. The set contains 6 calibrators: 0; 40; 150; 600; 2500; 10000 mU/l	6	pcs	blue (C1 – colourless)	2 months
3	CONTROL Control serum (0.6 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 11 ml	1	pcs	blue	until exp. date
5	SUBS TMB Substrate solution, 11 ml	1	pcs	colourless	until exp. date
6	BUF WASH 21X Washing solution concentrate 21x, 22 ml	1	pcs	colourless	Concentrate - until exp. date Diluted washing solution - 1 month at 2-8 °C or 5 days at RT
7	STOP Stop solution, 11 ml	1	pcs	colourless	until exp. date
8	N003 Plate sealing tape	2	pcs		N/A
9	K2381 Instruction PAPP-A EIA	1	pcs		N/A
10	K238Q QC data sheet PAPP-A EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100-250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25-250 µl;
- Dry thermostat for 37 °C± 0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0-3.0

5.3. Storage and stability of the Kit

Store the whole kit at 2 to 8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at 2-8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1.** Reagent Preparation

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18 to +25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 21X by 21 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Alternative units:

1 mU/l = 0.0045 µg/ml

7.5. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1-6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dispense 100 µl of CONJ HRP into the wells.
3	Pipet 10 µl of calibrators CAL 1-6, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
4	Incubate 60 minutes at 37 °C.
5	Prepare washing solution by 21x dilution of washing solution concentrate (BUF WASH 21X) by distilled water. Wash the strips 5 times.
6	Dispense 100 µl of SUBS TMB into the wells.
7	Incubate 10-20 minutes at 18-25 °C.
8	Dispense 100 µl of STOP into the wells.
9	Measure OD (optical density) at 450 nm.
10	Set photometer blank on first calibrator.
11	Apply lin-log method for data reduction.

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus pregnancy-associated plasma protein A concentration.

9.3. Determine the corresponding concentration of pregnancy-associated plasma protein A in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 mU/l	0.07
CAL 2	40 mU/l	0.12
CAL 3	150 mU/l	0.20
CAL 4	600 mU/l	0.49
CAL 5	2500 mU/l	1.55
CAL 6	10000 mU/l	2.35

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for PAPP-A. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units		Units alternative	
	mU/l		µg/ml	
	Lower limit	Upper limit	Lower limit	Upper limit
Males	-	150	-	0.68
Females	-	150	-	0.68

Pregnancy, week	Median, mU/l	SKO
9	969	2.9
10	1279	3.3
11	2153	3.4
12	3205	3.4
13	4250	3.6

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity.

The monoclonal antibody is specific for human PAPP-A. There is no cross-reactivity to other species

11.2. Analytical sensitivity

Sensitivity of the assay was assessed as being 20 mU/l.

11.3. Linearity

Linearity was checked by assaying dilution series of 5 samples with different pregnancy-associated plasma protein A concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery

Recovery was estimated by assaying 5 mixed samples with known pregnancy-associated plasma protein A concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

- Schaelike M, Kossakiewicz M, Kossakiewicz A, Schild RL Examination of a first-trimester Down syndrome screening concept on a mix of 11,107 high- and low-risk patients at a private center for prenatal medicine in Germany. // *Ultrasound Obstet Gynecol.* 2009 May;33(5):518-23
- Wortelboer EJ, Koster MP, Stoutenbeek P, Elvers LH, Loeber JG, Visser GH, Schielen PC. First-trimester Down syndrome screening performance in the Dutch population; how to achieve further improvement? // *Prenat Diagn.* 2009 Mar 17. [Epub ahead of print]
- Linskens IH, Levitus M, Frans A, Schielen PC, van Vugt JM, Blankenstein MA, Dijkstra Bloem HM. Performance of free beta-human chorionic gonadotrophin (free beta-hCG) and pregnancy associated plasma protein-A (PAPP-A) analysis between Delfia Xpress and AutoDelfia systems in The Netherlands. // *Clin Chem Lab Med.* 2009;47(2):222-6
- Schmidt P, Staboulidou I, Soergel P, Wüstemann M, Hillemanns P, Scharf A. Comparison of Nicolaides' risk evaluation for Down's syndrome with a novel software: an analysis of 1,463 cases. // *Arch Gynecol Obstet.* 2007 Jun;275(6):469-74. Epub 2007 Mar 1
- Spencer K. Accuracy of Down's syndrome risks produced in a prenatal screening program. // *Ann Clin Biochem.* 1999 Jan;36 (Pt 1):101-3.