



## **INSTRUCTION FOR USE**

### **TEST DEVICE 1000 TESTS**

**REF SA-TD1000**

### **TEST DEVICE 5000 TESTS**

**REF SA-TD5000**

### **TEST DEVICE 10000 TESTS**

**REF SA-TD10000**

**Biota Grup Saglik Sistemleri San. Ve Tic. Ltd. Őti.**  
Izzetpasa Mah. Vefa Poyraz Cad. No:5/5 Sisli, Istanbul- Turkey  
Tel: +90 212 248 52 64-65 Fax:+90 212 291 28 65  
e-mail: [biota@biota.com.tr](mailto:biota@biota.com.tr) web-site: [www.biota.com.tr](http://www.biota.com.tr)



## TEST DEVICE

<b>REF</b>	<b>SA-TD1000</b>	1.000 Tests
<b>REF</b>	<b>SA-TD5000</b>	5.000 Tests
<b>REF</b>	<b>SA-TD10000</b>	10.000 Tests

### 1. INTENDED USE

Dedicated smart card

Single use for the reloading of erythrocyte sedimentation rate tests performed on SABIO Series Automated ESR Analyzers.

### 2. INTRODUCTION AND PRINCIPLE OF THE METHOD

The insertion of Test Device allows the reloading of the instrument for the amount of tests indicated in the **REF**

### 3. CONTENTS

Test Device 1000 Tests (REF: SA-TD1000)	1 Smart Card containing 1000 tests
Test Device 5000 Tests (REF: SA-TD5000)	1 Smart Card containing 5000 tests
Test Device 10000 Tests (REF: SA-TD10000)	1 Smart Card containing 10000 tests

### 4. WARNING

**IMPORTANT!**

Smart card can be used only once. After being used for reloading, the instruments of the SABIO Series Automated ESR Analyzers will not recognise the card again.