



INSTRUCTION FOR USE

TEST DEVICE 1000 TESTS

REF SA-TD1000

TEST DEVICE 5000 TESTS

REF SA-TD5000

TEST DEVICE 10000 TESTS

REF SA-TD10000

Biota Grup Saglik Sistemleri San. Ve Tic. Ltd. Őti.
Izzetpasa Mah. Vefa Poyraz Cad. No:5/5 Sisli, Istanbul- Turkey
Tel: +90 212 248 52 64-65 Fax:+90 212 291 28 65
e-mail: biota@biota.com.tr web-site: www.biota.com.tr



TEST DEVICE

REF	SA-TD1000	1.000 Tests
REF	SA-TD5000	5.000 Tests
REF	SA-TD10000	10.000 Tests

1. INTENDED USE

Dedicated smart card

Single use for the reloading of erythrocyte sedimentation rate tests performed on SABIO Series Automated ESR Analyzers.

2. INTRODUCTION AND PRINCIPLE OF THE METHOD

The insertion of Test Device allows the reloading of the instrument for the amount of tests indicated in the **REF**

3. CONTENTS

Test Device 1000 Tests (REF: SA-TD1000)	1 Smart Card containing 1000 tests
Test Device 5000 Tests (REF: SA-TD5000)	1 Smart Card containing 5000 tests
Test Device 10000 Tests (REF: SA-TD10000)	1 Smart Card containing 10000 tests

4. WARNING

IMPORTANT!

Smart card can be used only once. After being used for reloading, the instruments of the SABIO Series Automated ESR Analyzers will not recognise the card again.