

Liquid Reagents – ready to use

G6PDH LYSE REAGENT

Single Reagent

For the preparation of samples used in the DIALAB G6PDH assay

REF	Cont.	
Y04550B	1 x 900 ml	G6PDH Lyse Reagent
Y04550LV	1 x 60 ml	G6PDH Lyse Reagent
Additionally offered:		
Y04500	10 x 18 ml	G6PDH, UV, kinetic
Y04560	6 x 0.5 ml	G6PDH Control Set
Y04701	5 x 100 ml	Hemoglobin Total Reagent
Y04702	5 x 50 ml	Hemoglobin Total Reagent
Y04705	1 x 2 ml	Hemoglobin Total Calibrator

COMPONENTS:

Sodium Azide
Detergents

CONCENTRATION

<0.1 %

REAGENT STABILITY AND STORAGE

Conditions: close immediately after use
Storage: at 2 - 8°C
Stability: up to the expiration date

TEST PROCEDURE

Please follow the procedure described in the instruction inserts for G6PDH (UV, kinetic) and G6PDH Deficiency Screen.

WARNINGS AND PRECAUTIONS

Take the necessary precautions for the use of laboratory reagents.

WASTE MANAGEMENT

Please refer to local legal requirements.



DIALAB Produktion und Vertrieb von chemisch – technischen Produkten und Laborinstrumenten Gesellschaft m.b.H.
A – 2351 Wiener Neudorf, Austria
IZ-NÖ Süd, Hondastrasse, Objekt M55
Phone: ++43 (0) 2236 660910-0
Fax: ++43 (0) 2236 660910-30 e-mail: office@dialab.at